



Manuals.plus /

› Mega Motion /

› Mega Motion Vance NM7001 Power Recliner Chair User Manual

Mega Motion NM7001

Mega Motion Vance NM7001 Power Recliner Chair User Manual

Model: NM7001 | Brand: Mega Motion

INTRODUCTION

This manual provides essential instructions for the safe and efficient use of your Mega Motion Vance NM7001 Power Recliner Chair. Please read all instructions carefully before assembly and operation. Retain this manual for future reference.

The Mega Motion Vance NM7001 is a heavy-duty, three-position power lift and recline chair designed to assist individuals weighing up to 500 lbs with sitting down and standing up. It features an electric motor for smooth transitions, a split back design, and a chaise seat for comfort.

SAFETY INFORMATION

WARNING: To reduce the risk of injury, ensure all safety precautions are followed.

- Always assemble the chair according to the instructions provided.
- Do not exceed the maximum weight capacity of 500 lbs (227 kg).
- Keep the area around the chair clear of obstacles to prevent tripping hazards.
- Ensure the power cord is not pinched or damaged.
- Keep children and pets away from moving parts.
- Do not allow anyone to sit on the footrest when it is extended.

SETUP AND ASSEMBLY

Your Mega Motion Vance NM7001 Power Recliner Chair requires some assembly. It is recommended that two people assemble the chair due to its weight and size.

1. **Unpacking:** Carefully remove all components from the shipping carton. Inspect for any damage. Report any damage to your dealer immediately.
2. **Base Placement:** Position the chair base in your desired location, ensuring sufficient clearance from

walls and other furniture for reclining and lifting functions.

3. **Backrest Attachment:** Slide the backrest onto the base until it locks securely into place. You may hear an audible click.
4. **Electrical Connections:** Connect the power transformer to the chair's motor and then plug the transformer into a standard wall outlet. Ensure all connections are secure.
5. **Handwand Connection:** Connect the handwand cable to the appropriate port on the chair.



Image: Dimensions of the Mega Motion Vance NM7001 Power Recliner Chair. This diagram illustrates the overall height, width, and depth, which are crucial for planning placement during setup.



Image: The Mega Motion Vance NM7001 Power Recliner Chair in a lifted position, showcasing its ability to assist users in standing.

OPERATING INSTRUCTIONS

The Vance NM7001 chair features an ergonomic push-button control handwand for easy operation of its three positions: upright, reclined, and lifted.

Using the Handwand

- **To Recline:** Press and hold the "Recline" button on the handwand. The chair will smoothly transition from an upright position to a partial or full recline. Release the button at your desired recline angle.
- **To Lift/Stand:** Press and hold the "Lift" button on the handwand. The chair will gently rise and tilt forward, assisting you to a standing position. Release the button when you are comfortably standing.
- **To Return to Upright:** From a reclined position, press and hold the "Lift" button. From a lifted position, press and hold the "Recline" button.
- **USB Charging Port:** A USB port is integrated into the handwand for convenient charging of small electronic devices.



Image: The Mega Motion Vance NM7001 Power Recliner Chair in an upright, seated position.

AS-7001 FEATURES AND BENEFITS

3-Position Reclining Lift Chair

- Three-position, all electric lift and recline
- Heavy duty chaise lounger
- Split back design
- 500 lb maximum weight capacity
- Heavy duty polyester knit fabric
- Furniture grade hardwood/plywood construction
- Counter-balanced recliner mechanism delivers smooth, quiet and easy operation
- Integrated emergency battery backup in case of power failure
- Four durable, easily adjustable steel leg levelers
- Convenient pocket on both sides of chair to store magazines, newspaper or a remote control
- Warranty: 1 year parts and service and lifetime on lift frame

Weight Capacity: 500 lbs

Outside Arm Width: 39.5"

Arm Height from Seat: 6.5"

Seat Width: 27"

Seat to Floor: 20"

Seat Depth: 20.5"

Top of Back to Seat: 27"

Top of Back to Floor: 43.5"

Overall Chair Weight: 143 lbs

Model AS-7001 Fabric Options (Terry)



Forest



Walnut



Bordeaux



Ocean

Image: The Mega Motion Vance NM7001 Power Recliner Chair in a fully reclined position, demonstrating its comfort features.



Image: The Mega Motion Vance NM7001 Power Recliner Chair in a lifted position, assisting a user to stand.

MAINTENANCE AND CARE

Proper maintenance ensures the longevity and optimal performance of your recliner chair.

- **Cleaning:** The polyester fabric upholstery should be spot washed only. Use a mild detergent and a damp cloth to clean spills or stains. Do not saturate the fabric.
- **Inspection:** Periodically check all electrical connections and cables for wear or damage. Ensure the chair's mechanism operates smoothly and quietly.
- **Clearance:** Regularly ensure that the area around the chair is free from obstructions that could interfere with its movement or pose a safety hazard.

TROUBLESHOOTING

If you encounter issues with your chair, refer to the following common troubleshooting steps:

- **Chair Not Responding:**
 - Check if the power cord is securely plugged into the wall outlet and the chair's transformer.
 - Ensure the handwand cable is firmly connected to the chair.
 - Verify that the wall outlet is functional by plugging in another device.

- **Power Failure:**

The chair is equipped with an integrated battery emergency backup system. In the event of a power outage, this system allows you to safely lower the chair to a seated position. The backup batteries are typically located near the motor and should be replaced periodically to ensure functionality.

- **Unusual Noises:**

If the chair makes unusual noises during operation, discontinue use and contact customer support. Do

not attempt to lubricate or repair the mechanism yourself.

SPECIFICATIONS

Model Number	NM7001 (MM-AS-7001-Bordeaux)
Brand	Mega Motion
Maximum Weight Capacity	500 lbs (227 kg)
Reclining Positions	3 (Upright, Partial Recline, Full Recline)
Overall Height	43.5 inches (110.5 cm)
Overall Width	39.5 inches (100.3 cm)
Overall Depth	44 inches (111.8 cm)
Seat Height from Floor	20 inches (50.8 cm)
Seat Depth	20.5 inches (52.1 cm)
Width Between Arms	27 inches (68.6 cm)
Chair Weight	143 lbs (64.9 kg)
Upholstery Material	Polyester Knit Fabric (Textile)
Frame Material	Furniture Grade Hardwood/Plywood, Heavy Duty Steel
Features	Dual Motor Lift, USB Port in Handwand, Integrated Battery Emergency Backup, Side Pockets



Image: Detailed list of AS-7001 features and benefits, including specifications.



Image: Overview of the NM-7001 product features and specifications.

WARRANTY INFORMATION

Your Mega Motion Vance NM7001 Power Recliner Chair comes with a comprehensive warranty:

- **1-Year In-Home Service:** Covers parts and labor for repairs performed in your home for one year from the date of purchase.
- **Lifetime Warranty on Steel Frame:** The heavy-duty steel frame is warranted for the lifetime of the original purchaser.

This warranty covers manufacturing defects and workmanship. It does not cover damage resulting from misuse, abuse, accidents, or unauthorized modifications. Please retain your proof of purchase for warranty claims.

CUSTOMER SUPPORT

For technical assistance, warranty claims, or replacement parts, please contact Mega Motion customer support. Refer to your purchase documentation or the Mega Motion official website for the most current

contact information.

When contacting support, please have your model number (NM7001) and proof of purchase readily available.