

Roscoe Medical BW1832

Roscoe Medical BW1832 Deluxe Wrist Brace User Manual

Model: BW1832

1. INTRODUCTION

The Roscoe Medical BW1832 Deluxe Wrist Brace is designed to provide comfortable support and relief for various wrist conditions, including sprains, strains, and post-cast support. Its ambidextrous design allows it to be worn on either the left or right hand, and it features an adjustable fit for most wrist sizes.

2. SAFETY INFORMATION

Please read all instructions carefully before using this product. If you experience increased pain, numbness, tingling, or skin irritation while wearing the brace, discontinue use and consult a healthcare professional. This product is not a substitute for medical treatment. Always follow your doctor's recommendations.

3. PACKAGE CONTENTS

- 1 x Roscoe Medical BW1832 Deluxe Wrist Brace (Ambidextrous, One Size)

4. SETUP AND APPLICATION

Follow these steps to properly apply the wrist brace:

1. **Identify Brace Orientation:** The Roscoe Medical BW1832 Deluxe Wrist Brace is ambidextrous. It can be worn on either the left or right hand. The removable aluminum support plate should be positioned along the palm side of your wrist.
2. **Open Straps:** Unfasten all hook-and-loop straps on the brace.
3. **Position Hand:** Slide your hand into the brace, ensuring your thumb passes through the designated thumb hole. The aluminum support plate should rest comfortably along the underside of your wrist and palm.
4. **Secure Straps:** Begin by securing the strap closest to your fingers, then proceed to the wrist straps. Adjust each

strap to achieve a snug, yet comfortable fit. The brace should provide firm support without restricting circulation or causing discomfort.



Figure 1: Inside view of the wrist brace, illustrating the soft interior and the pocket for the support plate.



Figure 2: Outside view of the wrist brace, showing the adjustable hook-and-loop straps for a customized fit.

5. OPERATING AND USAGE

Once properly applied, the Roscoe Medical BW1832 Deluxe Wrist Brace provides stability and support to the wrist. The lightweight neoprene and nylon construction offers comfort and flexibility during daily activities. The removable, contoured aluminum support plate can be adjusted or removed based on the desired level of support and mobility. Wear the brace for the duration recommended by your healthcare professional. Ensure the brace remains comfortably snug throughout the day, adjusting straps as needed.

6. MAINTENANCE AND CARE

Proper care will extend the life of your wrist brace:

- **Cleaning:** Hand wash the brace in cold water with mild soap. Do not use bleach or harsh detergents.
- **Drying:** Air dry completely. Do not machine dry or expose to direct heat.
- **Storage:** Store in a cool, dry place away from direct sunlight.
- **Aluminum Support Plate:** The aluminum support plate can be removed before washing if desired. Ensure it is reinserted correctly after the brace is dry.

7. TROUBLESHOOTING

- **Discomfort or Pain:** If the brace causes discomfort or pain, re-adjust the straps to ensure it is not too tight. Check the positioning of the aluminum support plate. If discomfort persists, discontinue use and consult a healthcare

professional.

- **Improper Fit:** Ensure the brace is correctly oriented on your hand and that your thumb is through the thumb hole. Adjust all straps evenly for a secure fit.
- **Loss of Support:** If the brace feels loose or no longer provides adequate support, check the condition of the hook-and-loop fasteners. Ensure the aluminum support plate is correctly inserted.

8. SPECIFICATIONS

Feature	Detail
Model Number	BW1832
Type	Deluxe Wrist Brace
Size	One Size
Hand Compatibility	Ambidextrous (Left or Right)
Material	Lightweight Neoprene & Nylon
Support	Removable, Contoured Aluminum Support Plate
Item Weight	3.2 ounces
Package Dimensions	8 x 4.3 x 1.1 inches

9. WARRANTY AND SUPPORT

For any questions, concerns, or support regarding your Roscoe Medical BW1832 Deluxe Wrist Brace, please contact Roscoe Medical customer service.

Manufacturer Contact Information: 800-376-7263

Please refer to the product packaging for specific warranty details, as terms may vary.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question