

[Manuals.plus](#) /

› [McKesson](#) /

› McKesson Lubricating Jelly User Manual

McKesson B00ZVU6JVI

McKesson Lubricating Jelly User Manual

Model: B00ZVU6JVI

1. PRODUCT OVERVIEW

McKesson Lubricating Jelly is a high-quality, water-based lubricant designed for a variety of medical and personal applications. It comes in sterile, individual 5-gram packets, ensuring hygiene and convenience. This non-greasy, non-staining formula is easy to apply and clean, providing optimal comfort and effectiveness.




Image 1.1: McKesson Lubricating Jelly 144-count box with a single packet.



Image 1.2: Top view of the McKesson Lubricating Jelly packaging.

- Clear • Latex-Free • Odorless
- Sterile • Greaseless



GLYCERINE & WATER BASED

Image 1.3: McKesson Lubricating Jelly box and packet, showing UPC 0612479225008.



Image 1.4: Key characteristics of the lubricating jelly.

2. KEY FEATURES

- **Versatile Use:** Ideal for lubricating rubber materials and synthetic surfaces, including a wide range of medical tools and home devices.
- **Easy Cleanup:** Water-soluble and non-staining formula ensures quick and effortless removal from skin, hair, and synthetic surfaces.
- **Smooth and Uniform Lubrication:** Enhanced viscosity provides consistent application on various medical instruments such as rectal thermometers and vaginal nozzles.
- **Comfortable and Clear Gel:** Non-greasy, fragrance-free, and non-sticky gel designed for maximum comfort without residue or odor.
- **Convenient Packaging:** 5-gram individual packet size is ideal for single use or on-the-go portability.

3. SAFETY INFORMATION

For medical use only. Follow instructions for proper handling and storage. If irritation or other adverse reactions occur, discontinue use and consult a healthcare professional.

Proposition 65 Warning: This product may contain chemicals known to the State of California to cause cancer and birth defects or other reproductive harm.

4. INGREDIENTS

The McKesson Lubricating Jelly is formulated with the following ingredients:

- Water
- Glycerin
- Sodium Hydroxide
- Carbomer 140G
- Polyethylene Glycol
- Propylparaben
- Methylparaben

A very small percentage of men and women may be sensitive to one or more of the ingredients in this lubricant. If irritation is experienced, discontinue use and notify a physician.

5. DIRECTIONS FOR USE (OPERATING)

Apply the lubricating jelly evenly on the medical instrument or surface before use. Use immediately after the product is opened. Do not use if the package is broken.



Image 5.1: Applying jelly from a tube (demonstration of consistency).



Image 5.2: Applying jelly from an individual packet.

6. PRODUCT SPECIFICATIONS

Specification	Detail
Brand	McKesson
Model Number	B00ZVU6JVI
Item Form	Jelly
Scent	Unscented

Material Feature	Greaseless
Number of Items	144 (individual 5g packets)
Product Dimensions	6.4 x 4.9 x 4 inches; 2.08 Pounds (packaging)
Manufacturer	McKesson Brand
Country of Origin	China
Date First Available	April 1, 2014

7. VIDEO DEMONSTRATION

Watch the official product video for an overview of McKesson Water Based Lubricating Jelly features and benefits.

Your browser does not support the video tag.

Video 7.1: McKesson Water Based Lubricating Jelly - Features & Benefits (Duration: 0:28).

8. MAINTENANCE

This product is designed for single-use and does not require maintenance. Ensure proper storage in a cool, dry, and well-ventilated area, away from direct sunlight/fluorescent lighting, at 59-86°F (15-30°C). Store in unopened, undamaged inner packaging to maintain sterility.

9. TROUBLESHOOTING

As a single-use medical lubricant, specific troubleshooting steps are generally not applicable. If you experience any issues related to product quality, packaging integrity, or adverse reactions, please refer to the Safety Information section and contact a healthcare professional or the manufacturer directly.

10. WARRANTY AND SUPPORT

Specific warranty information for this consumable product is not provided in the available data. For product support or inquiries, please contact McKesson Brand customer service. A satisfaction guarantee is indicated on the product packaging.

You can visit the official McKesson Store for more information:[McKesson Store on Amazon](#)

© 2024 McKesson Brand. All rights reserved.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question