

[Manuals.plus](#) /

› [MLAY](#) /

› MLAY Professional RF Beauty Machine, Radio Frequency Face Machine with 2 Probes, at Home Radio Frequency Skin Tightening for Lifting, Wrinkle Reducing, Anti Aging - Salon Effects Gel Included

MLAY MLAY-RF01

MLAY Professional RF Beauty Machine Instruction Manual

Model: MLAY-RF01

1. INTRODUCTION AND PRODUCT OVERVIEW

The MLAY Professional RF Beauty Machine (Model: MLAY-RF01) is designed to provide professional-grade radio frequency (RF) skin tightening treatments for home use. It utilizes bipolar RF technology (1.0MHz) to deliver energy deep into the dermis layers, aiming to stimulate collagen remodeling for visibly lifted, plumper skin and to reduce the appearance of wrinkles and cellulite.

This device is equipped with two probes for targeted treatment on various body and facial areas, ensuring a comprehensive skincare solution.



Figure 1: MLAY Professional RF Beauty Machine with included accessories.

The device generates warmth to promote collagen production. The maximum power output is 25W for the smaller probe and 50W for the larger probe, ensuring effective penetration for skin tightening and rejuvenation.



Figure 2: Illustration of RF energy penetrating the dermis to stimulate collagen.

2. SAFETY INFORMATION AND WARNINGS

Important: Read all instructions and warnings before use. Failure to follow these instructions may result in injury or device malfunction.

Contraindications:

- Do not use if you have a pacemaker or other implanted electronic devices.
- Do not use if you are pregnant or breastfeeding.
- Avoid use on areas with recent surgical procedures. It is recommended to wait at least three months after facial surgeries before using the device on the treated area.
- Do not use on broken, irritated, or infected skin.

- Ensure sufficient conductive gel is applied to avoid burns. People's sensation to temperature varies; choose a lower energy level or move the probe faster to prevent possible burns.

General Safety Precautions:

- Keep out of reach of children.
- Do not immerse the device in water or any other liquid.
- Disconnect the device from the charger after use and before cleaning.
- Clean the device with a soft cloth.

3. PACKAGE CONTENTS

Upon unboxing your MLAY Professional RF Beauty Machine, please ensure all the following items are present:

- MLAY RF Beauty Instrument (Main Unit)
- Small RF Probe (25W) for delicate areas like eyes, nose, and mouth.
- Large RF Probe (50W) for broader areas like cheeks, jawline, abdomen, arms, buttocks, and thighs.
- Power Adapter
- Conductive Gel (100g)
- Instruction Manual

4. SETUP

Follow these steps to set up your MLAY RF Beauty Machine:

1. Unpack all components from the box.
2. Connect the power adapter to the main unit and plug it into a power outlet.
3. Select the desired probe (small or large) and screw it tightly onto the handle of the main unit.



Figure 3: Visual guide for product setup and initial application.

5. OPERATING INSTRUCTIONS

For optimal results and safety, follow these operating guidelines:

1. Clean and dry your skin thoroughly before use.
2. Apply a thin, even layer of the provided conductive gel to the treatment area. Ensure enough gel is applied to allow the probe to glide smoothly.
3. Turn on the device. The display will show '00' minutes.
4. Press the 'TIME' button to select your desired treatment duration (10, 20, or 30 minutes).
5. Press the 'LEVEL' button to select the energy intensity (Level 1, 2, or 3). Start with Level 1 for your first few sessions, especially if you have sensitive skin, and gradually increase as comfortable.
6. Place the probe on the skin with full contact. Move the probe slowly in circular motions over the treatment area. For

facial areas, move upwards to lift the skin.

7. Reapply gel as needed during the treatment to maintain conductivity and prevent discomfort.
8. The device will automatically shut off once the selected time has elapsed.

Treatment Areas and Recommended Usage:

The MLAY RF01 is suitable for various areas of the face and body:

- **Face:** Use the smaller probe for delicate areas like around the eyes (to fade crow's feet), nose, and mouth (to reduce smile lines and nasolabial folds). Use the larger probe for cheeks and jawline to tighten sagging skin.
- **Body:** Use the larger probe for areas such as the abdomen, arms, buttocks, and thighs to improve skin firmness and reduce cellulite.



Figure 4: Visual results of fine line and wrinkle reduction.



Figure 5: Visual results of jawline tightening and lifting.

6. MAINTENANCE AND CARE

Proper maintenance ensures the longevity and effectiveness of your MLAY RF Beauty Machine:

- After each use, disconnect the device from the power outlet.
- Clean the probes and the main unit with a soft, damp cloth. Do not use abrasive cleaners or immerse the device in water.
- Store the device in a dry place, away from direct sunlight and extreme temperatures.

7. TROUBLESHOOTING

If you experience any issues with your MLAY RF Beauty Machine, please refer to the troubleshooting section in the full instruction manual or contact customer support for assistance.

8. SPECIFICATIONS

Specification	Value
Model Number	MLAY-RF01
Brand	MLAY
ASIN	B01IBSX35Q
UPC	756040275907
Item Weight	9.26 Pounds
Package Dimensions	11 x 10.5 x 8.5 inches
Color	Black
Skin Type	Normal

9. WARRANTY AND SUPPORT

MLAY products are designed for quality and durability. For warranty information, technical support, or any inquiries regarding your MLAY Professional RF Beauty Machine, please refer to the contact details provided in your product packaging or visit the official MLAY store on Amazon.
Visit the MLAY Store: [MLAY Amazon Store](#)

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question