

[Manuals.plus](#) /

› [Drive Medical](#) /

› Drive Medical Nitro Elite CF Foldable Rollator Walker with Seat, Black

Drive Medical RTL10266CF

Drive Medical Nitro Elite CF Foldable Rollator Walker with Seat

Model: RTL10266CF | Brand: Drive Medical

INTRODUCTION

Experience the ultra-lightweight elegance of Drive Medical's Nitro Carbon Fiber portable rollator walker with a seat. This exceptional rollator features hand craftsmanship, thoughtful features, and the finest materials, designed for ease of use and stability across various terrains.

SETUP INSTRUCTIONS

Unpacking and Initial Setup:

1. Carefully remove the Nitro Elite CF Rollator from its packaging. The unit comes fully assembled, requiring no tools for initial setup.
2. Unfold the rollator by gently pulling the seat grip up and allowing the frame to expand. Ensure the cross-brace design is fully extended for added stability.
3. Verify that the plastic clip holds the frame securely in the open position.

Navigate with confidence indoors & outdoors



Large 10" front wheels!

Figure 1: The rollator easily folds for transport and storage.

Comes fully assembled



Supports users up to 300 lb

Figure 2: Key dimensions of the assembled rollator.

Handle Height Adjustment:

1. Locate the push-button height adjustment mechanism on each handle.
2. Press the button and slide the handle up or down to your desired height. Ensure both handles are adjusted to the same height for balanced support.
3. Release the button and ensure the handle locks securely into place.

Personalized comfort with adjustable height handles



7-level adjustments!

Figure 3: Adjusting the handle height for personalized comfort.

OPERATING INSTRUCTIONS

Walking with the Rollator:

- Stand upright with the rollator in front of you, holding both handles firmly.
- Push the rollator forward a comfortable distance, then step forward with one foot, followed by the other. Maintain a steady pace.
- The large 10-inch front casters provide optimal steering and smooth navigation over various terrains, both indoors and outdoors.

Comes fully assembled



Supports users up to 300 lb

Figure 4: The rollator provides steady support for indoor and outdoor use.

Your browser does not support the video tag.

Video 1: Overview of the Drive Medical Nitro Elite CF Carbon Fiber Walker Rollator, demonstrating key features and ease of use.

Your browser does not support the video tag.

Video 2: Detailed demonstration of the Drive Medical Nitro Elite Carbon Fiber Rollator Walker's functionality.

Using the Brakes:

- To slow down or stop, gently pull the brake levers located on the handles towards you.
- To engage the parking brake (for sitting or stability), push the brake levers down until they lock into place. This will prevent the rollator from moving.
- To release the parking brake, pull the levers back up.
- The hidden brake cable ensures maximum safety and a sleek appearance.

Using the Seat:

- Ensure the parking brakes are engaged before attempting to sit.
- The luxury seat is made of easy-to-clean nylon and provides a comfortable resting spot.
- The removable back is also height-adjustable for added comfort.

Supportive & comfortable resting with foam seat & back



Figure 5: The padded seat and backrest provide supportive and comfortable resting.

Storage Bag:

- The rollator includes a spacious, removable zippered storage bag, highlighted with attractive contrast diamond stitching.
- A unique attachment keeps the storage bag securely in place when the rollator is open or folded.



Secure essentials with removable zippered storage bag

Figure 6: Securely store essentials in the removable zippered bag.

MAINTENANCE

Cleaning:

- Clean the frame with a damp cloth and mild detergent. Avoid abrasive cleaners.
- The nylon seat can be wiped clean with a damp cloth.
- Regularly clean the wheels to prevent dirt and debris buildup, which can affect rolling performance.

Inspection:

- Periodically check all screws, nuts, and bolts to ensure they are tight.
- Inspect the brake cables for any signs of wear or damage. Ensure they operate smoothly and effectively.
- Check the wheels for any cracks or excessive wear.

TROUBLESHOOTING

Rollator does not roll smoothly:

- Check wheels for debris (hair, lint, dirt) and clean thoroughly.
- Ensure all wheel axles are properly tightened.

Brakes are not engaging/disengaging properly:

- Inspect brake cables for kinks or damage.
- Ensure brake levers are not obstructed.
- If issues persist, consult a qualified technician.

SPECIFICATIONS

Feature	Detail
Model Number	RTL10266CF
Product Dimensions	27.75"D x 23.13"W x 38"H
Item Weight	14.5 Pounds
Material	Plastic, Carbon Fiber
Color	Black
Manufacturer	Drive Medical
Weight Capacity	Up to 300 lb

WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the official Drive Medical website or contact their customer service directly. Keep your purchase receipt for any warranty claims.

A detailed User Manual (PDF) is also available for download from the product's Amazon page for further reference.

IMPORTANT SAFETY WARNINGS

WARNING:

- This product is intended to provide mobility support. Do not use it for purposes other than its intended use.
- Always engage parking brakes before sitting on the seat or when the rollator is stationary.
- Ensure all parts are securely locked and adjusted before each use.
- Do not exceed the maximum weight capacity of 300 lbs.
- Exercise caution when navigating uneven surfaces or inclines.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B01JS6L8B6>

More documents and user manuals: <https://manuals.plus/>