



Manuals.plus /

- › Tommee Tippee /
- › Tommee Tippee Microwave Steam Sterilizer Instruction Manual, Model 52261030

## Tommee Tippee 52261030

# Tommee Tippee Microwave Steam Sterilizer Instruction Manual

Model: 52261030

## INTRODUCTION

---

Welcome to the Tommee Tippee Microwave Steam Sterilizer. This manual provides essential information for the safe and effective use of your sterilizer. Please read these instructions carefully before first use and retain them for future reference. This sterilizer is designed to eliminate 99.9% of harmful microbes and viruses from baby bottles and accessories using steam.

## IMPORTANT SAFETY INFORMATION

---

- Always keep this instruction manual for future reference.
- This appliance is intended for household use only. Do not use outdoors.
- Do not immerse the unit in water or other liquids.
- Ensure the sterilizer is placed on a stable, heat-resistant surface during use.
- Always add the recommended amount of water before operating. Operating without water can damage the unit.
- The contents will be extremely hot after the sterilization cycle. Allow a cooling period before handling. Use caution when opening the lid to avoid steam burns.
- Keep out of reach of children.
- Do not use abrasive cleaning agents or solvents.
- This product is designed for sterilizing baby bottles, pacifiers, breast pump parts, and other feeding accessories. Do not use for other purposes.
- Ensure your microwave is clean and free from food debris before use.
- Do not use metallic objects inside the sterilizer or microwave.

## PACKAGE CONTENTS

---

- 1 x Tommee Tippee Microwave Steam Sterilizer Unit
- 1 x Sterilizer Lid
- 1 x Accessory Tray

- 1 x Teat Tongs
- 1 x Instruction Manual (this document)



Image: The Tommy Tippee Microwave Steam Sterilizer, showing the main unit, clear lid, internal accessory tray, and a pair of white tongs, along with a pacifier.

## SETUP

1. Before first use, wash all components of the sterilizer and any items to be sterilized in warm soapy water, rinse thoroughly, and dry.
2. Place the sterilizer base on a flat, stable surface.
3. Add 200ml of cold tap water to the base of the sterilizer. Do not overfill or underfill.
4. Load bottles and accessories into the sterilizer. Ensure bottle openings face downwards. Place smaller items like teats and pacifiers in the accessory tray.

# Compatible

Suitable for most brands of bottles, teats, soothers and breast pump accessories.



Image: The sterilizer with bottles and various feeding accessories arranged inside, demonstrating how to load items for sterilization.

5. Secure the lid firmly onto the sterilizer base.

## OPERATING INSTRUCTIONS

1. After setup, place the loaded sterilizer into the microwave oven. Ensure the microwave turntable is clean and rotates freely.
2. Close the microwave door.
3. Set the microwave power and time according to the table below.

Microwave Wattage	Sterilization Time
800W - 1100W	4 minutes
500W - 750W	6 minutes

Note: Sterilization times are approximate and may vary based on microwave model and load. Always refer to your microwave's instruction manual for power settings.



Image: A graphic displaying "4 MINUTES" within a circular timer icon, representing the quick sterilization cycle.

4. Once the cycle is complete, allow the sterilizer to cool down for at least 3 minutes before removing it from the microwave. The unit and its contents will be extremely hot.
5. Carefully remove the sterilizer from the microwave using oven mitts or heat-resistant gloves.
6. Open the lid away from you to allow steam to escape safely.
7. Contents remain sterile for up to 24 hours if the lid is kept closed.

# Stays sterile

Contents stay sterile for 24 hours with the lid closed.



Image: A hand holding the closed sterilizer, illustrating that items inside remain sterile for an extended period when sealed.

8. Use the provided tongs to remove sterilized items.

## MAINTENANCE AND CLEANING

1. After each use, empty any remaining water from the sterilizer base.
2. Wash all components (base, lid, accessory tray, tongs) in warm soapy water. Rinse thoroughly.
3. For stubborn stains or limescale buildup, use a descaling solution suitable for kitchen appliances, following the manufacturer's instructions. Do not use harsh chemicals or abrasive pads.
4. Ensure all parts are completely dry before storing.
5. Store the sterilizer in a clean, dry place away from direct sunlight and heat.

## TROUBLESHOOTING

- **Sterilizer not working / No steam produced:** Ensure the correct amount of water (200ml) has been added to the base. Check that the microwave is functioning correctly and set to the appropriate power and time.
- **Items not sterile / Still wet:** Ensure items were thoroughly cleaned before sterilization. Verify that the lid

was securely closed during the cycle. Some moisture is normal after steam sterilization; allow items to air dry on a clean drying rack if desired.

- **Limescale buildup:** This is common in hard water areas. Descal the unit regularly using a suitable descaling solution. Follow the descaling product's instructions carefully.
- **Sterilizer too hot to handle:** Always allow a cooling period of at least 3 minutes after the microwave cycle. Use oven mitts or heat-resistant gloves when handling the hot sterilizer.

## SPECIFICATIONS

Feature	Detail
Model Number	52261030
Product Dimensions	10.79 x 10.79 x 5.91 inches
Item Weight	1.57 pounds
Material	BPA Free
Sterilization Cycle	4 minutes (with 1100W microwave)
Capacity	Up to 4 bottles

## WARRANTY INFORMATION

This Tommee Tippee product comes with a standard manufacturer's warranty. Please refer to the warranty card included with your purchase or visit the official Tommee Tippee website for detailed terms and conditions. Keep your proof of purchase for warranty claims.

## CUSTOMER SUPPORT

For further assistance, troubleshooting, or to inquire about replacement parts, please contact Tommee Tippee customer service.

- Website: [www.tommeetippee.com](http://www.tommeetippee.com)
- Contact information may vary by region. Please check the website for local support details.