

[Manuals.plus](#) /

› [Aroma Season](#) /

› Aroma Season Heated Eye Mask User Manual (Model B072QXQ9Y7)

Aroma Season B072QXQ9Y7

Aroma Season Heated Eye Mask User Manual

Model: B072QXQ9Y7

INTRODUCTION

Thank you for choosing the Aroma Season Heated Eye Mask. This manual provides essential information for the safe and effective use of your device. Please read it thoroughly before operation and retain it for future reference.

The Aroma Season Heated Eye Mask is designed to provide warm compression to the eyes, utilizing Graphene heating elements to generate far-infrared rays for quick and deep heating. This can assist in unclogging lacrimal glands and relieving discomfort associated with dry eyes, styes, chalazion, blepharitis, and MGD.



Figure 1: Aroma Season Heated Eye Mask (Violet)

Safety Information

- Do not use while sleeping.
- Stop heating immediately if you experience any discomfort.
- Exercise caution to prevent moderate temperature burns.

Product Overview

The Aroma Season Heated Eye Mask features a soft, breathable silk floss filling and a Graphene heating element. It includes a detachable controller for temperature and timer settings, and a natural lavender insert for enhanced relaxation.



Figure 2: Package Contents

SETUP

1. **Unpack:** Carefully remove all components from the packaging.
2. **Optional Lavender Insert:** If desired, place the natural dried lavender flower insert into the designated pocket within the eye mask.
3. **Connect Power:** Connect the eye mask's cable to the provided USB controller. Then, plug the USB end of the controller into a compatible USB power source (e.g., car charger, power bank, laptop, or a wall adapter).

Powered By USB

Plug in the eye mask for lasting soothing heat



Figure 3: Powering the Eye Mask via USB

OPERATING INSTRUCTIONS

1. **Power On:** Press and hold the power button on the controller to turn on the device. The display will illuminate.
2. **Adjust Temperature:** Use the temperature control button on the controller to select one of the 3 heating levels: 104°F (Low), 124°F (Medium), or 140°F (High).
3. **Set Timer:** Use the timer button on the controller to select a duration from 20 to 60 minutes (in 5-minute increments). A 20-minute warm therapy session is generally suggested for eye discomfort relief.
4. **Apply Mask:** Place the eye mask comfortably over your eyes, ensuring the adjustable strap provides a snug fit without being too tight.
5. **Relax:** Enjoy the soothing warmth. The mask will automatically shut off once the selected time has elapsed.
6. **Power Off:** To manually turn off the mask, press and hold the power button on the controller.

HELP

Soothing, effective relief
for Styes and Dry Eyes



✓ Dry Eyes



✓ Blepharitis



✓ MGD



✓ Styes




✓ Chalazion



✓ Ocular migraines

Figure 4: Controller for Temperature and Timer Settings

Upgraded Controller



5-Level Timing



Figure 5: Temperature and Timing Options

Video 1: Demonstration of the Heated Eye Mask for Dry Eye Relief (Seller Content)

MAINTENANCE

1. **Detachable Cover:** The outer cover of the eye mask is removable for cleaning.
2. **Washing:** Hand wash the cover with mild soap and water. Do not machine wash or tumble dry.
3. **Drying:** Allow the cover to air dry completely before reattaching it to the heating element.
4. **Storage:** Store the eye mask in a clean, dry place when not in use.

Detachable & Washable



Figure 6: Detachable and Washable Cover

TROUBLESHOOTING

If you encounter any issues with your Aroma Season Heated Eye Mask, please refer to the manufacturer's official support channels for troubleshooting assistance.

SPECIFICATIONS

Feature	Detail
Brand	Aroma Season
Model Number	B072QXQ9Y7
UPC	799355496728
Material	Satin (Cover), Silk Floss (Filling)
Color	Violet
Heating Element	Graphene

Temperature Settings	3 levels (104°F to 140°F)
Timer Settings	3 levels (20 to 60 minutes)
Special Features	Adjustable Strap, Lightweight, Detachable/Washable Cover, Lavender Insert
Power Source	USB
Package Dimensions	10.63 x 5.63 x 1.89 inches
Weight	6.35 ounces

WARRANTY & SUPPORT

For any product inquiries, warranty claims, or support needs, please contact Aroma Season directly. While specific warranty details are not provided in this manual, the product typically adheres to standard return policies. Please retain your proof of purchase for any claims.