

Respironics RE1102088

Respironics InnoSpire Replacement Filter User Manual

Model: RE1102088

1. PRODUCT OVERVIEW

The Respironics InnoSpire Disposable Replacement Filter (Model RE1102088) is designed for use with compatible Respironics InnoSpire nebulizer systems. This filter plays a crucial role in maintaining the cleanliness and efficiency of your nebulizer by preventing airborne particles from entering the compressor. Regular replacement ensures optimal performance and extends the life of your device.



Image 1.1: Respironics InnoSpire Disposable Replacement Filters. These small, white cylindrical filters are essential for maintaining the air quality and performance of your nebulizer system.

2. SAFETY INFORMATION

- Always ensure the nebulizer system is unplugged and turned off before attempting to replace the filter.
- Use only genuine Respironics replacement filters (Model RE1102088) to ensure compatibility and proper device function.
- Do not wash or reuse disposable filters. Attempting to clean or reuse a disposable filter can compromise its effectiveness and potentially damage the nebulizer.
- Keep filters out of reach of children and pets.
- Dispose of used filters according to local regulations for medical waste, if applicable, or general household waste.

3. PACKAGE CONTENTS

Your package contains:

- Respironics InnoSpire Disposable Replacement Filter(s) (Model RE1102088)

4. INSTALLATION AND REPLACEMENT INSTRUCTIONS

It is recommended to replace the air filter every 30 days or sooner if it appears discolored or dirty. Refer to your specific nebulizer system's user manual for exact filter housing location and removal instructions.

- 1. Prepare the Device:** Ensure your Respironics InnoSpire nebulizer system is turned off and unplugged from the power source.
- 2. Locate the Filter Housing:** The air filter housing is typically located on the side or back of the nebulizer compressor unit. It may have a small cap or cover.
- 3. Remove the Old Filter:** Carefully remove the filter housing cover. Gently pull out the old, used disposable filter. Note its orientation for correct installation of the new filter.
- 4. Insert the New Filter:** Take a new Respironics InnoSpire Disposable Replacement Filter (Model RE1102088) from its packaging. Insert the new filter into the filter housing, ensuring it is seated correctly and oriented as the old filter was. Do not force the filter into place.
- 5. Secure the Housing:** Replace the filter housing cover securely.
- 6. Resume Operation:** You may now plug in and turn on your nebulizer system.

5. CARE AND DISPOSAL

- Storage:** Store new, unused filters in their original packaging in a cool, dry place away from direct sunlight and extreme temperatures.
- Disposal:** Once a filter is removed, it should be immediately discarded. Do not attempt to clean, wash, or reuse the filter. Dispose of the used filter in regular household waste, unless local regulations require specific disposal methods for medical accessories.

6. TROUBLESHOOTING

Problem	Possible Cause	Solution
Filter does not fit into housing.	Incorrect filter model or improper orientation.	Ensure you are using Respironics InnoSpire Disposable Replacement Filter Model RE1102088. Verify the filter's orientation matches the original filter's placement. Do not force the filter.
Nebulizer performance is poor after filter replacement.	Filter housing not securely closed or filter not seated correctly.	Check that the filter is fully seated and the housing cover is securely closed. If issues persist, consult your nebulizer system's user manual or contact Respironics support.
New filter appears dirty quickly.	High levels of airborne particles in the environment.	This indicates the filter is working effectively. Replace the filter more frequently as needed. Ensure the nebulizer is used in a clean environment.

7. SPECIFICATIONS

Attribute	Detail
Model Number	RE1102088

Attribute	Detail
Brand	Respironics
Item Weight	0.634 ounces (approx. 18 grams)
Package Dimensions	3.35 x 1.42 x 0.43 inches (8.51 x 3.61 x 1.09 cm)
Type	Disposable Replacement Filter
Compatibility	Designed for compatible Respironics InnoSpire nebulizer systems. Refer to your device's manual.

8. WARRANTY AND SUPPORT

For specific warranty information regarding your Respironics InnoSpire Disposable Replacement Filter, please refer to the documentation provided with your original nebulizer system or contact Respironics customer support directly. As a disposable item, filters typically do not carry a separate warranty beyond initial manufacturing defects.

For technical assistance, product inquiries, or to purchase additional replacement filters, please visit the official Respironics website or contact their customer service department. Contact information can usually be found on the Respironics packaging or their official website.

Note: This manual provides general guidance. Always refer to your specific nebulizer system's user manual for detailed instructions and safety precautions.

© 2023 Respironics Inc. All rights reserved. Information subject to change without notice.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question