

[Manuals.plus](#) /

› [Starkey Products](#) /

› Starkey 6mm Comfort Closed/Occluded Ear Buds Instruction Manual

Starkey Products 6mm Comfort Closed/Occluded Ear Buds

Starkey 6mm Comfort Closed/Occluded Ear Buds Instruction Manual

1. PRODUCT OVERVIEW

These 6mm Starkey Comfort Closed/Occluded Ear Buds are designed for use with Starkey hearing aids. They provide a secure and comfortable fit, optimizing sound delivery and minimizing feedback. This pack contains 10 ear buds, ensuring you have replacements readily available.



This image displays multiple 6mm Starkey Comfort Closed/Occluded Ear Buds. These small, translucent dark gray silicone domes are designed to fit onto Starkey hearing aids, providing a comfortable and secure fit within the ear canal. Each ear bud features a closed design for optimal sound delivery and a small attachment point for the hearing aid receiver.

2. SETUP AND INSTALLATION

Installing the Ear Bud

1. Before inserting the ear bud, ensure your hearing aid is turned off.
2. To insert, gently pull up on the top of your ear to straighten the ear canal.
3. Carefully insert the ear bud until it is comfortable and sits snugly in the ear.

Removing the Ear Bud

1. To remove, first tilt your head toward the ear from which you are removing the ear bud.
2. Then, gently pull the ear bud out while still holding onto the hearing aid.

3. OPERATING INSTRUCTIONS

Once the ear bud is securely attached to your Starkey hearing aid and properly inserted into your ear, you can power on your hearing aid. The closed/occluded design of these ear buds helps to direct sound effectively into the ear canal, providing clear audio and reducing sound leakage. Ensure the ear bud remains comfortably seated for optimal performance throughout the day.

4. MAINTENANCE

Cleaning

- Clean the ear buds regularly with a soft, dry cloth to maintain hygiene and functionality.
- Avoid using harsh chemicals or abrasive materials, as these can damage the silicone.

Replacement

- Replace the ear buds when they become worn, damaged, or discolored.
- Always follow the manufacturer's instructions for replacement to ensure proper fit and function.

5. TROUBLESHOOTING

If you experience issues such as reduced sound quality, discomfort, or the ear bud not staying in place, consider the following:

- **Poor Sound Quality:** Ensure the ear bud is clean and free of debris. A clogged ear bud can significantly reduce sound clarity.
- **Discomfort:** Verify that the ear bud is correctly inserted and not causing pressure points. If discomfort persists, consult with your audiologist to ensure the 6mm size is appropriate for your ear canal.
- **Loose Fit:** If the ear bud feels loose or falls out, it may be worn out and require replacement. Ensure it is securely attached to the hearing aid receiver.
- **Feedback (Whistling):** This can occur if the ear bud is not properly sealed in the ear canal or if it is damaged. Re-insert the ear bud or replace it if necessary.

For persistent issues, contact your hearing aid specialist or the manufacturer for further assistance.

6. SPECIFICATIONS

Brand	Starkey Products
Size	6mm (10 count, Pack of 1)
Style	Uncorded, Closed/Occluded
Reusability	Disposable (intended for regular replacement)
Package Dimensions	6.18 x 4.69 x 0.43 inches
Weight	0.32 ounces
Manufacturer	Starkey

7. WARRANTY AND SUPPORT

Specific warranty information for these ear buds is not provided in this manual. For details regarding warranty coverage, product support, or to purchase additional ear buds, please refer to the official Starkey Products website or contact your authorized Starkey hearing aid provider.

You can visit the [Starkey Products Store](#) for more information.

8. IMPORTANT SAFETY INFORMATION

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is an accessory for hearing aids and should be used as directed by a healthcare professional or audiologist.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Email

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Details

Submit question