

Vicks VIH200

Vicks Sinus Inhaler User Manual

Model: VIH200 | Brand: Vicks

INTRODUCTION

The Vicks Sinus Steam Inhaler provides warm, non-medicated steam therapy designed for quick relief from sinus pain, congestion, cough, and cold symptoms. It targets the nasal passages, throat, and mouth to alleviate dryness and irritation, promoting easier breathing and more comfortable sleep. This portable device is compatible with Vicks VapoPads for an optional soothing menthol scent.

KEY FEATURES

- **Non-Medicated Steam Relief:** Provides temporary symptomatic relief from dryness and irritation caused by colds, seasonal allergies, and sinusitis.
- **Quick Treatment:** Delivers targeted steam therapy in just 5 to 15 minutes.
- **Easy to Use:** Operates with tap water and requires no filter.
- **VapoPads Compatibility:** Designed for use with Vicks VapoPads for enhanced comfort and soothing vapors (one sample pad included).
- **Adjustable Steam Control:** Allows users to customize steam intensity.
- **Auto Shut-off:** Ensures safety by automatically turning off when water runs out.
- **Easy to Clean:** Simple design for straightforward maintenance.

PRODUCT OVERVIEW



Figure 1: Vicks Sinus Inhaler and its retail packaging.



Figure 2: User demonstrating the targeted relief for sinus pressure, congestion, and seasonal discomfort.



Figure 3: The inhaler provides soothing comfort for congestion, cough, allergies, dryness, sinus pressure, and sore throat.



Figure 4: Fast-acting steam therapy penetrates nasal, sinus, and throat passages for quick soothing comfort in as little as 5 minutes.



Figure 5: Quick and easy to use features include a Comfort Flex Mask, VapoPad Tray, Works with Tap Water, No Filter Required, Two Adjustable Steam Levels, and Auto Shut-off Feature.

SETUP

1. **Unpack:** Carefully remove all components from the packaging.
2. **Place on Flat Surface:** Ensure the inhaler is placed on a stable, flat surface.
3. **Unlock Top Housing:** Twist and lift the top mask assembly to access the water reservoir.
4. **Fill with Water:** Fill the water reservoir with up to 2.5 tablespoons (40 ml) of clean tap water. Do not exceed the "MAX" fill line indicated inside the reservoir.
5. **Insert VapoPad (Optional):** If desired, slide a Vicks VapoPad into the designated tray located on the side of the mask assembly.
6. **Reattach Top Housing:** Securely place the mask assembly back onto the base and twist to lock it in place.
7. **Connect to Power:** Plug the power cord into a standard 110-volt electrical outlet.

OPERATING INSTRUCTIONS

1. **Power On:** Press the power button located on the front of the unit. The indicator light will illuminate.
2. **Wait for Steam:** Allow approximately 5-6 minutes for the unit to heat up and begin producing steam.
3. **Adjust Steam Intensity:** Use the adjustable steam control dial on the side of the mask assembly to set your preferred steam intensity.
4. **Inhale Steam:** Place your face comfortably into the soft mask, ensuring a good seal around your nose and mouth. Breathe normally through your nose and mouth, inhaling the warm steam.
5. **Treatment Duration:** Use for 5 to 15 minutes per session, or as needed for relief.
6. **Auto Shut-off:** The unit will automatically turn off when the water reservoir is empty.

MAINTENANCE

Regular cleaning of your Vicks Sinus Inhaler is essential to ensure optimal performance and hygiene.

- **Disconnect Power:** Always unplug the unit from the electrical outlet before cleaning.
- **Allow to Cool:** Let the unit cool down completely before handling, as components can be hot after use.
- **Disassemble:** Remove the mask assembly from the base. The mask can be detached for thorough cleaning.
- **Clean Components:** Wash the mask assembly and water reservoir with warm, soapy water. Rinse thoroughly to remove all soap residue.
- **Wipe Base:** Wipe the exterior of the base unit with a damp cloth. Do not immerse the base in water.
- **Dry Completely:** Ensure all parts are completely dry before reassembling or storing the unit.

TROUBLESHOOTING

Problem	Possible Cause	Solution
No steam produced.	Unit not plugged in; water reservoir empty; not enough time for heating.	Ensure unit is plugged in and power button is on. Fill water reservoir to MAX line. Allow 5-6 minutes for heating.
Steam output is too low.	Steam control dial set to low; mineral buildup.	Adjust steam control dial to a higher setting. Clean the unit to remove mineral buildup.
Unit does not turn on.	No power; faulty outlet.	Check power connection and try a different outlet.

SPECIFICATIONS

- **Model Number:** VIH200
- **Product Dimensions:** 5.16 x 5.16 x 8.58 inches
- **Weight:** 1.2 Pounds
- **Manufacturer:** Kaz, Incorporated
- **Country of Origin:** China
- **Voltage:** 110 volts

WARRANTY & SUPPORT

For detailed warranty information and customer support, please refer to the official Vicks website or the user manual provided with your product.

- Official User Manual (PDF): [Download Here](#)
- Visit the Vicks Store: [Vicks Store on Amazon](#)

