

[Manuals.plus](#) /

› [UpStart Components](#) /

› User Manual: 12 CPAP Filter Replacement for ResMed S9 VPAP ST

UpStart Components ResMed S9 VPAP ST CPAP Filter Replacement

User Manual: 12 CPAP Filter Replacement for ResMed S9 VPAP ST

Brand: UpStart Components

INTRODUCTION

This manual provides essential information for the proper use and maintenance of your UpStart Components CPAP Filter Replacements. These filters are designed to be compatible with ResMed S9 VPAP ST CPAP machines.

The filters are hypo-allergenic and play a crucial role in removing dust and pollen particles, ensuring you breathe clean air during therapy.

PRODUCT OVERVIEW



Image: A pack of twelve individual CPAP filters, showing their rectangular shape and white material.



Image: A close-up view of a single CPAP filter, highlighting its texture and thickness.



Image: Top-down view of a single CPAP filter, showing its flat surface.



Image: A slightly angled view of a single CPAP filter, emphasizing its thin profile.



Image: The product packaging for UpStart Components filters, displaying the brand logo and the twelve filters inside.

SETUP AND INSTALLATION

These replacement filters are designed for use with ResMed S9 VPAP ST CPAP machines. Refer to your CPAP machine's original user manual for specific instructions on accessing and replacing the air filter compartment.

1. Ensure your CPAP machine is turned off and unplugged from the power source before beginning installation.
2. Locate the air filter compartment on your ResMed S9 VPAP ST machine. This is typically found at the back or side of the device.
3. Open the filter compartment cover.
4. Carefully remove the old, dirty filter. Dispose of the old filter responsibly.

5. Insert a new UpStart Components replacement filter into the compartment, ensuring it fits snugly and correctly. The filter should lie flat without being bent or creased.
6. Close the filter compartment cover securely.
7. Plug your CPAP machine back into the power source.

It is recommended to wash your hands before and after handling filters to maintain hygiene.

OPERATION

Once installed, the UpStart Components CPAP filter will automatically function as part of your ResMed S9 VPAP ST machine's air filtration system. The filter works to remove airborne particles such as dust and pollen from the air before it reaches your respiratory system.

No special operation is required for the filter itself; its function is passive filtration as air passes through your CPAP device.

MAINTENANCE AND REPLACEMENT

Regular replacement of your CPAP filter is essential for maintaining optimal air quality and the efficiency of your CPAP machine.

- **Frequency:** These disposable filters should be replaced as soon as there is visible dirt or debris. For best filter efficiency, it is recommended to replace the filter every month or sooner, depending on usage and environmental conditions.
- **Inspection:** Periodically inspect the filter for signs of discoloration, dust accumulation, or damage.
- **Benefits of Replacement:** Regular replacement ensures continued filtration of airborne allergens and particles, contributing to healthier therapy.

Do not attempt to wash or reuse disposable filters, as this can compromise their filtration capabilities and potentially introduce contaminants.

TROUBLESHOOTING

While CPAP filters are simple components, issues can arise if they are not properly maintained or installed.

- **Reduced Airflow or Unusual Noise:**
 - Check if the filter is clogged with excessive dust or debris. Replace the filter if it appears dirty.
 - Ensure the filter is correctly seated in its compartment and the cover is securely closed. An improperly installed filter can restrict airflow or cause air leaks.
- **Visible Dirt or Particles in Air Path:**
 - This indicates the filter is no longer effectively trapping particles. Replace the filter immediately.
 - Ensure you are using the correct type and size of filter for your ResMed S9 VPAP ST machine.
- **Allergy Symptoms Worsening:**
 - If you experience increased allergy symptoms, your filter may be saturated or overdue for replacement. Replace the filter.

If problems persist after replacing the filter and ensuring proper installation, consult your CPAP machine's user manual or contact a qualified service technician.

SPECIFICATIONS

Attribute	Detail
Product Type	CPAP Filter Replacement
Compatibility	ResMed S9 VPAP ST CPAP Machines
Quantity	12 Filters (6 Pack - 2 Filters per Pack)
Material	Hypo-allergenic filtration material
Dimensions	2.19 x 1.5 x 0.13 inches
Brand	UpStart Components
Manufacturer	UpStart Components
ASIN	B075X1KSTM

WARRANTY AND SUPPORT

Please note that these are UpStart Components brand replacement parts and are not OEM products. This product is not affiliated with any OEM brands and is not covered under any warranties offered by the original manufacturers. Any warranties for this product are offered solely by **Upstart Battery**. For warranty claims or support inquiries regarding your UpStart Components CPAP filters, please contact Upstart Battery directly. For general product information or assistance, you may also refer to the [UpStart Components brand page on Amazon](#).

© 2024 UpStart Components. All rights reserved.
This manual is for informational purposes only. Always refer to your CPAP machine's original manufacturer manual for specific device instructions.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question