

[Manuals.plus](#) /

> [T POWER](#) /

> T POWER AC Adapter Charger for Spectra M1 Electric Breast Pump Instruction Manual

## T POWER AC Adapter Charger for Spectra M1

# T POWER AC Adapter Charger for Spectra M1 Electric Breast Pump Instruction Manual

## 1. INTRODUCTION

---

This manual provides essential instructions for the safe and effective use of the T POWER AC Adapter Charger. This adapter is specifically designed to provide power to the Spectra M1 Electric Breast Pump, including models M1 Plus REF IFU SP01-L1 UE09WCP-084090SPC.

Please read this manual thoroughly before using the product and retain it for future reference.



Image 1.1: T POWER AC Adapter Charger with its 6.6-foot cable, ready for connection.

## 2. SETUP INSTRUCTIONS

1. **Unpack the Adapter:** Carefully remove the T POWER AC Adapter Charger from its packaging.
2. **Inspect for Damage:** Before use, visually inspect the adapter and its cable for any signs of damage. Do not use if damaged.
3. **Connect to Power Outlet:** Plug the AC input prongs of the adapter into a standard wall power outlet (AC 100V - 240V).
4. **Connect to Device:** Insert the DC output connector of the adapter into the power input port of your Spectra M1 Electric Breast Pump. Ensure a secure connection.

The adapter features an extra-long 6.6 ft cable for convenient placement.

## 3. OPERATING INSTRUCTIONS

Once properly connected, the T POWER AC Adapter Charger will supply the necessary power to your Spectra M1 Electric Breast Pump.

- **Power On:** Operate your Spectra M1 Electric Breast Pump according to its manufacturer's instructions. The adapter provides continuous power.
- **Continuous Use:** The adapter is designed for reliable, continuous operation within its specified voltage range.



Image 3.1: Side view of the T POWER AC Adapter Charger, showing product details and safety warnings.

## 4. MAINTENANCE AND SAFETY

To ensure the longevity and safe operation of your T POWER AC Adapter Charger, follow these guidelines:

- **Cleaning:** Disconnect the adapter from the power outlet and the device before cleaning. Use a dry, soft cloth to wipe the surface. Do not use liquid cleaners.
- **Environment:** Avoid exposing the adapter to moisture, high humidity, extreme temperatures, or direct sunlight. Ensure adequate ventilation around the adapter during use.
- **Cable Care:** Do not bend, twist, or pull the cable excessively. Regularly inspect the cable for cuts, fraying, or other damage. Replace the adapter if the cable is damaged.
- **Disassembly:** Do not attempt to open, modify, or repair the adapter. This can void the warranty and pose safety risks.

- **Safety Features:** The adapter is manufactured with multiple smart features safeguarding against incorrect voltage (IV), short circuit (SC), and internal overheating (IO). It is CE/FCC/RoHS certified.



**power**

Over-Current Protection    Over-Temp Protection    Over-Voltage Protection    Short Circuit Protection

**SUPERIOR SAFETY POWER ADAPTER**



**Warm Tips**  
Please check our compatible Model List and Notice. Make sure the Model# or Part# match your device.  
**please feel free to contact us for any inquiries**



**Professional of Power Adapters**  
T-Power adapters are tested on the original device to ensure compatibility. Certified by T-power technicians 100% compatible.



**Safety Features with Safety certified**  
Safeguards against incorrect voltage, short circuit, and internal overheating. (CE/FCC/RoHS certified)



**Reinforced Power Cord**  
Long lasting DC cord for your convenience and charge your device.

**T-Power is a globally trusted provider of premier power supplies and power supply products.**



**30-Day Return for Any Reason**

You may return your undamaged product and packaging within 30 days of purchase to receive a FULL REFUND for any reason.



**1-Year Limited Manufacturer Warranty**

We will take care of all quality-related issues , **Easy Exchang** with a REPLACEMENT.\*



**Technical Support**

Received the wrong item?  
Looking for Charger for your device?  
Adapter tip doesn't fit your device?  
We can help, just let us know.

Image 4.1: Overview of the superior safety features integrated into the T POWER adapter design.

## 5. TROUBLESHOOTING

If you encounter issues with your T POWER AC Adapter Charger, refer to the following common problems and solutions:

- **No Power to Device:**
  - Ensure the adapter is securely plugged into both the wall outlet and the Spectra M1 Electric Breast Pump.
  - Verify that the wall outlet is functional by testing it with another device.
  - Check the adapter cable for any visible damage. If damaged, discontinue use.
- **Adapter Overheating:**
  - If the adapter feels excessively hot to the touch, immediately disconnect it from the power source and the device.
  - Ensure there is adequate air circulation around the adapter. Do not cover it during operation.
  - If overheating persists, contact customer support.
- **Device Not Charging/Operating Correctly:**

- Confirm that the adapter's specifications match the power requirements of your Spectra M1 Electric Breast Pump.
- Verify that your specific Spectra M1 model (M1 Plus REF IFU SP01-L1 UE09WCP-084090SPC) is listed as compatible.
- Refer to your breast pump's own instruction manual for device-specific troubleshooting.

## 6. SPECIFICATIONS

|                            |   |
|----------------------------|---|
| <b>Input Voltage Range</b> | AC 100V - 240V  |
| <b>Cable Length</b>        | 6.6 feet (approximately 2 meters)   |
| <b>Product Dimensions</b>  | 4 x 4 x 2 inches  |
| <b>Item Weight</b>         | 4 ounces  |
| <b>Manufacturer</b>        | T Power   |
| <b>Certifications</b>      | CE, FCC, RoHS   |
| <b>Compatibility</b>       | Spectra M1 Electric Breast Pump (M1 Plus REF IFU SP01-L1 UE09WCP-084090SPC) |

# Longer Lifespan

Tested to withstand  
10,000 bends




Image 6.1: Illustration of the adapter's durable construction and long lifespan features.

## 7. WARRANTY AND SUPPORT

T POWER stands behind the quality of its products. This AC Adapter Charger comes with:

- **30-Day Return Policy:** You may return undamaged products within 30 days of purchase for a full refund.
- **1-Year Limited Manufacturer Warranty:** This warranty covers quality-related issues and allows for easy exchange with a replacement.

For technical assistance, warranty claims, or any inquiries, please contact T POWER customer support. When contacting support, please have your product model number and purchase date available.

**Contact Information:** Please refer to the seller's information on your purchase platform or the official T POWER website for the most current contact details.

