

[Manuals.plus](#) /

› [LKPower](#) /

› LKPower AC Adapter for Walgreens BA-616W Blood Pressure Monitor User Manual

LKPower 180130454

LKPower AC Adapter for Walgreens BA-616W Blood Pressure Monitor User Manual

1. INTRODUCTION

This user manual provides essential information for the safe and effective operation of your LKPower AC Adapter, model 180130454. This adapter is specifically designed to provide power to the Walgreens BA-616W Healthy Living Blood Pressure Monitor. Please read this manual thoroughly before use and retain it for future reference.

2. PRODUCT FEATURES

- High-quality replacement AC wall power charger adapter.
- Universal input voltage: AC 100-240V 50/60Hz, suitable for worldwide use.
- Integrated over-voltage protection to safeguard your device.
- Integrated over-heat protection for safe operation.
- Tested and validated to ensure compatibility and reliable performance with the Walgreens BA-616W Blood Pressure Monitor.

3. SAFETY INFORMATION

To ensure safe operation and prevent damage to the adapter or your device, please observe the following safety precautions:

- Use the adapter only with compatible devices, specifically the Walgreens BA-616W Blood Pressure Monitor.
- Do not expose the adapter to water, moisture, or extreme temperatures.
- Avoid dropping or subjecting the adapter to severe impact.
- Do not attempt to open, modify, or repair the adapter. Refer servicing to qualified personnel.
- Ensure the power outlet is easily accessible and unplug the adapter when not in use for extended periods.
- Keep out of reach of children.

4. SETUP

1. **Inspect the Adapter:** Before first use, visually inspect the AC adapter and its cable for any signs of damage. Do not use if damaged.



Image: Front and side view of the LKPower AC Adapter. The adapter features a standard two-prong US plug on one side and an integrated cable with a barrel connector on the other. The adapter casing is black with ventilation grooves.

2. **Connect to Device:** Insert the barrel connector end of the AC adapter cable into the power input port of your Walgreens BA-616W Blood Pressure Monitor. Ensure a snug fit.



Image: The LKPower AC Adapter's cable is connected to a white electronic device, which appears to be a blood pressure monitor, resting on a white desk. The adapter itself is plugged into a power strip behind the device.

3. **Plug into Power Outlet:** Insert the two-prong plug of the AC adapter into a standard AC 100-240V wall power outlet.

5. OPERATING INSTRUCTIONS

Once the adapter is properly connected to both the blood pressure monitor and a power outlet, your device should receive power. The blood pressure monitor can now be operated according to its own user manual. The adapter converts the AC power from the wall outlet to the DC power required by your device.

6. MAINTENANCE

- **Cleaning:** Disconnect the adapter from the power outlet and the device before cleaning. Wipe the adapter with a soft, dry cloth. Do not use liquid cleaners or solvents.
- **Storage:** When not in use, store the adapter in a cool, dry place, away from direct sunlight and extreme temperatures. Avoid tightly coiling the cable to prevent damage.
- **Cable Care:** Do not bend, twist, or pull the cable excessively. Ensure the cable is not pinched or walked on.




Image: The LKPower AC Adapter is shown with its physical dimensions indicated. The main body of the adapter measures approximately 1.6 inches in width, 2.6 inches in length, and 1.2 inches in height. The attached power cable has a length of 76 inches.

7. TROUBLESHOOTING

If your device does not power on when connected to the LKPower AC Adapter, consider the following:

- **Check Connections:** Ensure the barrel connector is fully inserted into the device's power port and the AC plug is firmly seated in the wall outlet.
- **Test Power Outlet:** Plug another working appliance into the same power outlet to confirm the outlet is functional.
- **Inspect Adapter:** Look for any visible damage to the adapter casing or cable. If damaged, discontinue use.
- **Device Compatibility:** Confirm that your device is indeed the Walgreens BA-616W Blood Pressure Monitor, as this adapter is designed for specific compatibility.

If the issue persists after performing these checks, please contact customer support.

8. SPECIFICATIONS

Model Number	180130454
--------------	-----------

Brand	LKPower
Input Voltage	AC 100-240V, 50/60Hz
Compatible Devices	Walgreens BA-616W Healthy Living Blood Pressure Monitor
Connector Type	Barrel Connector
Special Feature	Travel-friendly design
Power Source	Corded Electric
Manufacturer	LKPower

9. WARRANTY AND SUPPORT

For warranty information or technical support, please refer to the documentation provided with your purchase or visit the official LKPower website. Please have your product model number (180130454) available when contacting support.

Manufacturer: LKPower

Date First Available: January 30, 2018

© 2024 LKPower. All rights reserved. This manual is for informational purposes only.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B079FQW4LT>

More documents and user manuals: <https://manuals.plus/>