

Pic Solution 02038204100000

Pic Solution AIRProjet Plus Ultrasonic Aerosol Nebulizer Instruction Manual

1. INTRODUCTION

The Pic Solution AIRProjet Plus is an advanced ultrasonic aerosol nebulizer designed for effective respiratory therapy. This device utilizes high-frequency sound waves to transform liquid medication into a fine mist for inhalation. Its compact design and quiet operation make it suitable for both adults and children, promoting respiratory well-being efficiently and comfortably.

This manual provides essential information for the safe and correct use of your AIRProjet Plus nebulizer. Please read all instructions carefully before first use and retain this manual for future reference.

2. SAFETY INFORMATION

Always observe the following safety precautions to prevent injury or damage to the device:

- Read all instructions thoroughly before operating the device.
- Do not immerse the main unit in water or any other liquid.
- Keep the device and its accessories out of reach of children.
- Use only medications prescribed by a healthcare professional for aerosol therapy.
- Ensure the device is placed on a stable, flat surface during operation.
- Unplug the device from the power outlet before cleaning or when not in use.
- Do not block the ventilation openings during operation.
- If the device is damaged or malfunctions, discontinue use immediately and contact customer support.

3. PACKAGE CONTENTS

Verify that all components are present in the package:

- AIRProjet Plus Ultrasonic Aerosol Nebulizer Main Unit
- 2-in-1 Soft Mask (suitable for adults and children)
- Connecting Tube
- Mouthpiece
- Power Adapter (not explicitly listed but implied for operation)
- Instruction Manual

4. PRODUCT OVERVIEW

Familiarize yourself with the components of your AIRProjet Plus nebulizer.



Figure 4.1: The AIRProjet Plus main unit shown alongside its retail packaging, highlighting its compact design.



Figure 4.2: Side view of the AIRProjet Plus nebulizer, showcasing its sleek and ergonomic form factor.



Figure 4.3: Top view of the AIRProjet Plus, illustrating the medication chamber where liquid medicine is placed for nebulization.



Figure 4.4: Essential accessories for the AIRProjet Plus, including the mask, connecting tube, and various connectors for personalized use.

5. SETUP

Follow these steps for initial setup before starting your therapy session:

1. Place the AIRProjet Plus main unit on a clean, stable, and flat surface.
2. Open the medication chamber cover (refer to Figure 4.3).
3. Carefully pour the prescribed liquid medication directly into the medication chamber. Do not overfill. The device is designed for direct medication loading, eliminating the need for additional water.
4. Close the medication chamber cover securely.
5. Connect one end of the connecting tube to the main inhalation opening on the device.
6. Attach the desired accessory (mask or mouthpiece) to the other end of the connecting tube. For children, use the '2-in-1' mask, ensuring it covers both the nose and mouth for optimal therapy efficiency. For adults, the mask can be used, or a mouthpiece if preferred.
7. Plug the power adapter into the device and then into a suitable electrical outlet.

6. OPERATING INSTRUCTIONS

Once the device is set up, proceed with the following steps to begin your aerosol therapy:

1. Press the power button to switch on the appliance. The device will begin nebulizing the medication silently and quickly.
2. Position the mask or mouthpiece comfortably. For masks, ensure a snug fit over the nose and mouth to prevent mist leakage.
3. Breathe normally and deeply through the accessory to inhale the medication mist.
4. Continue the treatment session until the medication flow becomes intermittent or stops, indicating the medication chamber is empty.
5. Press the power button again to turn off the device.
6. Unplug the device from the power outlet.

The AIRProjet Plus is compatible with most common drugs and saline solutions used for aerosol therapy.

7. CLEANING AND MAINTENANCE

Proper cleaning and maintenance are crucial for the longevity and hygienic operation of your nebulizer:

1. **After Each Use:** Disassemble all accessories (mask, mouthpiece, connecting tube). Wash them thoroughly with warm water and a mild detergent. Rinse them well under running water and allow them to air dry completely before storing. Do not sterilize accessories in boiling water or an autoclave.
2. **Main Unit:** Wipe the exterior of the main unit with a soft, damp cloth. Do not use abrasive cleaners or immerse the unit in water. Ensure the medication chamber is clean and dry after each use.
3. **Storage:** Store the clean and dry nebulizer and its accessories in a cool, dry place, away from direct sunlight and extreme temperatures.

8. TROUBLESHOOTING

If you encounter any issues with your AIRProjet Plus, refer to the table below for common problems and solutions:

Problem	Possible Cause	Solution
Device does not turn on.	Power adapter not connected properly or power outlet issue.	Ensure the power adapter is securely plugged into the device and a working electrical outlet.
No mist is produced.	No medication in the chamber, insufficient medication, or chamber cover not closed properly.	Add prescribed medication to the chamber. Ensure the chamber cover is securely closed.
Mist production is weak or intermittent.	Low medication level, dirty medication chamber, or clogged accessories.	Refill medication. Clean the medication chamber and all accessories thoroughly.
Unusual noise during operation.	Device not on a flat surface or internal issue.	Place the device on a stable, flat surface. If noise persists, contact customer support.

If the problem persists after attempting these solutions, please contact Pic Solution customer support for assistance.

9. SPECIFICATIONS

- **Model:** AIRProjet Plus

- **Model Number:** 02038204100000
- **Technology:** Ultrasonic Nebulization (High Frequency Sound Waves)
- **Product Dimensions:** 9.1 x 6.9 x 8.6 inches (23.1 x 17.5 x 21.8 cm)
- **Item Weight:** 2.86 pounds (1.3 kg)
- **Manufacturer:** PIC Solution

10. WARRANTY AND SERVICE LIFE

The Pic Solution AIRProjet Plus nebulizer is designed for durability and reliable performance.

- **Service Life of Appliance:** The estimated service life for the main AIRProjet Plus appliance is 3 years.
- **Service Life of Accessories:** The estimated service life for accessories (mask, tube, mouthpiece) is 1 year. Regular inspection and replacement of accessories are recommended to ensure optimal performance and hygiene.

For warranty claims or further support, please refer to the warranty information provided at the time of purchase or contact Pic Solution customer service.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Email

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Details

Submit question