

[Manuals.plus](#) /

- › [Soehnle Professional](#) /
- › SOEHNLE PROFESSIONAL Multina Comfort 8352 Baby Scale User Manual

Soehnle Professional 8352.01

SOEHNLE PROFESSIONAL Multina Comfort 8352 Baby Scale

User Manual

Model: 8352.01 | Brand: Soehnle Professional

1. INTRODUCTION

Thank you for choosing the SOEHNLE PROFESSIONAL Multina Comfort 8352 Baby Scale. This precision instrument is designed for accurate and comfortable weighing of infants and toddlers up to 20 kg. Its user-friendly design, including one-button operation and a removable baby tray, makes it an ideal tool for monitoring your child's growth. Please read this manual carefully before first use to ensure proper operation and to maintain the longevity of your scale.

2. SAFETY INSTRUCTIONS

- Always place the scale on a stable, flat, and hard surface to ensure accurate readings.
- Do not expose the scale to extreme temperatures, direct sunlight, or high humidity.
- Keep the scale away from strong electromagnetic fields.
- Do not immerse the scale in water or other liquids. Clean with a damp cloth only.
- Ensure the baby is placed securely and centrally on the weighing tray. Never leave a child unattended on the scale.
- Do not exceed the maximum weight capacity of 20 kg.
- If the scale is not used for an extended period, remove the batteries to prevent leakage.

3. PACKAGE CONTENTS

Please check the package contents upon unpacking:

- SOEHNLE PROFESSIONAL Multina Comfort 8352 Baby Scale Base Unit
- Removable Baby Weighing Tray
- 1 x 9V Battery
- User Manual (this document)

Note: An optional baby length measuring rod may be available separately.

4. SETUP

4.1 Battery Installation

1. Locate the battery compartment on the underside of the scale base unit.
2. Open the battery compartment cover.
3. Insert the 9V battery, ensuring correct polarity (+/-).
4. Close the battery compartment cover securely.

4.2 Attaching the Weighing Tray

The baby weighing tray is designed for easy attachment and removal.

1. Align the slots on the underside of the baby weighing tray with the corresponding pegs on the scale base unit.
2. Gently press down until the tray clicks securely into place.
3. To remove, lift the tray firmly upwards from the base unit.



Image: The SOEHNLE PROFESSIONAL Multina Comfort 8352 Baby Scale with its removable weighing tray securely attached, ready for use.

4.3 Optional Length Measuring Rod Attachment

If you have the optional length measuring rod, it can be attached to the scale for simultaneous weight and length measurements.



Image: The SOEHNLE PROFESSIONAL Multina Comfort 8352 Baby Scale with the optional length measuring rod installed, allowing for both weight and height measurements.

5. OPERATING INSTRUCTIONS

5.1 Powering On/Off

- To turn the scale ON, press the single operating button. The display will show "0.00" indicating it is ready for use.
- The scale will automatically turn OFF after a period of inactivity to conserve battery life. You can also press and hold the button to manually turn it off.

5.2 Weighing a Baby

1. Ensure the scale is on a stable surface and the display shows "0.00".
2. Carefully place the baby onto the center of the weighing tray.
3. Wait for the display to stabilize. The weight will be shown in kilograms (kg) or pounds (lb) depending on the selected unit.
4. Record the weight reading.
5. Carefully remove the baby from the scale.



Image: A child standing on the SOEHNLE PROFESSIONAL Multina Comfort 8352 scale base, demonstrating its versatility for weighing toddlers without the tray.

5.3 Tare Function (Zeroing)

The Tare function allows you to subtract the weight of a blanket or diaper before weighing the baby, ensuring only the baby's weight is measured.

1. Place the blanket or diaper on the weighing tray.
2. Press the operating button briefly. The display will reset to "0.00".
3. Now, place the baby on the blanket/diaper. The display will show only the baby's weight.

5.4 Unit Selection (kg/lb)

The scale supports both kilograms (kg) and pounds (lb). Refer to the specific instructions on the scale's display or a dedicated button (if present) to switch between units. Typically, this is done by a quick press of the power button when the scale is off, or a dedicated unit button.

6. MAINTENANCE AND CARE

- Clean the scale with a soft, damp cloth. Do not use abrasive cleaners or chemical solvents.
- Ensure the scale is completely dry before storage or next use.
- Store the scale in a cool, dry place, away from direct sunlight and extreme temperatures.
- Regularly check the battery compartment for any signs of corrosion.

7. TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Scale does not turn on.	Dead or incorrectly installed battery.	Check battery polarity. Replace the 9V battery.
Inaccurate readings.	Scale not on a flat, stable surface; weight limit exceeded; external interference.	Place scale on a hard, level surface. Ensure weight is within 20 kg. Move away from electronic devices.
Display shows "Lo".	Low battery power.	Replace the 9V battery.
Display shows "Err" or "O-Ld".	Overload (weight exceeds 20 kg).	Remove the excess weight immediately.

8. SPECIFICATIONS

Feature	Detail
Model	Multina Comfort 8352
Manufacturer Reference	8352.01
Brand	Soehnle Professional
Maximum Weight Capacity	20 kg (44 lbs)
Accuracy	10g
Display Type	LCD
Power Source	1 x 9V Battery (included)
Product Dimensions	60.1 x 34.6 x 9.2 cm (approx.)
Item Weight	3.5 kg (approx.)
Material	Plastic
Color	White
Measurement Type	Kilogram, Pound

9. WARRANTY AND SUPPORT

SOEHNLE PROFESSIONAL products are manufactured to high-quality standards and undergo strict quality controls. For information regarding warranty terms and conditions, please refer to the warranty card included with your product or visit the official Soehnle Professional website. For technical support, spare parts, or service inquiries, please contact Soehnle Professional customer service through their official channels.

Official Website: www.soehnle-professional.com

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B07JZW26Q2>

More documents and user manuals: <https://manuals.plus/>