

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

> [Medline](#) /

> [Medline MediGuard ES Nitrile Exam Gloves User Manual](#)

Medline MG3003

Medline MediGuard ES Nitrile Exam Gloves User Manual

Model: MG3003

INTRODUCTION

This user manual provides essential information for the safe and effective use of Medline MediGuard ES Nitrile Exam Gloves. These gloves are designed for medical examination purposes, offering tactile sensitivity and durable protection. Please read this manual thoroughly before use.

PRODUCT OVERVIEW

The Medline MediGuard ES Nitrile Exam Gloves are powder-free and offer a combination of tactile feel and durable nitrile protection. Key features include:

- **Enhanced Strength:** 77% stronger than standard stretch vinyl exam gloves.
- **Improved Grip:** Textured fingertips for enhanced grip and dexterity.
- **Easy Identification:** Color-coded sizing on each box helps quickly identify the correct size.
- **Convenient Packaging:** Available in various box counts (e.g., 300 gloves per box for Large size, 250 for XL) to maximize space and reduce storage needs.



Image description: A box of Medline MediGuard ES Nitrile Exam Gloves, displaying the product name, size (Large), and key features. The box is predominantly blue and white.

SETUP AND PREPARATION

Selecting the Correct Size

Proper fit is crucial for comfort, dexterity, and protection. Refer to the size guide on the product packaging to ensure you select the appropriate glove size. Gloves that are too small may tear easily, while gloves that are too large may compromise grip and dexterity.

Inspection Before Use

Before donning, visually inspect each glove for any tears, punctures, or defects. Do not use damaged gloves.

OPERATING AND USAGE INSTRUCTIONS

Donning Gloves

1. Ensure hands are clean and dry.
2. Remove one glove from the dispenser box.

3. Carefully insert one hand into the glove, ensuring fingers are fully seated.
4. Repeat for the other hand.
5. Adjust gloves for a comfortable and secure fit.

Doffing Gloves (Removal)

1. Grasp the outside of one glove near the wrist.
2. Peel the glove off, turning it inside out, and hold it in the gloved hand.
3. With the ungloved hand, slide fingers under the wrist of the remaining glove.
4. Peel the second glove off, turning it inside out, enclosing the first glove inside.
5. Dispose of the gloves immediately in an appropriate waste receptacle.

MAINTENANCE AND STORAGE

Store gloves in their original packaging in a cool, dry place, away from direct sunlight, fluorescent lighting, and excessive heat. Avoid storing near ozone-generating equipment. Ensure the box is sealed to prevent contamination.

TROUBLESHOOTING

Gloves Tear Easily

If gloves tear frequently, consider the following:

- **Incorrect Size:** The gloves may be too small. Try a larger size.
- **Improper Donning:** Ensure hands are dry and free of sharp objects (e.g., rings, long fingernails).
- **Material Stress:** Avoid excessive stretching or contact with sharp edges.

Poor Grip

The gloves feature textured fingertips for improved grip. If grip is still an issue, ensure hands are dry and the gloves are not contaminated with oils or lubricants.

SPECIFICATIONS

Feature	Detail
Model Number	MG3003
Material	Nitrile
Powder	Powder-Free
Size	Large
Package Dimensions	13.54 x 9.21 x 3.39 inches
Item Weight	2.25 pounds
Manufacturer	Medline Industries Healthcare
ASIN	B07N4PKBSW

Note: Specifications are subject to change without notice.

SAFETY INFORMATION

- For single use only. Do not reuse gloves.
- Not made with natural rubber latex.
- If irritation or allergic reaction occurs, discontinue use immediately and consult a medical professional.
- Keep out of reach of children.
- These gloves are intended for medical examination purposes and are not suitable for all chemical exposures or industrial applications. Consult specific chemical resistance charts if needed.

DISPOSAL

Dispose of used gloves according to local regulations for medical waste or general waste, depending on the nature of their use and any contaminants.

SUPPORT AND CONTACT INFORMATION

For further assistance or inquiries regarding Medline MediGuard ES Nitrile Exam Gloves, please contact Medline customer service.

Website: www.medline.com

Manufacturer: Medline Industries Healthcare