

Homedics CELL-501-EU

Homedics Replacement Filter for Anti-Cellulite Beauty Vacuum Smoother (Models CELL-500, CELL-600)

This manual provides instructions for the Homedics CELL-501-EU Replacement Filter.

INTRODUCTION

This document provides essential information for the proper use and maintenance of your Homedics Replacement Filter, model CELL-501-EU. This filter is designed for use with the Homedics Beauty Anti Cellulite Vacuum Smoother, specifically models CELL-500 and CELL-600, to ensure optimal performance and hygiene.

IMPORTANT SAFETY INFORMATION

Please read and understand all safety instructions before using the Homedics Beauty Anti Cellulite Vacuum Smoother with the replacement filter.

- Do not use the device or filter near water or in wet environments.
- Discontinue use immediately if the device causes any pain, discomfort, or adverse skin reactions.
- Consult a physician before use if you have any medical conditions, skin sensitivities, or are pregnant.
- Keep out of reach of children.

PRODUCT OVERVIEW



Figure 1: Homedics Beauty Anti Cellulite Vacuum Smoother. This image shows the main device with its controls and transparent suction cup attachment.



Figure 2: Homedics Replacement Filters (CELL-501-EU) shown alongside the main device. The image highlights the circular filter pads.

The Homedics Replacement Filter (CELL-501-EU) is an essential component for maintaining the hygiene and efficiency of your Homedics Beauty Anti Cellulite Vacuum Smoother. These filters prevent debris and impurities from entering the device's internal mechanisms, ensuring consistent suction and performance.

INSTALLATION AND SETUP

Follow these steps to install the replacement filter:

1. Ensure the Homedics Beauty Anti Cellulite Vacuum Smoother is turned off and disconnected from any power source.
2. Carefully remove the used filter from the device. Refer to your main device's instruction manual for specific filter compartment location.
3. Insert a new Homedics Replacement Filter (CELL-501-EU) into the designated filter slot. Ensure it is securely seated.
4. Reattach any covers or components that were removed during the filter replacement process.

OPERATING INSTRUCTIONS (WITH FILTER INSTALLED)

Once the filter is installed, you can operate your Homedics Beauty Anti Cellulite Vacuum Smoother as follows:

1. **Prepare the Device:** Attach the filter to the HoMedics Beauty Anti Cellulite Vacuum Smoother.
2. **Select Settings:** Choose the appropriate EMS (Electrical Muscle Stimulation) and vacuum massage settings based on the specific body area you intend to treat and your comfort level.
3. **Customize Treatment:** Utilize the two included cup sizes (with your main device) to personalize your treatment for different body contours.
4. **Apply to Skin:** Gently glide the device over the target skin area. It is recommended to use a suitable massage oil or cream to facilitate smooth movement and prevent skin irritation.
5. **Duration:** Follow the recommended treatment duration as specified in your main device's instruction manual.

MAINTENANCE AND FILTER REPLACEMENT

Regular filter replacement is crucial for the longevity and hygienic operation of your device.

- **Filter Lifespan:** Replace the filter after approximately 20 uses. This ensures optimal filtration and prevents clogging.
- **Cleaning:** The filters are disposable and not designed for cleaning or reuse.
- **Storage:** Store unused filters in a cool, dry place away from direct sunlight.

TROUBLESHOOTING

This section addresses common issues related to the filter. For issues with the main device, refer to its dedicated instruction manual.

- **Reduced Suction:** If you experience reduced suction, first check if the filter is correctly installed. If it is, the filter may be clogged and requires replacement.
- **Filter Not Fitting:** Ensure you are using the correct replacement filter (CELL-501-EU) for your Homedics Beauty Anti Cellulite Vacuum Smoother (CELL-500 or CELL-600).

SPECIFICATIONS

Product Name	Homedics Replacement Filter
--------------	-----------------------------

Model Number	CELL-501-EU
Compatible Devices	Homedics Beauty Anti Cellulite Vacuum Smoother (CELL-500, CELL-600)
Material	Charcoal (as per product specifications, likely for filtration)
Units	1 Total (referring to the pack quantity, typically)
Dimensions (Package)	7.2 x 6.8 x 2.6 cm
Weight (Package)	300 grams
Manufacturer	FKA Brands
Country of Origin	China

WARRANTY AND SUPPORT

For warranty information and customer support regarding your Homedics products, please refer to the documentation included with your original Homedics Beauty Anti Cellulite Vacuum Smoother device or visit the official Homedics website. Keep your proof of purchase for any warranty claims.