

Omron HEM-6232T

OMRON HEM-6232T Automatic Wrist Blood Pressure Monitor

User Instruction Manual

[Information](#)

[Setup](#)

[Operation](#)

[Product Overview](#)

[Safety](#)

[Maintenance](#)

[Troubleshooting](#)

[Specifications](#)

[Warranty &](#)

[Support](#)

1. PRODUCT OVERVIEW

The OMRON HEM-6232T Automatic Wrist Blood Pressure Monitor is a compact, portable, and accurate device designed for convenient blood pressure measurement. It incorporates OMRON's renowned technology to provide reliable readings and features for comprehensive health monitoring.

Key features include:

- **Memory Storage:** Stores up to 100 measurements for 2 different users, allowing for long-term tracking.
- **Morning Hypertension Indicator:** Alerts users to potential morning hypertension.
- **Position Sensor:** Guides the user to the correct wrist position for accurate readings.
- **Movement Error Detector:** Identifies and alerts if body movement affects measurement accuracy.
- **Irregular Heartbeat Detector:** Indicates the presence of an irregular heartbeat during measurement.
- **Bluetooth Connectivity:** Seamlessly connects with the OMRON connect app for data synchronization and analysis.
- **Average Value Function:** Calculates the average of previous readings.
- **Date and Time Display:** Records the date and time of each measurement.
- **Compact Design:** Easy to transport and use at home, work, or while traveling.



Figure 1.1: Front view of the OMRON HEM-6232T Wrist Blood Pressure Monitor, displaying a typical blood pressure reading and pulse rate.



Figure 1.2: The OMRON HEM-6232T monitor alongside its retail packaging, highlighting its compact size and design.



Figure 1.3: The monitor display indicating compatibility with the OMRON connect mobile application, enabling data synchronization.



Figure 1.4: A detailed diagram labeling the key components of the OMRON HEM-6232T, including the cuff, display, and control buttons.



Figure 1.5: An angled perspective of the OMRON HEM-6232T, showcasing its sleek black design and wrist cuff.

2. SAFETY INFORMATION

Before using the OMRON HEM-6232T, please read all instructions thoroughly. This device is intended for general health monitoring and should not be used for self-diagnosis or treatment of medical conditions. Always consult with a healthcare professional for interpretation of blood pressure readings and before making any medical decisions.

- **Consult your doctor:** Always consult your physician regarding your blood pressure readings and any health concerns.
- **Read instructions for use:** Ensure you fully understand the operation and safety guidelines provided in this manual.
- Do not adjust medication based on measurements from this device without consulting your doctor.
- Keep the device out of reach of children and pets.
- Use only specified batteries and dispose of them properly.

3. SETUP

3.1. Battery Installation

The OMRON HEM-6232T operates on 2 AAA alkaline batteries. Ensure the batteries are inserted with the correct polarity as indicated inside the battery compartment.

1. Open the battery cover on the back of the monitor.
2. Insert two new AAA alkaline batteries, matching the (+) and (-) terminals.
3. Close the battery cover securely.

3.2. Initial Settings (Date, Time, User)

Upon first use or after battery replacement, you will need to set the date and time. The device supports two user profiles (User 1 and User 2).

1. Press the **START/STOP** button to turn on the device.
2. Follow the on-screen prompts to set the year, month, day, hour, and minute using the adjustment buttons (refer to

your full instruction manual for specific button functions).

3. Select the desired user profile (User 1 or User 2) before taking a measurement. This ensures readings are stored correctly for each individual.

3.3. Pairing with OMRON connect App

The OMRON HEM-6232T can connect to the OMRON connect app via Bluetooth for easy data management.

1. Download the "OMRON connect" app from your smartphone's app store (available for iOS and Android).
2. Enable Bluetooth on your smartphone.
3. Open the OMRON connect app and follow the on-screen instructions to add a new device.
4. On the blood pressure monitor, activate Bluetooth pairing mode (refer to your full instruction manual for specific steps, usually involves holding a button).
5. Confirm the pairing on both the app and the monitor.

4. OPERATION

4.1. Proper Cuff Placement

Accurate readings depend on correct cuff placement. The cuff should be placed directly on your wrist, approximately 1-2 cm above your hand.

1. Remove any tight clothing or accessories from your wrist.
2. Wrap the cuff firmly around your left wrist (or right, if preferred, but be consistent), ensuring the monitor's display is on the inside of your wrist.
3. The bottom edge of the cuff should be about 1 to 2 cm from the bottom of your palm.
4. Secure the cuff snugly, but not too tightly.

4.2. Correct Body Position for Measurement

For accurate results, maintain a relaxed and still posture during measurement.

- Sit comfortably in a chair with your back supported and feet flat on the floor.
- Place your elbow on a table so that the cuff is at the same level as your heart. The monitor has a position sensor (blue light) to help guide you.
- Remain still and silent during the measurement. Avoid talking or moving.
- Take measurements at the same time each day, ideally in the morning and evening.

4.3. Taking a Measurement

1. Ensure the cuff is correctly placed and you are in the proper position.
2. Press the **START/STOP** button. The cuff will begin to inflate.
3. Remain still until the measurement is complete and the readings are displayed.
4. The monitor will display your Systolic Pressure (SYS), Diastolic Pressure (DIA), and Pulse Rate.
5. The readings will automatically be saved to the selected user profile.

4.4. Understanding Your Readings

The display shows three main values:

- **SYS (Systolic Pressure):** The pressure when your heart beats.
- **DIA (Diastolic Pressure):** The pressure when your heart rests between beats.

- **PULSE:** Your heart rate per minute.

The monitor also features an **Hypertension Indicator** which will illuminate if your systolic or diastolic pressure is above recommended levels. An **Irregular Heartbeat Symbol** will appear if an irregular heartbeat is detected during the measurement.

5. MAINTENANCE

5.1. Cleaning

To ensure the longevity and hygiene of your blood pressure monitor:

- Wipe the monitor and cuff with a soft, dry cloth.
- If necessary, use a cloth lightly dampened with water or a mild detergent, then wipe dry.
- Do not use abrasive cleaners, solvents, or immerse the device in water.

5.2. Storage

Store the device in its provided carrying case when not in use. Keep it in a clean, dry place away from direct sunlight, extreme temperatures, and humidity.

5.3. Battery Replacement

When the low battery indicator appears on the display, replace both AAA batteries simultaneously. Do not mix old and new batteries or different types of batteries.

6. TROUBLESHOOTING

If you encounter issues with your OMRON HEM-6232T, refer to the common problems and solutions below. If the problem persists, contact customer support.

Problem	Possible Cause	Solution
Inaccurate Readings	Incorrect cuff placement, movement during measurement, improper body position (wrist not at heart level).	Ensure cuff is snug and correctly positioned. Remain still and silent. Use the position sensor to keep the monitor at heart level.
Display is dim or blank	Low batteries, poor lighting conditions (as noted in reviews).	Replace batteries. Ensure adequate lighting to view the display. Adjust viewing angle if necessary.
Bluetooth connection issues	Bluetooth off on phone, app not open, device not in pairing mode.	Ensure Bluetooth is enabled on your smartphone. Open the OMRON connect app. Re-initiate pairing mode on the monitor.
Error message on display	Cuff not inflated properly, excessive movement, irregular heartbeat detected.	Re-apply cuff and try again. Remain still. If irregular heartbeat symbol persists, consult a doctor.

7. SPECIFICATIONS

Product Name	Wrist Blood Pressure Monitor
Model Name	HEM-6232T-LA
Brand	Omron
Power Source	Batteries (2 x AAA alkaline)
Usage	Wrist
Display Type	Digital
Weight	217 grams
Band Size	19.5 cm
Product Dimensions (L x W x H)	11 cm x 9 cm x 9 cm
Included Components	Monitor, Wrist Cuff, 2 AAA Alkaline Batteries, Carrying Case, OMRON connect Guide, Instruction Manual
Memory	100 readings for 2 users
Connectivity	Bluetooth

8. WARRANTY AND SUPPORT

8.1. Warranty Information

The OMRON HEM-6232T Automatic Wrist Blood Pressure Monitor comes with a **5-year warranty**. This warranty covers defects in materials and workmanship under normal use. Please retain your proof of purchase for warranty claims.

8.2. Customer Support

For technical assistance, warranty claims, or further inquiries, please contact OMRON customer support. Refer to the contact information provided in the packaging or on the official OMRON website for your region.

For more information, visit the official OMRON website: www.omronhealthcare.com

© 2024 OMRON Healthcare. All rights reserved.

This manual is for informational purposes only. Consult a healthcare professional for medical advice.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question