

foreverLily 1

foreverLily 7-Color LED Facial Mask Phototherapy Device User Manual

1. INTRODUCTION

Thank you for choosing the foreverLily 7-Color LED Facial Mask Phototherapy Device. This manual provides essential information for the safe and effective use of your device. Please read it thoroughly before first use and keep it for future reference. This device is designed to support skin rejuvenation, reduce acne, and improve overall skin health through LED light therapy.

2. SAFETY INFORMATION

- Read all instructions carefully before using the device.
- Keep the device out of reach of children.
- Do not use the device on broken, irritated, or infected skin.
- Consult a healthcare professional before use if you are pregnant, breastfeeding, have epilepsy, thyroid conditions, light sensitivity, or are taking photosensitizing medications.
- Avoid direct eye contact with the LED lights. Protective eyewear is recommended.
- Do not immerse the device in water or use it in wet environments.
- Use only the provided power adapter.
- Discontinue use if you experience any discomfort or adverse reactions.

3. PRODUCT OVERVIEW

The foreverLily 7-Color LED Facial Mask is a phototherapy device designed for home use. It features 256 LED lamps and offers 7 different light colors, each with specific benefits for various skin concerns. The mask is made from environmental ABS material, ensuring it is non-toxic, odorless, and durable. It conforms to the facial contour for comfortable wear.



Image 3.1: The foreverLily 7-Color LED Facial Mask with its remote control, power adapter, and cables.

BEAUTY TREATMENT



Image 3.2: The LED Facial Mask shown with its remote control, indicating the various light colors available.

4. PACKAGE CONTENTS

Please check the package to ensure all items are present:

- 1 x LED Facial Mask
- 1 x Remote Control
- 1 x Power Cord
- 1 x Head Retention Strap
- 1 x Power Adapter
- 1 x User Manual (English)

5. SETUP

1. Unpack all components from the box.
2. Connect the power cord to the LED facial mask.
3. Connect the power cord to the power adapter.
4. Plug the power adapter into a suitable electrical outlet.
5. Attach the head retention strap to the mask for a secure and comfortable fit during use.

6. OPERATING INSTRUCTIONS

Follow these steps to operate your LED facial mask:

1. **Preparation:** Cleanse your face thoroughly and ensure it is dry. You may apply a high-quality moisturizing skincare product before use for enhanced effects.
2. **Wear the Mask:** Place the mask comfortably on your face and secure it with the head retention strap. Ensure your eyes are protected from direct light.
3. **Power On:** Press the power button on the remote control to turn on the device.
4. **Select Light Color:** Use the remote control to cycle through the 7 available light colors (Red, Blue, Green, Yellow, Purple, Cyan, White) and select the desired color for your treatment.
5. **Adjust Intensity:** The device offers 10 speed/intensity levels. Use the intensity buttons on the remote to adjust the light intensity to a comfortable level.
6. **Set Timer:** Select a treatment duration using the timer function on the remote control. Available options are 5, 10, 15, 20, 25, or 30 minutes.
7. **Treatment:** Relax and allow the mask to perform the light therapy for the selected duration. The device will automatically turn off once the timer expires.
8. **Power Off:** After the session, you can manually turn off the device by pressing the power button on the remote control if it hasn't already turned off automatically.



Image 6.1: A user wearing the foreverLily LED Facial Mask during a treatment session.

7. LED LIGHT FUNCTIONS AND BENEFITS

Each of the 7 LED light colors offers unique benefits for skin treatment:

- **Red Light (630 nm):** Increases cell activity, accelerates metabolism, and promotes collagen formation. Helps smooth skin and reduce fine lines.
- **Green Light (520 nm):** Neutralizes, balances, and stabilizes skin. Can reduce skin oil secretion, smooth wrinkles, and fine lines.
- **Blue Light (470 nm):** Bactericidal and anti-inflammatory. Inhibits inflammation, repairs skin, and improves sensitive skin.
- **Yellow Light (590 nm):** Improves cellular oxygen exchange, replenishes energy to skin cells. Gentle treatment for acne.
- **Purple Light:** Dual-frequency light. Repairs acne, relaxes skin, and improves lymphatic metabolism.

- **Cyan Light:** Enhances cell energy and promotes metabolism. Prevents allergies.
- **White Light:** Improves fine lines and loose skin. Enhances lymphatic metabolism.



RED(630nm)

Increase cell activity, speed up metabolism, and promote collagen formation

(520nm)GREEN

Neutralization, balance and stability, can reduce skin oil secretion



BLUE(470nm)

Bactericidal anti-inflammatory, inhibit inflammation, repair skin

(590nm)YELLOW

Improve cell oxygen exchange, replenish energy to skin cells



PURPLE

Dual frequency light, repairing acne.

CYAN

Enhances cell energy and promotes metabolism



WHITE

Improve fine lines and loose skin.

Image 7.1: Visual representation of each LED light color and its corresponding wavelength and skin benefits.



Image 7.2: Overview of the comprehensive skin benefits achievable with the LED facial mask, including tightening, whitening, and anti-aging effects.

TIGHTING

WHITENING

REJUVENATION

BLEMISH

ANTI - AGING

ACNE



Tips: The effect will better if used with a high quality moisturizing skin care product

Image 7.3: Examples of skin improvements for various conditions such as sun spots, blemishes, and wrinkles, illustrating the device's effectiveness.

8. MAINTENANCE

- **Cleaning:** After each use, gently wipe the mask with a soft, dry cloth. Do not use abrasive cleaners or immerse the mask in water.
- **Storage:** Store the device in a cool, dry place, away from direct sunlight and extreme temperatures. Keep it in its original packaging or a protective case to prevent damage.
- **Power Cord:** Regularly inspect the power cord and adapter for any signs of damage. If damaged, discontinue use and contact customer support.

9. TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on	Power cord not properly connected; Power adapter not plugged in; No power from outlet.	Ensure all connections are secure. Plug the adapter into a working outlet.
LED lights are dim or not working	Low power supply; Device malfunction.	Check power connection. If problem persists, contact customer support.
Remote control not responding	Remote control batteries are low or depleted.	Replace batteries in the remote control.
Discomfort during use	Light intensity too high; Skin sensitivity.	Reduce light intensity. Discontinue use if discomfort persists and consult a doctor.

10. SPECIFICATIONS

- **Model Number:** 1
- **Material:** ABS
- **Input:** 100-240 V, 50/60 Hz
- **Output:** 12 V, 1 A
- **LED Quantity:** 256
- **Weight:** 260 g
- **Dimensions (approx.):** 19 cm x 14 cm (Mask)

11. WARRANTY AND SUPPORT

The foreverLily 7-Color LED Facial Mask Phototherapy Device comes with a standard manufacturer's warranty. For specific warranty details, please refer to the documentation included with your purchase or contact your retailer.

If you encounter any issues or have questions regarding your device, please contact foreverLily customer support through the retailer's platform or the contact information provided in your purchase documentation. Our support team is available to assist you with troubleshooting, parts replacement, and general inquiries.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question