

[Manuals.plus](#) /

> [Beurer](#) /

> Beurer BF720 Body Analysis Scale User Manual

Beurer BF720

Beurer BF720 Body Analysis Scale User Manual

Model: BF720 (Item No. 74938)

INTRODUCTION

Welcome to your Beurer BF720 Body Analysis Scale. This manual provides essential information for the safe and effective use of your new scale. The BF720 is designed to provide a comprehensive analysis of your body values, including body weight, body fat, body water percentage, muscle percentage, bone mass, BMI, and calorie consumption (AMR/BMR). It features Bluetooth connectivity for seamless data transfer to the Beurer HealthManager Pro app, eight user memories with automatic recognition, and five activity levels.

IMPORTANT SAFETY INFORMATION

Please read these instructions carefully before use and keep them for future reference.

- This scale is intended for personal use only and not for medical or commercial purposes.
- Do not step on the scale with wet feet or if the surface is wet, as this may cause slipping.
- Keep the scale away from strong electromagnetic fields.
- Do not attempt to disassemble or repair the scale yourself. Contact authorized service personnel.
- Batteries should be handled carefully. Ensure correct polarity. Remove batteries if the scale will not be used for an extended period.
- This device is not suitable for persons with medical implants (e.g., pacemakers) as the electrical current may affect their function.
- Pregnant women should not use this scale for body analysis measurements, as the results may be inaccurate.

PRODUCT OVERVIEW

The Beurer BF720 Body Analysis Scale features a sleek design with a clear digital display and integrated electrodes for body analysis.



Figure 1: Top view of the Beurer BF720 Body Analysis Scale, showing the digital display and four silver electrodes.



Figure 2: A person interacting with the Beurer BF720 Body Analysis Scale, demonstrating its use in a home environment.

SETUP

Unpacking and Initial Inspection

Carefully remove the scale from its packaging. Inspect the device for any signs of damage. If damaged, do not use and contact customer support.

Battery Installation

The Beurer BF720 requires 3 AAA batteries (included).

1. Open the battery compartment on the underside of the scale.
2. Insert the 3 AAA batteries, ensuring correct polarity (+/-).
3. Close the battery compartment securely.

Placement

Place the scale on a firm, flat surface. Avoid carpets or uneven flooring, as this can affect measurement accuracy.

App Download and Pairing (Beurer HealthManager Pro)

To utilize the full features of your BF720 scale, download the free Beurer HealthManager Pro app.

1. Search for "Beurer HealthManager Pro" in your smartphone's app store (iOS or Android).
2. Install the app.
3. Open the app and follow the on-screen instructions to create a user profile.
4. Ensure Bluetooth is enabled on your smartphone.
5. Step onto the scale briefly to activate it. The scale will display "BT" or a Bluetooth symbol when ready to pair.
6. Within the app, navigate to "Devices" or "Add Device" and select the BF720 scale to pair.
7. Follow any further prompts in the app to complete the pairing process and set up your personal data (age, height, gender, activity level).

OPERATING INSTRUCTIONS

Taking a Weight Measurement

To take a simple weight measurement:

1. Place the scale on a hard, flat surface.
2. Step onto the scale with bare feet. Ensure your feet are positioned on the electrodes.
3. Stand still until your weight is displayed and stabilizes.
4. The scale will automatically turn off after a few seconds.

Taking a Body Analysis Measurement

For a comprehensive body analysis (weight, fat, water, muscle, bone mass, BMI, BMR/AMR):

1. Ensure your personal profile is set up in the Beurer HealthManager Pro app and paired with the scale.
2. Step onto the scale with bare feet, ensuring good contact with the electrodes.
3. Stand still. The scale will first display your weight, then cycle through your body analysis values.
4. Your measurements will automatically transfer to the Beurer HealthManager Pro app if Bluetooth is active and the app is open or running in the background.

User Recognition

The BF720 features automatic user recognition for up to eight users. Ensure your profile is correctly set up in the app. The scale will attempt to identify you based on your weight and other stored data. If multiple users have similar weights, the scale may prompt you to select your user number.

Activity Levels

The scale allows for five different activity levels to ensure accurate calorie requirement calculations. Adjust your activity level in the Beurer HealthManager Pro app based on your physical activity.

- **Level 1:** No physical activity.
- **Level 2:** Light physical activity (e.g., walking, light housework).
- **Level 3:** Moderate physical activity (e.g., 1-2 times per week sports).
- **Level 4:** Active physical activity (e.g., 3-5 times per week sports).
- **Level 5:** Very active physical activity (e.g., daily intensive sports).

MAINTENANCE AND CARE

Cleaning

Clean the scale regularly with a damp cloth. Do not use abrasive cleaning agents or immerse the scale in water.

Storage

Store the scale in a dry place, away from extreme temperatures and direct sunlight.

Battery Replacement

When the battery indicator appears on the display, replace all 3 AAA batteries. Dispose of old batteries according to local regulations.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Scale does not turn on.	Batteries are dead or incorrectly inserted.	Check battery polarity or replace batteries.
Inaccurate measurements.	Scale on uneven surface; not standing still; dirty electrodes; clothing worn.	Place on firm, flat surface; stand still; clean electrodes; measure with bare feet and minimal clothing.
Bluetooth connection issues.	Bluetooth off on phone; app not open; scale not activated; interference.	Ensure Bluetooth is on; open app; step on scale to activate; move closer to phone; restart phone/scale.
"Err" message on display.	Overload; measurement error.	Ensure weight is within limit (720 lbs); re-attempt measurement.

TECHNICAL SPECIFICATIONS

- **Model:** BF720 (Item No. 74938)
- **Dimensions:** 11.81 x 11.81 x 0.98 inches (30 x 30 x 2.5 cm)
- **Weight:** Approximately 3.53 Pounds (1.6 Kilograms)
- **Weight Limit:** 720 Pounds (326.5 Kilograms)
- **Material:** Alloy Steel
- **Power Supply:** 3 x 1.5V AAA batteries (included)
- **Measurement Units:** kg / lb / st
- **Connectivity:** Bluetooth
- **App Compatibility:** Beurer HealthManager Pro (iOS/Android)
- **User Memories:** 8
- **Activity Levels:** 5
- **Functions:** Weight, Body Fat, Body Water, Muscle Percentage, Bone Mass, BMI, AMR/BMR

WARRANTY AND SUPPORT

For warranty information, technical support, or service inquiries, please refer to the official Beurer website or contact Beurer customer service. Keep your proof of purchase.

Manufacturer: Beurer

Website: www.beurer.com

Date First Available: August 15, 2023