

[Manuals.plus](#) /

› [XSOUL](#) /

› XSOUL At-Home IPL Hair Removal Device D-1103 User Manual

XSOUL D-1103

XSOUL At-Home IPL Hair Removal Device

MODEL: D-1103

Thank you for choosing the XSOUL IPL Hair Removal Device. Please read this manual carefully before use and keep it for future reference.

1. Introduction and Overview

The XSOUL IPL Hair Removal Device utilizes advanced Intense Pulsed Light (IPL) technology to target hair follicles and disrupt the hair regrowth cycle. Designed for safe and effective at-home use, this device offers a convenient solution for permanent hair reduction on various body areas.

- **Advanced IPL Technology:** Breaks the cycle of hair regrowth for lasting results.
- **High Flash Capacity:** Equipped with 999,999 flashes, sufficient for lifelong use.
- **Adjustable Energy Levels:** 5 energy levels to customize treatment intensity based on skin sensitivity.
- **Dual Operation Modes:** Features Auto Mode for continuous flashing on large areas and Manual Mode for precise treatment on smaller, sensitive areas.

XSOUL

IPL HAIR REMOVAL

LEAVES YOU WITH SILKY SMOOTH SKIN


efficient
hair removal


safe &
painless


eyes
protection


5 energy
levels


years of
usage

AUTO

manual & auto
modes



Figure 1: XSOUL IPL Hair Removal Device highlighting its main features.

2. Safety Information

Your safety is paramount. Please adhere to all safety guidelines provided in this manual.

- **Eye Protection:** Always wear the provided safety goggles during operation to protect your eyes from the intense pulsed light.

- **Skin Contact:** The device will only activate when in full contact with the skin, ensuring safe and targeted treatment.
- **Skin Sensitivity:** It is recommended to start treatment at the lowest energy level (Level 1) and gradually increase it based on your skin's adaptability.
- **Contraindications:** The device is **not suitable** for:
 - Areas with tattoos
 - Sunburned skin
 - Men's facial hair areas (beard, mustache)
 - Pregnant women
 - Individuals with very light hair colors (white/grey, light blonde, red) or very dark skin tones. Refer to the compatibility chart below.



Figure 2: Skin Tone & Body Hair Type Chart. Green check indicates applicability, red cross indicates not applicable.




Figure 3: Important areas and conditions where the XSOUL IPL device should not be used.

3. Product Components

Your XSOUL IPL Hair Removal Device package includes the following items:

- XSOUL IPL Hair Removal Device
- Power Adapter
- Safety Goggles
- Disposable Razors (2 units)
- User Manual



Figure 4: All components included with the XSOUL IPL Hair Removal Device.

4. Setup

1. **Prepare the Treatment Area:** Shave any visible hair from the area you intend to treat. Ensure the skin is clean, dry, and free from any lotions, creams, or deodorants.
2. **Wear Safety Goggles:** Put on the provided safety goggles to protect your eyes from the intense light flashes.
3. **Connect Power:** Plug the power adapter into the device and then into an electrical outlet. The device is corded and not rechargeable.



Figure 5: Visual guide for preparing and using the XSOUL IPL device.

5. Operating Instructions

1. **Power On:** Press and hold the power button (usually the main button on the device) for a few seconds until the display lights up and the fan starts running.
2. **Select Energy Level:** The device has 5 energy levels. Short press the power button to cycle through the levels. For first-time users, begin with Level 1 and gradually increase as your skin tolerates it.
3. **Choose Operation Mode:** The device offers two modes:
 - **Manual Mode:** Short press the flash button for a single flash. Ideal for precise treatment on smaller, sensitive areas like the bikini line, face, and armpits.
 - **Auto Mode:** Press and hold the flash button for continuous flashing. This mode is suitable for larger areas such as the back, arms, and legs.
4. **Begin Treatment:** Place the device's flash window firmly against the skin. The device will only flash when it detects full skin contact. Press the flash button (or hold for Auto Mode) to emit a pulse of light. Move the device to the next untreated area and repeat. Avoid treating the same area twice in one session.

IS THE XSOUL IPL HAIR REMOVAL RIGHT FOR YOU?

SKIN TONE & BODY HAIR TYPE CHART

Skin Tone							
Hair Color							
White/Grey		✗	✗	✗	✗	✗	✗
Light Blonde		✗	✗	✗	✗	✗	✗
Red		✗	✗	✗	✗	✗	✗
Dark Blonde Light Brown		✓	✓	✓	✓	✗	✗
Brown		✓	✓	✓	✓	✗	✗
Dark Brown		✓	✓	✓	✓	✗	✗
Black		✓	✓	✓	✓	✗	✗

 = It's match!
  = Try a different hair removal method

* Not applicable for tatooed skin and tanned skin

If you are not sure of your skin/hair type, please refer to the Skin Tone & Body Hair Type Chart prior to usage.

Note: different skin sensitivity and exposure settings may lead to different user experience. Start from Power Level 1 and work your way up as needed.



Figure 6: Manual and Auto modes for targeted and broad area hair removal.

Safe, Fast & Effective

EFFICIENT REDUCTION IN HAIR REGROWTH



Before



After




Week 1-3

1st-5th week: 2-3 times/week



Week 4-5



Week 6-7



After Week 8

6th-8th week: 1 time/week

Figure 7: Device interface displaying energy levels, smart sensor, and mode selection.

Your browser does not support the video tag.

Video 1: Demonstration of the XSOUL IPL device in use for hair removal, highlighting its ease of use and effectiveness.

Your browser does not support the video tag.

Video 2: User demonstrating the application of the XSOUL IPL device on the skin, showing the flashing mechanism.

6. Usage Schedule and Hair Removal Plan

For optimal results, follow the recommended treatment schedule:

- **Weeks 1-5:** Use the device 2-3 times per week.
- **Weeks 6-8:** Use the device 1 time per week.
- **After Week 8:** Continue with maintenance treatments as needed, typically once a month or as hair regrowth slows.

Consistency is key to achieving significant hair reduction. Individual results may vary.



Figure 8: Visual representation of hair reduction over an 8-week treatment period.

7. Maintenance

To ensure the longevity and optimal performance of your device:

- **Cleaning:** After each use, gently wipe the flash window with a soft, dry cloth to remove any hair residue or debris. Ensure the device is unplugged before cleaning.

- **Storage:** Store the device in a cool, dry place away from direct sunlight and moisture.

8. Troubleshooting

If you encounter any issues with your XSOUL IPL Hair Removal Device, please refer to the following common troubleshooting tips:

- **Device Not Flashing:** Ensure the device is firmly pressed against the skin, making full contact. Check if the power adapter is securely plugged into both the device and the electrical outlet.
- **Overheating:** If the device feels excessively warm, allow it to cool down for a few minutes before resuming use. This is a normal protective measure.
- **Skin Irritation:** If you experience redness or discomfort, reduce the energy level or discontinue use. Ensure your skin is clean and dry before treatment.

For further assistance, please contact customer support.

9. Specifications

Feature	Detail
Product Dimensions	8.27 x 5.51 x 3.66 inches
Item Model Number	D-1103
Weight	8.96 ounces
UPC	194358027205
Manufacturer	XSOUL

10. Warranty and Support

XSOUL provides comprehensive after-sales service and product guarantee to ensure your satisfaction.

- **Return Policy:** This product is eligible for a refund/replacement for up to 2026 days from purchase.
- **Customer Support:** For any questions, concerns, or assistance, please contact us at help@xsoulofficial.com.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0828JD1WF>

More documents and user manuals: <https://manuals.plus/>