

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

- › [ROSA RUGOSA](#) /
- › [ROSA RUGOSA Menstrual Cup Sterilizer: User Manual](#)

ROSA RUGOSA Menstrual Cup Sterilizer

ROSA RUGOSA Menstrual Cup Sterilizer: User Manual

Model: Menstrual Cup Sterilizer

INTRODUCTION

The ROSA RUGOSA Menstrual Cup Sterilizer is a 3-in-1 device designed for efficient and hygienic cleaning, drying, and storage of menstrual cups. This device provides a chemical-free sterilization process, eliminating up to 99.9% of germs, and features an automatic shut-off for convenience and safety.

SAFETY INFORMATION

- **Anti-Scalding Design:** The sterilizer is constructed from superior PP material, offering high-temperature resistance. Its double-walled design helps prevent burns by keeping the exterior cool during operation.
- **Automatic Shut-Off:** The device is equipped with an auto shut-off feature that activates after 8-10 minutes or when the water dries, preventing dry boiling and ensuring safety.
- Always ensure the lid is properly closed before starting the sterilization cycle.
- Use only clean water for sterilization. Distilled water is recommended to prevent mineral buildup.
- Keep out of reach of children.

PRODUCT OVERVIEW

The ROSA RUGOSA Menstrual Cup Sterilizer combines sterilizing, secondary cleaning, and storage functions. Its compact design ensures privacy and convenience.



Image: Product features including automatic sterilization, container function, privacy protection, and double heat insulation.

The device is made from Bisphenol-A (BPA) Free, Latex Free, Chlorine Free, and Phthalate Free materials, ensuring safety for your menstrual cup.

SETUP

1. Unpack the sterilizer and ensure all components are present.
2. Place the sterilizer on a stable, flat surface near a power outlet.
3. Ensure the device is clean before first use.

OPERATING INSTRUCTIONS

Follow these steps for effective sterilization of your menstrual cup:

1. **Add Water:** Remove the lid and add water to the water level mark inside the heating plate.
2. **Place Menstrual Cup:** Place your clean menstrual cup upside down onto the central stand within the sterilizer.
3. **Close Lid:** Securely place the lid back onto the sterilizer.
4. **Start Sterilization:** Plug in the device and press the power button. The indicator light will illuminate.
5. **Wait:** The device will automatically sterilize for 8-10 minutes. Steam will be visible from the top vent.
6. **Automatic Shut-Off:** The sterilizer will automatically turn off once the cycle is complete or the water has evaporated. The indicator light will turn off.

7. **Remove Cup:** Carefully remove the lid and your sterilized menstrual cup. Be cautious as the cup and internal components will be hot.



Image: Step-by-step guide for using the sterilizer.

Your browser does not support the video tag.

Video: Official demonstration of the Rosa Rugosa Menstrual Cup Sterilizer in use.

MAINTENANCE

- **Cleaning:** After each use, allow the device to cool down. Wipe the heating plate and interior with a soft, damp cloth. Avoid abrasive cleaners.
- **Descaling:** If mineral deposits accumulate on the heating plate (especially with hard water), use a descaling solution or a mixture of water and vinegar to clean it. Rinse thoroughly.
- **Storage:** The sterilizer doubles as a storage container for your menstrual cup, offering privacy and protection. Ensure the cup is dry before storing.

TROUBLESHOOTING

- **Device not turning on:** Ensure the device is properly plugged into a working power outlet. Check if the power button is pressed.
- **No steam or insufficient steam:** Verify that enough water has been added to the heating plate, up to the indicated level. Ensure the lid is securely closed.
- **Mineral buildup:** If you notice white residue on the heating plate, descale the unit as described in the Maintenance section. Using distilled water can prevent this.

SPECIFICATIONS

Product Dimensions	6.34 x 3.98 x 3.94 inches
Weight	9.91 ounces
Manufacturer	ROSA RUGOSA
Material Type Free	Bisphenol-A (BPA) Free, Latex Free, Chlorine Free, Phthalate Free
Color	Pink

SUITABLE FOR ALL MENSTRUAL CUPS ON THE MARKET:



Image: Product dimensions and compatibility with various menstrual cup designs.

WARRANTY AND SUPPORT

The ROSA RUGOSA Menstrual Cup Sterilizer comes with a 1-year hassle-free warranty and a 30-day money-back guarantee. For any support or inquiries, please contact our friendly service team.

