



[Manuals.plus](#) /

> [Mia&Coco](#) /

> Mia&Coco Electric Heating Pad User Manual

Mia&Coco HP05-001B

Mia&Coco Electric Heating Pad User Manual

Model: HP05-001B | Brand: Mia&Coco

1. INTRODUCTION

Thank you for choosing the Mia&Coco Electric Heating Pad. This product is designed to provide soothing heat therapy for various body parts, including the back, neck, shoulders, and legs. Its soft polyester material and user-friendly controls ensure a comfortable and effective experience. Please read this manual thoroughly before use to ensure safe and proper operation.

2. IMPORTANT SAFETY INSTRUCTIONS

When using electrical appliances, basic precautions should always be followed to reduce the risk of fire, electric shock, and personal injury. These heating pads are manufactured from high-quality materials and are designed for safe use when operated according to these instructions.

- Do not use while sleeping.
- Do not use on an infant, invalid, or a person insensitive to heat.
- Do not use with liniment, salve, or ointment preparations, as skin burns could result.
- Never pull this pad by the supply cord.
- Do not use if the cord or plug is damaged.
- Do not use if the pad is wet.
- Always disconnect the pad from the electrical outlet when not in use.
- The heating pad features an automatic shut-off function, turning off after approximately 90 minutes for safety.

Sicher mit ETL-Zertifikat zu verwenden

Mia & Coco Heizkissen sind aus hochwertigem Material gefertigt und von hoher Qualität



Close-up view of the soft polyester fabric and the ETL safety certification label on the heating pad, indicating compliance with safety standards.

3. SETUP

Before first use, unpack the heating pad and its controller. Ensure all components are present and undamaged.

1. Unroll the heating pad completely. Do not fold or bunch the pad during use.
2. Connect the controller cable securely to the heating pad's connector. Ensure it clicks into place.
3. Plug the power cord into a standard electrical outlet.



The Mia&Coco Electric Heating Pad, showing the pad, controller, and product packaging.

4. OPERATING INSTRUCTIONS

The Mia&Coco Electric Heating Pad features an easy-to-use controller with multiple heat settings and an automatic shut-off function.

Powering On/Off

- To turn on the heating pad, press the power button on the controller. The display will light up.
- To turn off the heating pad, press the power button again until the display turns off.

Adjusting Temperature Settings

The controller allows you to select from 3 heat levels to customize your comfort.

- Press the temperature adjustment button (usually indicated by an up/down arrow or a 'level' icon) to cycle through the available heat settings.

- The current heat level will be displayed on the controller's screen.

Automatic Shut-Off

For your safety and energy conservation, the heating pad is equipped with an automatic shut-off feature. The pad will automatically turn off after approximately 90 minutes of continuous use.



Drei Heizstufen & Auto-OFF

Aus Sicherheitsgründen wird es

in **90 Minuten** automatisch ausgeschaltet

User holding the controller, illustrating the three heat level settings and power button for easy operation.



The heating pad can be comfortably applied to the lower back for targeted heat therapy.



The heating pad can also be used on the shoulder area for soothing relief.

Weiche und komfortable Berührung

weich und warm zum Berühren und Bedecken von Rücken
und Nacken, Schultern, Beinen oder Bauch



For abdominal or hip comfort, the heating pad can be placed while relaxing.



Versatile application: The heating pad can be used on various body parts including the back, shoulders, and knees, and can be easily folded for storage.

5. MAINTENANCE

Proper care and maintenance will extend the life of your Mia&Coco Electric Heating Pad.

Cleaning Instructions

The heating pad is machine washable for easy cleaning.

1. Always disconnect the power cord from the electrical outlet before cleaning.
2. Detach the controller cable from the heating pad.
3. Machine wash the pad on a gentle cycle with cold water (max 30°C / 86°F).
4. Do not bleach or iron the heating pad.
5. Do not dry clean.

Drying and Storage

- Tumble dry on low heat or line dry. Ensure the pad is completely dry before reconnecting the controller or storing.
- Store the heating pad in a cool, dry place, preferably in its original packaging or a protective bag.
- Do not wrap the cord tightly around the controller or pad, as this can damage the cord.

6. TROUBLESHOOTING

If you encounter any issues with your heating pad, please refer to the following common solutions:

Problem	Possible Cause	Solution
Pad not heating	Not plugged in; Controller not connected properly; Power button not pressed.	Ensure the power cord is securely plugged into a working outlet. Check that the controller cable is firmly connected to the pad. Press the power button to turn on.
Pad heats unevenly	Pad is folded or bunched during use.	Ensure the heating pad is laid flat and not folded or bunched during operation.
Controller display is off	Automatic shut-off activated; Power disconnected.	The pad may have automatically turned off after 90 minutes. Press the power button to restart. Check power connection.

If the problem persists after trying these solutions, please contact customer support.

7. SPECIFICATIONS

Feature	Detail
Model Number	HP05-001B
Dimensions	Approx. 30 x 60 cm (12 x 24 inches)
Material	100% Polyester
Voltage	230 Volts
Frequency	60 Hz
Power	60 Watts
Heat Settings	3 Levels
Automatic Shut-off	Approx. 90 minutes

8. WARRANTY AND SUPPORT

Mia&Coco products are designed and manufactured to high-quality standards. For any questions, concerns, or warranty claims, please contact Mia&Coco customer support through the retailer's platform or the official brand website. Please retain your proof of purchase for warranty purposes.

