

[Manuals.plus](#) /

› [Travel Buggy](#) /

› Travel Buggy Lightweight Adjustable Telescopic Folding Retractable Wheelchair Ramp User Manual

Travel Buggy Telescopic Folding Retractable Wheelchair Ramp

Travel Buggy Lightweight Adjustable Telescopic Folding Retractable Wheelchair Ramp User Manual

Model: Telescopic Folding Retractable Wheelchair Ramp

1. INTRODUCTION

Thank you for choosing the Travel Buggy Lightweight Adjustable Telescopic Folding Retractable Wheelchair Ramp. This manual provides essential information for the safe and effective use, setup, maintenance, and troubleshooting of your new ramp. Please read these instructions carefully before initial use and retain them for future reference.

2. SAFETY INFORMATION

Your safety is paramount. Adhere to all safety guidelines to prevent injury or damage.

- **Maximum Weight Capacity:** Do not exceed the maximum load capacity of **600 lbs (272 kg)**. Overloading can cause structural failure and severe injury.
- **Proper Angle:** Ensure the ramp is placed at a safe and gradual incline. Avoid steep angles. Consult your mobility device's manual for recommended ramp angles.
- **Secure Placement:** Always ensure both ends of the ramp are securely resting on stable, level surfaces before use. The top lip should be firmly seated against the vehicle or elevated surface.
- **Non-Skid Surface:** The ramp features a non-skid traction surface. Keep this surface clean and free from debris, ice, or water to maintain optimal grip.
- **Side Guardrails:** The integrated side guardrails are designed to help prevent wheels from slipping off the edge. Always center the mobility device on the ramp.
- **Supervision:** Always have an attendant present when a user is operating a manual or power wheelchair on the ramp.
- **Inspection:** Before each use, inspect the ramp for any signs of damage, wear, or loose components. Do not use if damaged.
- **Storage:** Store the ramp in a dry, protected area when not in use.

3. PACKAGE CONTENTS

Verify that all items are present in your package:

- Telescopic Wheelchair Ramps (2 units)
- Carry Bag (1 unit)



Image: The included black carry bag for convenient transport and storage of the ramps.

4. SETUP INSTRUCTIONS

Follow these steps to properly set up your telescopic wheelchair ramps:

1. **Unpack:** Carefully remove the ramps from the carry bag.
2. **Position:** Place the top lip of each ramp securely onto the edge of the elevated surface (e.g., vehicle trunk, doorway threshold). Ensure the ramps are parallel and spaced appropriately for your mobility device's wheel width.
3. **Extend:** While holding the ramp firmly, gently pull the inner section outwards to extend the ramp to the desired length. The ramp has two locking positions: 4 feet and 7 feet. Ensure the locking pins engage fully at your chosen length. You will hear a click when the pins are engaged.
4. **Verify Stability:** Before use, apply light pressure to each ramp to confirm it is stable and will not shift or slide.



Image: The telescopic ramps shown in both their folded and fully extended configurations.



Image: A wheelchair being loaded onto the ramps, demonstrating proper positioning for vehicle access.

5. OPERATING INSTRUCTIONS

Once the ramps are securely set up, follow these steps for safe operation:

1. **Approach:** Position the mobility device directly in front of the ramps, aligning the wheels with the center of each ramp.
2. **Ascending:** Slowly and steadily push or drive the mobility device up the ramps. Maintain a consistent speed and avoid sudden movements. Ensure the attendant maintains control and stability.
3. **Descending:** When descending, back the mobility device down the ramps slowly and carefully. An attendant should guide the device, ensuring it remains centered and stable.
4. **Clearance:** Be mindful of head clearance and any obstacles when entering or exiting a vehicle or doorway.
5. **Removal:** After use, retract the ramps by pressing the locking pins and gently pushing the inner section back into the outer section. Store them in the provided carry bag.



Image: A wheelchair being loaded into a car trunk using the ramps, showing the gradual incline.



Image: A wheelchair being loaded into a minivan, illustrating the ramp's versatility for different vehicle types.



Image: The wheelchair fully inside the minivan, with the ramps providing a stable path.

6. MAINTENANCE

Regular maintenance ensures the longevity and safe operation of your ramps.

- **Cleaning:** Clean the ramps regularly with mild soap and water. Rinse thoroughly and dry completely to prevent corrosion. Do not use abrasive cleaners or solvents.
- **Inspection:** Periodically check all components for wear, damage, or loose parts. Pay close attention to the locking pins, non-skid surface, and structural integrity of the aluminum.
- **Lubrication:** If the extension mechanism becomes stiff, apply a small amount of silicone-based lubricant to the sliding surfaces. Avoid oil-based lubricants that can attract dirt.
- **Storage:** Always store the ramps in their carry bag in a dry, protected environment away from extreme temperatures and direct sunlight.

7. TROUBLESHOOTING

If you encounter issues with your ramps, refer to the following common problems and solutions:

- **Problem:** Ramp extension is difficult or stiff.
Solution: Check for any debris or obstructions in the sliding mechanism. Clean the ramps thoroughly. If still stiff, apply a silicone-based lubricant to the sliding surfaces as described in the Maintenance section.
- **Problem:** Locking pins do not engage properly.
Solution: Ensure the ramp is fully extended to one of the two designated locking positions (4 ft or 7 ft). Check if the locking pins are bent or obstructed. Do not force the extension; if pins are damaged, discontinue use and contact customer support.
- **Problem:** Ramp feels unstable during use.
Solution: Re-verify that both ends of the ramp are securely placed on stable, level surfaces. Ensure the top lip is firmly seated. Confirm the ramps are parallel and correctly spaced. Do not use if instability persists.

8. SPECIFICATIONS

Feature	Detail
Brand	Travel Buggy
Material	Aluminum
Color	Silver

Item Weight	8.62 kg (19 lbs)
Product Dimensions (Extended)	2.13L x 0.15W x 0.08H Meters (7 ft x 5.9 in x 3.1 in)
Load Capacity	600 lbs
Foldable	Yes
Inner Surface Width	5.2 inches
UPC	638126208217

9. WARRANTY AND SUPPORT

For warranty information or technical support, please refer to the documentation included with your purchase or contact Travel Buggy customer service directly. You can often find contact details on the manufacturer's official website or through your retailer. Visit the [Travel Buggy Store](#) for more products and information.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B08HVVH1NF>

More documents and user manuals: <https://manuals.plus/>