

[Manuals.plus](#) /

› [Breo](#) /

› Breo WOWOS Hand Massager User Manual

Breo WOWOS

Breo WOWOS Hand Massager User Manual

Model: WOWOS

INTRODUCTION

The Breo WOWOS Hand Massager is designed to provide targeted massage and heat therapy for hands and fingers. This manual provides essential information for the safe and effective operation, maintenance, and troubleshooting of your device.



Image: The Breo WOWOS Hand Massager, a white, oval-shaped device with a central display, shown alongside a smartphone illustrating the companion app's hand massage interface.

SAFETY INFORMATION

Adherence to the following safety guidelines is crucial for preventing injury and ensuring proper device function:

- Do not use on broken skin or open wounds.
- Keep away from heat sources.
- Do not submerge in water.
- Allow the massager to cool down before storing.
- Always follow the instructions provided in the user manual for safe and effective use.
- Women who are pregnant or in menstrual period should not use the massager.
- Do not use the product if the massage position has skin injury, redness and inflammation, congestion or other discomfort symptoms.
- Persons with heart disease or history of heart disease, or with acute disease, skin disease, malignant tumor or physical abnormality should not use the product.
- Persons with bleeding tendency or other blood diseases, or with osteoporosis, soft-tissue disorders or apoplexia should not use the product.
- Stop using the product in case of any discomfort during use.

WHAT'S IN THE BOX

Upon unpacking, verify that all components are present:

- WOWOS Massager
- Micro-USB Cable
- Plastic Gloves
- Charger
- User Manual (this document)



Image: A flat lay showing the Breo WOWOS Hand Massager, its user manual, a micro-USB cable, a charger, and a pair of plastic gloves.

PRODUCT OVERVIEW

The Breo WOWOS Hand Massager features an ergonomic design with intuitive controls and a clear display.

LCD SCREEN



6 Modes

- Default Mode 1-5
- APP DIY Mode 6



Force 1-3



Power Button



Start/Pause



Safety Button

Image: A close-up of the Breo WOWOS Hand Massager's control panel, highlighting the LCD screen and the function of each button: Mode, Force, Power, Start/Pause, and Safety Button.

Key Features:

- **Pressure and Heat Technology:** Combines pressure and heat to relax acupoints with a waving rhythm, relieving hand and finger pain and coldness.
- **Reflexology Massage:** Applies pressure to 36 reflective areas of the palm through rhythmic massage, relaxing the palm and fingers.
- **Fingertip Heating:** Effective fingertip heating for muscle relaxation and relief from hand fatigue, discomfort, and cold fingers.
- **APP Control:** Personalize settings using the Breo APP with 5 default massaging modes, 3 different intensities, and a customizable mode.
- **Cordless Design:** Allows for convenient use anywhere once charged.
- **Non-slip Mat:** Integrated texture at the bottom of the hand insertion area helps keep the hand stable during massage.

HEAT & POINT MASSAGE

Acupuncture
Point Massage



Non-slip mat

Image: A diagram illustrating the heat and acupuncture point massage function of the Breo WOWOS Hand Massager, showing internal components and a non-slip mat.

SETUP

1. Charging the Device:

Connect the Micro-USB cable to the charging port located on the device (refer to the diagram in the Product Overview section for port location). Connect the other end to the provided charger or a compatible USB power source. The battery indicator on the display will show charging status. Ensure the device is unplugged before turning it on for use.



Image: A graphic highlighting the long-lasting battery feature of the Breo WOWOS Hand Massager, showing a battery icon.

2. Initial Power On:

Press and hold the power button for approximately 3 seconds to turn on the device. The display will illuminate, showing the current mode and intensity settings.

OPERATING INSTRUCTIONS

The Breo WOWOS Hand Massager offers various modes and intensities for a personalized experience.

1. Placing Your Hand:

Insert your hand into the massager, ensuring your palm rests comfortably on the internal texture. The device is designed to accommodate various hand sizes.

2. Selecting a Mode:

Use the 'Mode' button (circular arrow icon) to cycle through the available massage modes. The display will illuminate, showing a hand icon with highlighted areas indicating which part of the hand is being massaged in the selected mode.

3. Adjusting Intensity:

Use the '+' and '-' buttons to adjust the massage intensity. There are 3 intensity levels available.

4. Starting and Pausing:

Press the 'Start/Pause' button (play/pause icon) to begin or pause the massage cycle.

5. Air Release:

To quickly release air pressure and remove your hand, press the 'Air Release' button (wavy lines icon).

6. APP Control for Personalized Settings:


Download the Breo APP to unlock customizable massage settings. The app allows you to:

- Tap on the hand icons to switch modes.
- Select from "soft", "smoothing" or "strong" for different intensities.
- Enter customization mode by tapping "custom". Select "left hand" or "right hand" accordingly. Tap on the numbers (1-5) to customize massage areas.
- Tap "Heating" to turn warm compression ON/OFF.
- Tap on the icon on upper right corner to customize the working duration (15min, 10min, 5min).

INTELLIGENT APP OPERATED

Multi-Modes Options & Personal Customized Mode

Multi-Mode Options



Personal Customized Mode



Image: Three smartphone screens displaying the Breo WOWOS app interface, showing multi-mode options and a personal customized mode for hand massage.



Image: A series of five smartphone screens illustrating the different massage modes and customization options available through the Breo WOWOS app, including soft, soothing, strong, and heating settings.

MAINTENANCE

Proper maintenance ensures the longevity and optimal performance of your massager:

- **Cleaning:** Wipe the exterior of the device with a soft, damp cloth. Do not use abrasive cleaners or submerge the device in water.
- **Storage:** Store the massager in a cool, dry place away from direct sunlight and extreme temperatures.
- **Hand Care with Gloves:** For enhanced hand care, apply hand cream and wear the provided plastic gloves before using the massager. This helps with absorption and keeps the device interior clean.



01. Tap on the hand icons to switch mode.

02. Select from "soft", "soothing" or "strong" for different intensities.

03. Enter customization mode by tapping "custom".

Select "left hand" or "right hand" accordingly.

Tap on the numbers (1-5) to customize massage areas.

04. Tap "Heating" to turn warm compression ON/OFF.

05. Tap on the icon on upper right corner to customize the working duration (15min, 10min, 5min).

Image: A three-step visual guide for hand care using the massager: 1) Slathering Hand Cream, 2) Wearing Gloves, and 3) Applying Heat with the massager.

TROUBLESHOOTING

If you encounter any issues with your Breo WOWOS Hand Massager, refer to the following common solutions:

Problem	Possible Cause	Solution
Device does not turn on.	Low battery or still charging.	Ensure the device is fully charged and unplugged from the power source. Press and hold the power button for 3 seconds.
Heating function is not working.	Heating takes time to activate.	The heating function may take 8-12 minutes to heat up. Ensure heating is enabled via the device controls or the Breo APP.
Massage feels too weak or too strong.	Incorrect intensity setting.	Adjust the intensity using the '+' and '-' buttons on the device or through the Breo APP.
Device stops unexpectedly.	Automatic shut-off or low battery.	The massager has an automatic shut-off after 30 minutes. Recharge the device if the battery is low.

If issues persist, please contact Breo Service Team for further assistance.

SPECIFICATIONS

- **Product Dimensions:** 11.42 x 4.72 x 7.09 inches
- **Item Model Number:** WOWOS
- **Batteries:** 1 Lithium Ion battery required (included)
- **Date First Available:** December 10, 2020
- **Manufacturer:** breo
- **ASIN:** B08Q87CMLV
- **Use for:** Hands
- **Power Source:** Battery Powered
- **Material:** Plastic
- **Item Weight:** 1 Kilograms (2.2 Pounds)
- **Brand:** Breo



Image: A diagram illustrating the dimensions of the Breo WOWOS Hand Massager, showing its height, width, and depth in inches, with a hand inserted for scale.

WARRANTY & SUPPORT

Breo offers comprehensive support for your WOWOS Hand Massager:

- **Satisfaction Guarantee:** 100% satisfaction guarantee within 30 days for any issues.
- **Warranty:** 12-months warranty from the date of purchase.
- **Customer Service:** Friendly customer service is available for assistance. Please refer to the contact information provided in your product packaging or on the official Breo website for support.

© 2023 Breo. All rights reserved.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B08Q87CMLV>

More documents and user manuals: <https://manuals.plus/>