

Cubitt Tracking Technologies B08N6YYJSJT

CUBITT Smart Digital Weight Scale Instruction Manual

Model: B08N6YYJSJT

INTRODUCTION

Thank you for choosing the CUBITT Smart Digital Weight Scale. This device is designed to provide comprehensive body composition analysis, helping you track your health and fitness progress. It measures 13 essential body metrics and syncs data with your smartphone via Bluetooth for easy monitoring.

SAFETY PRECAUTIONS

- Do not use if you have a pacemaker or other internal medical device.
- Do not use if you are pregnant.
- The scale should not be used to diagnose or treat any medical condition. Always consult a healthcare professional for medical advice.
- Place the scale on a hard, flat surface for accurate measurements. Avoid carpets or uneven floors.
- Keep the scale dry. Do not immerse it in water or use harsh chemical cleaners.
- Stand still on the scale with bare feet for best results.

PACKAGE CONTENTS

- CUBITT Smart Digital Weight Scale
- 2 x AAA Batteries
- Instruction Manual




Image: The retail packaging for the CUBITT Smart Digital Weight Scale, showing the scale and a smartphone displaying the app interface.

PRODUCT OVERVIEW

The CUBITT Smart Digital Weight Scale features a sleek design with a clear LED display and four high-precision sensors. The scale uses Bio-electrical Impedance Analysis (BIA) technology to provide detailed body composition measurements.



Image: A top-down view of the black CUBITT Smart Digital Weight Scale, showing the LED display and four circular sensor pads.



Image: An angled view of the black CUBITT Smart Digital Weight Scale, highlighting its slim profile and sensor pads.

INITIAL SETUP

1. Battery Installation

1. Open the battery compartment on the underside of the scale.
2. Insert the two AAA batteries, ensuring correct polarity (+/-).
3. Close the battery compartment securely.

2. App Download and Installation

The CUBITT Smart Digital Weight Scale requires the Cubitt Health App for full functionality and tracking of body composition data. The app is compatible with Android and iOS devices.

- Search for "Cubitt Health" in the [Apple App Store](#) (for iOS devices) or [Google Play Store](#) (for Android devices).

- Download and install the app on your smartphone.
- Register for a new account or log in if you already have one.

3. Pairing with Your Device

1. Ensure Bluetooth is enabled on your smartphone.
2. Open the Cubitt Health App.
3. Step onto the scale briefly to activate it. The display will show "0.00".
4. Follow the in-app instructions to connect the scale. The app will automatically detect the scale and prompt you to pair.
5. Once paired, your weight and other body composition data will sync automatically to the app after each measurement.

13 essential measures



body water



skeletal muscle



muscle mass



fat-free body weight



metabolic age

Image: A person smiling while looking at their smartphone, which displays the Cubitt Health App interface, demonstrating the connection process.

OPERATING THE SCALE

1. Taking a Measurement

1. Place the scale on a hard, flat surface.
2. Open the Cubitt Health App on your smartphone (optional, but recommended for full data tracking).
3. Step onto the scale with bare, dry feet, ensuring your feet are positioned on the four sensor pads.
4. Stand still until the weight reading stabilizes on the LED display.
5. If the app is open and connected, your weight and body composition data will automatically sync and be displayed in the app. If the app is not open, only your weight will be shown on the scale's display. Data will upload to the app when you open it and measure again.

2. Understanding Your Measurements

The Cubitt Health App provides 13 essential body composition measurements:

- **Weight:** Your total body mass.
- **BMI (Body Mass Index):** A measure of body fat based on height and weight.
- **Body Fat:** The percentage of fat in your body.
- **Visceral Fat:** Fat stored around your organs.
- **Body Water:** The percentage of water in your body.
- **Muscle Mass:** The total amount of muscle in your body.
- **Bone Mass:** The amount of bone tissue in your body.
- **BMR (Basal Metabolic Rate):** The number of calories your body needs at rest.
- **Skeletal Muscle:** Muscle attached to bones, responsible for movement.
- **Fat-Free Body Weight:** Your total body weight minus body fat.
- **Metabolic Age:** An indicator of your body's internal age compared to your chronological age.
- **Protein:** The percentage of protein in your body.
- **Subcutaneous Fat:** Fat stored directly under the skin.

Smart Scale



Measure
your
achievements

Image: An infographic illustrating the 13 essential body composition measurements provided by the scale, including body water, skeletal muscle, muscle mass, fat-free body weight, and metabolic age.

*Follow up
with your
body*



Body fat



Body weight



Muscle mass



Bone mass



Skeletal muscle



Body water

Image: A screenshot of the Cubitt Health App interface displaying various body composition measurements and progress tracking.

3. Multiple User Profiles

The Cubitt Health App supports an unlimited number of user profiles. Each user can create their own account within the app to track their individual progress. To add a new user:

1. Open the Cubitt Health App.
2. Navigate to the profile or user management section (usually found in settings).
3. Select "Add User" or "Create New Profile" and follow the prompts to enter their details.
4. When a new user steps on the scale, ensure their profile is selected in the app to record their data correctly.

CARE AND MAINTENANCE

- **Cleaning:** Wipe the scale surface with a soft, damp cloth. Do not use abrasive cleaners or immerse the scale in water.
- **Storage:** Store the scale in a cool, dry place. Avoid extreme temperatures and direct sunlight.
- **Battery Replacement:** Replace batteries when the low battery indicator appears on the display.
- **Avoid Impact:** Do not drop the scale or subject it to heavy impact, as this may damage the sensors.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Scale does not turn on.	Batteries are dead or incorrectly installed.	Check battery polarity or replace with new AAA batteries.
Inaccurate weight readings.	Scale is on an uneven surface; not standing still; scale needs calibration.	Place scale on a hard, flat surface. Stand still. Step on and off the scale quickly to re-calibrate to 0.00.
App not syncing data.	Bluetooth is off; app is not open; scale not paired; app permissions.	Ensure Bluetooth is on. Open the Cubitt Health App. Re-pair the scale if necessary. Check app permissions for Bluetooth access.
Body composition data not showing.	Not standing with bare feet; feet not on sensors; app not connected.	Ensure bare, dry feet are on all four sensor pads. Confirm app is open and connected during measurement.

PRODUCT SPECIFICATIONS

Feature	Specification
Brand	Cubitt Tracking Technologies
Model	B08N6YYJSJT
Color	Black
Display Type	LED
Weight Limit	400 Pounds
Readout Accuracy	High (BIA Technology)

Feature	Specification
Measurements	13 Body Composition Metrics
Connectivity	Bluetooth
App Compatibility	Android, iOS (Cubitt Health App)
Power Source	2 x AAA Batteries (included)
Product Dimensions	0.98 x 10.63 x 10.63 inches
Item Weight	2.9 pounds

WARRANTY INFORMATION

CUBITT Tracking Technologies provides a limited warranty for this product against defects in materials and workmanship under normal use. Please refer to the official Cubitt website or contact customer support for detailed warranty terms and conditions.

CUSTOMER SUPPORT

For technical assistance, troubleshooting, or any questions regarding your CUBITT Smart Digital Weight Scale, please contact our customer support team:

- **Website:** Visit the [Cubitt Tracking Technologies Store on Amazon](#) for product information and support resources.
- **Email:** Refer to your product packaging or the Cubitt website for specific support email addresses.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B08R6L1CG6>

More documents and user manuals: <https://manuals.plus/>