

[Manuals.plus](#) /

› [Fisher & Paykel](#) /

› Fisher & Paykel Brevida Nasal Pillows Mask Elbow Instruction Manual

Fisher & Paykel Brevida Nasal Pillows Mask Elbow

Fisher & Paykel Brevida Nasal Pillows Mask Elbow Instruction Manual

Essential guide for the proper use and care of your Brevida Nasal Pillows Mask Elbow.

INTRODUCTION

This manual provides instructions for the Fisher & Paykel Brevida Nasal Pillows Mask Elbow. This component is designed to connect the breathing tube to the frame of the Brevida Nasal Pillows Mask, facilitating airflow during CPAP therapy. Please read these instructions carefully before use to ensure proper assembly, operation, and maintenance.

PRODUCT OVERVIEW



Figure 1: The Fisher & Paykel Brevida Nasal Pillows Mask Elbow. This image shows the elbow component, which serves as the connection point between the CPAP breathing tube and the mask frame.

The Brevida Nasal Pillows Mask Elbow is a crucial part of your CPAP mask system, ensuring a secure and comfortable connection for effective therapy. It is specifically designed for compatibility with the Fisher & Paykel Brevida Nasal Pillows Mask.

SETUP AND ASSEMBLY

1. **Inspect Components:** Before assembly, carefully inspect the elbow for any signs of damage, wear, or foreign objects. Do not use if damaged.
2. **Connect to Mask Frame:** Align the elbow's smaller end with the opening on the Brevida Nasal Pillows Mask frame. Gently push and twist until it is securely seated. Ensure a snug fit to prevent air leaks.
3. **Connect to Breathing Tube:** Attach the larger end of the elbow to your CPAP machine's breathing tube. Push the tube firmly onto the elbow connector until it is secure.

4. **Verify Connection:** Once assembled, gently tug on both the mask frame and the breathing tube to confirm that the elbow is firmly attached at both ends.

Note: Ensure all connections are secure to maintain optimal CPAP therapy pressure and prevent air leakage.

OPERATING INSTRUCTIONS

The Brevida Nasal Pillows Mask Elbow functions as a passive connector within your CPAP system. Once properly assembled, it requires no further operational steps. Its primary role is to provide a flexible and secure conduit for pressurized air from your CPAP machine to the nasal pillows mask.

- Ensure the elbow is clean and free from obstructions before each use.
- Confirm that the elbow is securely attached to both the mask frame and the CPAP breathing tube.
- Periodically check for any signs of wear or damage during regular use.

MAINTENANCE AND CLEANING

Regular cleaning of your Brevida Nasal Pillows Mask Elbow is essential for hygiene and optimal performance. Follow these steps:

1. **Daily Cleaning:** Disconnect the elbow from the mask frame and breathing tube. Wash the elbow in warm, soapy water using a mild detergent.
2. **Rinse Thoroughly:** Rinse the elbow completely under running water to remove all soap residue.
3. **Air Dry:** Allow the elbow to air dry completely out of direct sunlight before reassembling. Ensure it is thoroughly dry to prevent bacterial growth.
4. **Inspection:** During cleaning, inspect the elbow for any cracks, tears, or discoloration. Replace the elbow if any damage is observed.

Do not use: Bleach, alcohol, or harsh cleaning solutions, as these can damage the material and reduce the lifespan of the elbow.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Air Leakage at Elbow Connection	Elbow not fully seated; elbow or mask frame damaged; elbow worn out.	Ensure elbow is pushed firmly onto both the mask frame and breathing tube. Inspect for cracks or tears; replace if damaged or worn.
Difficulty Attaching/Detaching Elbow	Misalignment; residue buildup; new component stiffness.	Ensure correct alignment. Clean connections thoroughly. Apply gentle, consistent pressure.
Whistling Noise from Elbow	Small air leak; foreign object partially obstructing airflow.	Check all connections for tightness. Inspect the elbow for any debris or damage.

SPECIFICATIONS

- **Product Name:** Fisher & Paykel Brevida Nasal Pillows Mask Elbow

- **Compatibility:** Specifically designed for Fisher & Paykel Brevida Nasal Pillows Mask
- **Quantity:** 1 (Elbow only, other mask parts sold separately)
- **Material:** Medical-grade plastic (assumed, based on typical CPAP components)
- **Manufacturer:** Fisher & Paykel Healthcare
- **ASIN:** B08SCL9RYJ

WARRANTY AND SUPPORT

For warranty information or technical support regarding your Fisher & Paykel Brevida Nasal Pillows Mask Elbow, please refer to the documentation provided with your complete Brevida mask system or contact Fisher & Paykel Healthcare directly. You can typically find contact information on the manufacturer's official website.

Manufacturer Website: www.fphcare.com

© 2024 Fisher & Paykel Healthcare. All rights reserved. Information subject to change without notice.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B08SCL9RYJ>

More documents and user manuals: <https://manuals.plus/>