



Manuals.plus /

› VnioLife /

› VnioLife Ionic Foot Bath Detox Machine HK-802FS User Manual

VnioLife HK-802FS

VnioLife Ionic Foot Bath Detox Machine (Model HK-802FS) User Manual

Comprehensive Instructions for Safe and Effective Use

1. INTRODUCTION

This manual provides detailed instructions for the VnioLife Ionic Foot Bath Detox Machine, Model HK-802FS. Please read this manual thoroughly before operating the device to ensure proper use, safety, and optimal performance. Keep this manual for future reference.

2. SAFETY INFORMATION

WARNING: Failure to follow these safety instructions may result in injury or damage to the device.

- Use with caution around water. Ensure the main unit and power adapter do not come into contact with water.
- Keep out of reach of children. This device is not a toy.
- Avoid contact with skin and eyes when handling the array or water during operation.
- Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is for personal care use only.
- Always unplug the device from the power outlet before cleaning or when not in use.
- Do not operate the device if the power cord or plug is damaged.

3. PACKAGE CONTENTS

Verify that all components are present in the package:

- Foot Bath Basin
- Control Unit

- Detox Arrays (typically 2-3 units)
- Power Adapter
- Disposable Liners (100 count)
- Zipper Carrying Bag
- User Manual (this document)



Image 3.1: All components included in the VnioLife Ionic Foot Bath Detox Machine package, including the foot bath basin, control unit, detox arrays, power adapter, disposable liners, and carrying bag.

4. PRODUCT OVERVIEW

Familiarize yourself with the main parts of the device:

- **Control Unit:** Features an ON/OFF button and LED indicators for operation.
- **Detox Array:** The submersible component that facilitates the electrolysis process in the water.
- **DC Power Input:** Port on the control unit for connecting the power adapter.

- **Array Connection:** Port on the control unit for connecting the detox array.

DETOX FOOT SPA

**LED Working Indicator
Programmed
30 Min session**



ON/OFF



DC Power Input To Array



Image 4.1: Close-up view of the control unit, showing the ON/OFF button, LED indicators, DC power input, and array connection port. A foot bath with a liner and array is also shown.

5. SETUP

Follow these steps to set up your VnioLife Ionic Foot Bath Detox Machine:

1. Place the foot bath basin on a flat, stable surface.
2. Insert a disposable liner into the foot bath basin. This helps maintain hygiene and simplifies cleaning.
3. Fill the basin with warm water, ensuring it covers the feet up to the ankles. Do not overfill.
4. Add a small amount of non-iodized salt (e.g., sea salt or Himalayan salt) to the water. Approximately 1/2 teaspoon is a common starting point, but adjust based on desired conductivity. Stir to dissolve.
5. Place the detox array into the water, ensuring it is fully submerged.
6. Connect the array cable to the 'ARRAY' port on the control unit.
7. Connect the power adapter to the 'DC INPUT' port on the control unit, then plug the adapter into a standard AC 110V power outlet.

8. Carefully place your feet into the water.



✓ **Happiness is so simple**

✓ **Foot bath detoxification machine**

✓ **The choice of more sub-healthy population**

Image 5.1: A user's feet placed inside the foot bath basin, which is lined with a disposable plastic bag and contains the detox array submerged in water.

6. OPERATING INSTRUCTIONS

Once the setup is complete, you can begin a session:

1. Press the ON/OFF button on the control unit to start the device. The LED indicators will illuminate.
2. The device is programmed for an automatic 30-minute session. The timer cannot be adjusted during operation.
3. During the session, the detox array will initiate an electrolysis process, which causes the water to change color. This color change is a result of the interaction between the array, the water, and any dissolved salts or minerals.
4. After 30 minutes, the device will automatically shut off.
5. Carefully remove your feet from the basin.
6. Press the ON/OFF button again to ensure the unit is completely off, then unplug the power adapter from

the wall outlet.

Understanding Water Color Changes

The device operates on the principle of electrolysis, a chemical reaction that breaks down water into oxygen and hydrogen. These gases are attracted by the metal array parts, leading to an oxidation reaction. The final color of the water is a combination of this oxidation reaction and any substances present in the water. It is important to note that the water color will change even if feet are not present in the water, due to the interaction between the array and the water itself.




Image 6.1: Visual representation of water color changes over 5, 15, and 30 minutes during an ionic foot bath session, illustrating the progression of the electrolysis and oxidation process.



Image 6.2: Diagram illustrating the fixed 30-minute session duration of the VnioLife Foot SPA, with indicators for 5, 10, 15, 20, 25, and 30 minutes.

7. MAINTENANCE AND CARE

Proper maintenance ensures the longevity and hygiene of your device:

- **After Each Use:** Dispose of the used liner and water.
- **Cleaning the Basin:** Wipe the foot bath basin with a damp cloth and mild soap, then rinse thoroughly. Ensure no water enters the control unit.
- **Cleaning the Array:** The detox array should be cleaned after each use. Rinse it under running water to remove any residue. For stubborn buildup, a soft brush can be used. Ensure the array is completely dry

before storage.

- **Storage:** Store the device and its components in the provided zipper carrying bag in a cool, dry place, away from direct sunlight and moisture.



Image 7.1: The VnioLife Ionic Foot Bath Detox Machine with its customized blue carrying bag, demonstrating its portability and ease of storage.

8. TROUBLESHOOTING

If you encounter issues with your VnioLife Ionic Foot Bath Detox Machine, refer to the following common concerns:

- **Device does not power on:**
 - Ensure the power adapter is securely plugged into both the control unit and a working wall outlet.
 - Check if the wall outlet is functional by plugging in another appliance.
- **Water changes color but feet are not in the basin:**
 - This is normal. The water color change is primarily due to the electrolysis process between the detox array, the water, and any dissolved minerals or salts. It is not solely indicative of substances being drawn from the body.
- **No water color change:**
 - Ensure the detox array is fully submerged in water.
 - Verify that a sufficient amount of salt has been added to the water to increase conductivity.
 - Check that the array cable is securely connected to the control unit.
 - Inspect the detox array for visible damage or excessive buildup that might impede its function.

9. PRODUCT SPECIFICATIONS

Feature	Specification
Model Number	HK-802FS
Product Dimensions	22.52 x 13 x 10.25 inches
Item Weight	5 Pounds
Power Source	Corded Electric (Input: AC 110V, Output: DC 12V)
Material	Acrylonitrile Butadiene Styrene
Manufacturer	WL
First Available Date	January 15, 2021

10. WARRANTY AND SUPPORT

For warranty information, technical support, or any questions not covered in this manual, please contact the manufacturer directly. You can often find contact details on the product packaging or by visiting the official VnioLife store online.

Visit the VnioLife Store for more information.