

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

› [Bellababy](#) /

› [Bellababy USB Cable for Electric Breast Pump 8015 User Manual](#)

## Bellababy Cable

# Bellababy USB Cable for Electric Breast Pump 8015 User Manual

Model: Cable



Official Bellababy Product

---


## 1. INTRODUCTION

This user manual provides essential information for the safe and effective use of your Bellababy USB Cable, specifically designed for the Bellababy Electric Breast Pump model 8015. Please read this manual thoroughly before using the product and retain it for future reference.

The Bellababy USB Cable serves as a replacement or spare power source cable for your compatible electric breast pump, ensuring continuous operation and charging capabilities.

## 2. PRODUCT OVERVIEW

The Bellababy USB Cable is a dedicated charging and power supply cable for the Bellababy Electric Breast Pump 8015. It features a standard USB-A connector on one end for power input and a proprietary connector on the other end designed to fit the breast pump's power port.



**Image 2.1:** Bellababy USB Cable. This image displays the USB cable, showing both the standard USB-A connector and the specialized connector for the breast pump.

**Components Included:**

- 1 x Bellababy USB Charging Cable

### 3. SETUP INSTRUCTIONS

---

Follow these steps to properly connect the USB cable to your Bellababy Electric Breast Pump 8015:

1. **Identify Ports:** Locate the USB power input port on your Bellababy Electric Breast Pump 8015. This port is typically found on the side or back of the main pump unit.
2. **Connect to Pump:** Insert the specialized connector end of the Bellababy USB Cable firmly into the power input port of the breast pump. Ensure it is fully seated and secure.
3. **Connect to Power Source:** Insert the standard USB-A connector end of the cable into a compatible USB power source. This can be a USB wall adapter (not included), a power bank, or a computer USB port. Ensure the power source provides adequate power (typically 5V/1A or 5V/2A).
4. **Verify Connection:** Once connected, the breast pump should indicate that it is receiving power or

charging (e.g., a charging indicator light on the pump may illuminate).

**Note:** Always use a power source that meets the specifications of your Bellababy Electric Breast Pump 8015 to prevent damage to the device or the cable.

## 4. OPERATING INSTRUCTIONS

---

The Bellababy USB Cable functions as a power conduit for your electric breast pump. Its operation is straightforward:

- **Powering the Pump:** When connected to a suitable power source, the cable supplies electrical power to the breast pump, allowing it to operate even if its internal battery is depleted or removed.
- **Charging the Battery:** If your Bellababy Electric Breast Pump 8015 has a rechargeable internal battery, connecting the USB cable to a power source will initiate the charging process. The charging status will typically be indicated by the breast pump's display or indicator lights.

**Important:** Disconnect the cable from the power source and the breast pump when not in use or when the battery is fully charged to conserve energy and prolong battery life.

## 5. MAINTENANCE AND CARE

---

Proper care of your Bellababy USB Cable will ensure its longevity and reliable performance:

- **Cleaning:** Wipe the cable with a soft, dry cloth. Do not use abrasive cleaners, solvents, or immerse the cable in water. Ensure the connectors are free from dust and debris.
- **Storage:** Store the cable in a cool, dry place, away from direct sunlight and extreme temperatures. Avoid tightly coiling or bending the cable sharply, as this can damage the internal wires.
- **Handling:** Always grasp the plug ends when connecting or disconnecting the cable. Avoid pulling on the cable itself, as this can strain the connections and lead to damage.
- **Inspection:** Periodically inspect the cable for any signs of wear, fraying, cuts, or damage to the connectors. If any damage is observed, discontinue use immediately and replace the cable.

## 6. TROUBLESHOOTING

---

If you encounter issues while using your Bellababy USB Cable, refer to the following common troubleshooting steps:

- **Pump Not Powering On/Charging:**
  - Ensure both ends of the USB cable are securely connected to the breast pump and the power source.
  - Try a different USB power source (e.g., a different wall adapter or USB port) to rule out a faulty power source.
  - Check if the power source provides sufficient power (e.g., 5V/1A or 5V/2A).
  - Inspect the cable for visible damage (frays, cuts, bent pins in connectors). If damaged, replace the cable.
  - Ensure the breast pump itself is functioning correctly with its original cable or another known working cable.
- **Loose Connection:**
  - If the cable feels loose in the pump's port, gently try to reseat it. Avoid forcing it.
  - Check for any debris in the pump's port or the cable's connector that might prevent a secure

connection.

If the issue persists after following these steps, please contact Bellababy customer support for further assistance.

## 7. SPECIFICATIONS

Feature	Detail
Product Name	Bellababy USB Cable
Model Number	Cable (for Electric Breast Pump 8015)
Brand	Bellababy
Power Source	USB Cable
Operating Mode	Electric
Special Feature	Lightweight, Portable
Color	White
Dimensions (Package)	11.1 x 10.6 x 4.9 cm
Weight (Package)	30 grams
ASIN	B08XJWSGRT

## 8. WARRANTY AND SUPPORT

For warranty information, technical support, or any inquiries regarding your Bellababy USB Cable or Electric Breast Pump 8015, please refer to the documentation included with your original breast pump purchase or visit the official Bellababy website.

You may also contact Bellababy customer service directly for assistance.