

Engage Engage Rechargeable

User Manual for Engage Premium Rechargeable OTC Hearing Aids

Brand: Engage | Model: Engage Rechargeable

PRODUCT OVERVIEW

The Engage Premium Rechargeable OTC Hearing Aids are designed to provide clear, natural sound for adults with mild to moderate hearing loss. These devices offer advanced features including app-controlled sound settings, Bluetooth streaming, and a comfortable, discreet fit for all-day wear.

- **OTC Hearing Aids:** A rechargeable solution that delivers clear, natural sound to help you meet all your hearing needs. Wirelessly stream music from your phone, take calls, and listen to the world around you.
- **Features:** Award-winning OTC hearing aids with nano coating for moisture resistance, a speaker in the ear for greater clarity, directional microphones, a rechargeable battery, app-controlled sound settings & streaming from your Android device.
- **High-Definition Natural Sound:** Designed for adults with signs of mild to moderate hearing loss, Lucid Hearing technology helps produce more natural sound than traditional compression-based hearing aids, meaning clearer sound without distortion.
- **Lucid Hearing App:** Personalize your sound with the Lucid Hearing App & select from four pre-programmed hearing profiles. Adjust volume, control hearing aids together or independently, get support, and check the battery status.
- **In the Box:** Includes one pair of hearing aids, a rechargeable case, nine pairs of assorted ear tips, 10 wax stops, a cleaning tool, USB-C charging cable, power adapter, and a user manual.

SETUP

1. Remove the hearing aids from the charging case.
2. Insert into your ears using appropriate ear tips.
3. Turn on by pressing the button on each device.
4. Adjust volume via the app or directly on the aid.
5. To connect to your phone for streaming, enable Bluetooth and select 'Engage' in your phone's settings.



Figure 1: Engage Hearing Aids and Charging Case. The aids are ready to wear directly from the case.

HIGH DEFINITION SOUND

- Customize your hearing experience with the Lucid Hearing App® & LucidShape®
- Explore and select your preferred hearing program
- Monitor your battery life

Bluetooth
streaming from
your device




Figure 2: Discreet Fit of the Engage Hearing Aid. Designed for comfortable all-day wear.

OPERATING INSTRUCTIONS

App Control and Personalization

The Lucid Hearing App allows you to fine-tune your listening experience. You can select from four pre-programmed hearing profiles or customize your own settings. The app also provides volume control for both hearing aids, either together or independently, and allows you to monitor battery status.

DISCREET, ALL DAY COMFORT



Available in 3 colors

Figure 3: Lucid Hearing App Interface. Customize your sound settings and monitor your devices.

Bluetooth Streaming

Engage hearing aids support Bluetooth streaming from your Android device. This feature allows you to wirelessly stream music, take phone calls, and listen to other audio directly through your hearing aids, enhancing your daily interactions.

Your browser does not support the video tag.

Video 1: Engage Hearing Aids Overview. This video demonstrates the key features and benefits of the Engage Bluetooth Hearing Aids, including their rechargeable nature, customizability via the app, and discreet fit. It also highlights their suitability for mild to moderate hearing loss and the ability to stream calls and music.

MAINTENANCE

Regular maintenance ensures optimal performance and longevity of your Engage hearing aids. Use the provided cleaning tool to remove any earwax or debris. Store the aids in their charging case when not in use. The aids feature a nano coating for moisture resistance, but it is important to keep them dry and clean.

- **Cleaning:** Use the provided cleaning tool to remove earwax and debris.
- **Storage:** Always store the hearing aids in their rechargeable case when not in use.
- **Charging:** Charge the aids using the included USB-C cable and power adapter.
- **Wax Stops:** Replace wax stops as needed to prevent blockages and maintain sound quality.

WHAT'S INCLUDED

Includes everything needed to get started!

Images not to scale

Rechargeable batteries *allow up to 16 hours of listening*

Figure 4: What's Included. This image displays all components provided with the Engage hearing aids, emphasizing the cleaning tool and wax stops for maintenance.

SPECIFICATIONS

Feature	Detail
Item Weight	10.5 ounces
Package Dimensions	7.52 x 4.57 x 1.93 inches

Item Model Number	Engage Rechargeable
Brand	Engage
Color	Grey
Material	Plastic
Power Source	Battery Powered
First Available Date	July 28, 2021




WARRANTY AND SUPPORT




Information regarding specific warranty details and direct support contact for the Engage Premium Rechargeable OTC Hearing Aids is not available in this manual. Please refer to the packaging or the manufacturer's official website for the most up-to-date warranty information and customer support contacts.



© 2023 Engage. All rights reserved.

Related Documents - Engage Rechargeable

	<p>Lucid Hearing TALA™ Premium Rechargeable OTC Hearing Aids User Manual</p> <p>Comprehensive user manual for Lucid Hearing's TALA™ Premium Rechargeable OTC Hearing Aids. Learn about setup, features, app usage, care, troubleshooting, and technical specifications.</p>
	<p>Lucid fio™ Premium Rechargeable OTC Hearing Aids User Manual</p> <p>This user manual provides comprehensive instructions for the Lucid fio™ premium rechargeable OTC hearing aids, covering setup, operation, care, troubleshooting, and important safety information.</p>
	<p>TTQ OTC Hearing Aids for Seniors: Smart Noise Cancelling & Bluetooth</p> <p>User manual for TTQ OTC Hearing Aids, featuring smart noise cancellation, Bluetooth connectivity, 3-level volume adjustment, and a rechargeable, fashionable design for seniors with hearing loss.</p>

 <p>ENHANCE SELECT 50R User guide Rechargeable Receiver-in-the-Ear Hearing Aids Let's hear what you can hear.™</p>	<p>Jabra Enhance Select 50R User Guide: Rechargeable Hearing Aids</p> <p>Comprehensive user guide for Jabra Enhance Select 50R rechargeable receiver-in-the-ear hearing aids. Learn about setup, usage, safety warnings, maintenance, and troubleshooting.</p>
 <p>Love at first sound Phonak Audio Marvel Your Guide to Marketing</p>	<p>Phonak Audio Marvel: Your Comprehensive Marketing Guide</p> <p>This marketing toolkit provides hearing care professionals with a comprehensive guide to promoting Phonak Audio Marvel hearing aids. It details product features, target audiences, campaign strategies, and available marketing assets to drive business growth and client satisfaction.</p>
 <p>2021 HEARING AID Purchasing Guide What you need to know before making a purchase</p>	<p>Davidson Hearing Aid Centres: 2021 Purchasing Guide for Better Hearing</p> <p>A comprehensive guide from Davidson Hearing Aid Centres detailing the journey to better hearing, types of hearing loss, the impact of hearing loss on health, modern hearing aid technology, available styles, costs, and essential accessories. Includes product spotlights and contact information.</p>