



[Manuals.plus](#) /

> [ALLJOY](#) /

> ALLJOY 24"x33" Weighted Heating Pad Instruction Manual

ALLJOY Heat

ALLJOY 24"x33" Weighted Heating Pad Instruction Manual

Brand: ALLJOY | Model: Heat

1. INTRODUCTION

This manual provides essential information for the safe and effective use of your ALLJOY 24"x33" Weighted Heating Pad. Designed to deliver soothing heat therapy, this extra-large heating pad offers targeted relief for the neck, shoulders, and back, and can also be applied to other areas such as the waist, abdomen, and legs. Please read all instructions carefully before use and retain this manual for future reference.

2. IMPORTANT SAFETY INFORMATION

- Avoid use on broken or damaged skin.
- Do not use while sleeping.
- Keep away from flammable materials.
- Utilize the automatic shut-off feature to prevent excessive heating.
- Do not use pins or other metallic objects to fasten the heating pad.
- Never leave the heating pad unattended, especially when children or individuals with reduced physical, sensory, or mental capabilities are present.
- Do not fold, sit on, or crush the heating pad during operation.

3. PACKAGE CONTENTS

Upon unpacking, please ensure all items are present and in good condition:

- ALLJOY 24"x33" Weighted Heating Pad
- Detachable Controller with Power Cord
- Instruction Manual



Image: The ALLJOY 24"x33" Weighted Heating Pad, its detachable controller, and product packaging.

4. SETUP INSTRUCTIONS

1. **Unpack:** Carefully remove the heating pad and all accessories from the packaging.
2. **Connect Controller:** Locate the connector on the heating pad and firmly plug in the detachable controller. Ensure a secure connection.
3. **Power Connection:** Plug the power cord from the controller into a standard electrical outlet.
4. **Positioning:** Place the heating pad on the desired area of your body (neck, shoulders, back, etc.). Use the snap fasteners at the neck and the tie straps around the waist to secure the pad for hands-free use. The weighted edges help keep the pad in place.

Wrap Design & Shoulder Weighted Edges

Fits naturally and comfortably around the neck and shoulders for maximum warmth



Neck Closure
Snap Fasteners



Weighted Design




Image: A user securing the heating pad with snap fasteners around the neck and tie straps around the waist for a snug fit.

5. OPERATING INSTRUCTIONS

The ALLJOY heating pad features an easy-to-use controller for personalized heat therapy.

1. **Power On:** Press the power button on the controller to turn on the heating pad. The display will illuminate.
2. **Select Heat Setting:** Use the '+' and '-' buttons to adjust the heat level. The pad offers 6 heat settings, ranging from 100°F to 150°F. Choose a comfortable temperature for your needs.
3. **Set Timer:** Press the timer button to cycle through the available auto shut-off times: 30, 60, or 90 minutes. This feature ensures safety and prevents excessive heating.
4. **Enjoy Heat Therapy:** Allow the pad to heat up. It is designed for fast heating, penetrating muscles within seconds.
5. **Power Off:** To turn off the heating pad manually, press the power button. The pad will also automatically shut off once the set timer expires.



Image: The digital controller showing temperature and timer settings for customized heat therapy.

Official Product Video: Heating Pad for Back Pain Relief

Your browser does not support the video tag.

Video: An official ALLJOY video demonstrating the features and use of the heating pad for back pain relief, including how to wear it and adjust settings.

6. CARE AND MAINTENANCE

To ensure the longevity and hygiene of your heating pad, follow these cleaning instructions:

- **Disconnect Power:** Always unplug the heating pad from the electrical outlet before cleaning.
- **Remove Controller:** Detach the controller from the heating pad. The pad is machine washable, but the controller is not.
- **Machine Wash:** Machine wash the heating pad on a gentle cycle with cold water. Do not bleach or iron.
- **Tumble Dry:** Tumble dry on low heat or air dry. Ensure the pad is completely dry before reattaching the controller and using it again.

High Quality Material




Image: The heating pad being placed into a washing machine, highlighting its machine-washable feature after removing the controller.

7. TROUBLESHOOTING

If you encounter any issues with your ALLJOY heating pad, please refer to the following common solutions:

- **Pad Not Heating:** Ensure the power cord is securely plugged into both the heating pad and the electrical outlet. Check if the outlet is functional.
- **Inconsistent Heat:** Verify that the heating pad is laid flat and not folded or bunched up during use, as this can affect heat distribution.
- **Controller Display Not Working:** Unplug the unit, wait a few minutes, and then plug it back in. If the issue persists, contact customer support.

8. PRODUCT SPECIFICATIONS

Feature	Detail
Product Dimensions	24 x 1 x 33 inches
Item Weight	2.7 Pounds
Material	Microplush Flannel Fabric
Heat Settings	6 settings (100°F-150°F)
Timer Settings	30, 60, 90 minutes (Auto Shut-Off)
Special Features	Machine Washable, Automatic Shut Off, Weighted Edges, Snap Fasteners, Fixture Straps
Safety Certification	ETL-approved

9. WARRANTY AND CUSTOMER SUPPORT

The ALLJOY heating pad is ETL-approved for safe use. For any issues or inquiries regarding your product, please feel free to contact ALLJOY customer support. We are committed to providing a satisfactory solution.

Additional protection plans may be available for purchase to extend coverage beyond the standard warranty period. Please refer to your purchase documentation for details on available protection plans.

Fast Heating, Heating Evenly, Warm Enough



Machine Washable



Warms Up in Seconds



Timing Setting



6/10 Heat Setting




Image: The heating pad in use, highlighting its fast heating capability and customizable timing settings.

Effective Full Body Pain Relief

Neck, Shoulder, Back, Abdomen, Leg, Waist



Image: Various applications of the heating pad for effective full-body pain relief, including neck, shoulder, back, abdomen, and legs.



**Skin-Friendly
Flannel**

Image: An illustration detailing the high-quality, multi-layer material construction of the heating pad.

Official Product Video: Extra Large Heating Pad for Back Pain and Cramps Relief

Your browser does not support the video tag.

Video: An official ALLJOY video showcasing the extra-large size and various ways to use the heating pad for comprehensive pain and cramp relief.