

toberto HP4076

toberto Electric Heating Pad HP4076 User Manual

Model: HP4076

INTRODUCTION

The toberto Electric Heating Pad is designed to provide targeted heat therapy for pain and cramp relief. Made from a high-quality, ultra-soft micro-plush material, it offers comfort and is gentle on the skin. This heating pad features 6 adjustable heating settings and 4 timer modes, along with an automatic shut-off function for safety and convenience. Its generous 16" x 30" size makes it suitable for various body areas, including the neck, back, shoulders, abdomen, and legs.



The tobeto Electric Heating Pad provides soothing heat for various body parts.

Larger Size for Back Pain Relief



Effective heat therapy for abdominal pain, knee pain, shoulder pain, and backache.

6

Heating Modes

4

Timer Modes



0.5-1-1.5-2H Auto Shut Off

The 16" x 30" size offers ample coverage for larger areas.

IMPORTANT SAFETY INFORMATION

Please read all instructions and warnings before using this heating pad to ensure safe operation and prevent injury or damage to the product.

- Do not sit, stand, or lie directly on the heating pad, as this can increase the risk of scalding and damage the internal components.
- Ensure the warm-up time and temperature settings are within a healthy and controllable range for your comfort. Prolonged use at high temperatures may cause burns, especially for individuals with reduced skin sensitivity.
- Avoid using the heating pad in an unforeseen or unsupervised state.
- If you are using heat therapy for a medical condition, consult your doctor before use.
- The heating pad automatically turns off after a set period (up to 2 hours) for safety considerations. This is a protective feature.

Super Soft Plush Fabric



Internal NTC heating wire and PTC alloy wire provide fast heating and prevent overheating.

SETUP


1. Unpack the heating pad and controller from its packaging.
2. Connect the power cord securely to the heating pad's input port. Ensure the connection is firm.
3. Plug the power cord into a standard electrical outlet. The 7-foot power cord provides ample reach.

OPERATING INSTRUCTIONS

The heating pad features an LCD controller for easy operation.

Fast and Safe

Built in NTC heating wire provides fast heating,
also prevents overheating



The controller allows selection of 6 heating modes and 4 timer modes.

Power On/Off

- Press the power button on the controller to turn the heating pad on or off.

Adjusting Heating Settings

The heating pad offers 6 temperature settings, heating up quickly in approximately 30 seconds:

- **Setting 1:** Approximately 98°F (36.7°C)
- **Setting 2:** Approximately 104°F (40°C)
- **Setting 3:** Approximately 115°F (46.1°C)
- **Setting 4:** Approximately 125°F (51.7°C)
- **Setting 5:** Approximately 132°F (55.6°C)
- **Setting 6:** Approximately 140°F (60°C)

Use the temperature control buttons on the LCD controller to cycle through these settings.

Adjusting Timer Modes

The heating pad includes 4 timer modes with automatic shut-off:

- 0.5 hours (30 minutes)
- 1 hour
- 1.5 hours
- 2 hours

Use the timer control buttons on the LCD controller to select your desired duration. The pad will automatically turn off after the selected time.


Wet Heat Therapy

For wet heat therapy, lightly spray water on one side of the heating pad. This can help to loosen muscles and prepare them for exercise. Do not immerse the pad in water or use it when excessively wet.

MAINTENANCE

The tobeto Electric Heating Pad is designed for easy cleaning and long-term use.

- **Machine Washable:** The micro-plush fabric is machine washable. Disconnect the controller and power cord completely before washing.
- Wash on a gentle cycle with cold water.
- Do not bleach or iron.
- Tumble dry on low heat or air dry. Ensure the heating pad is completely dry before reconnecting the controller and using it again.
- Proper care helps maintain the softness and comfort of the fabric for years.



The super soft plush fabric is skin-friendly and easy to maintain.

TROUBLESHOOTING

- **Controller Blinking:** If the controller continues to blink, it indicates that the heating pad is not firmly connected to the power cord or the power outlet. Ensure all connections are secure.
- **Automatic Shut-Off:** If the heating pad automatically turns off during use, it may be due to the overheat protection being activated or the timer reaching its set duration. Press the power button again to turn it back on.
- If you experience any other problems or concerns with your heating pad, please refer to the Warranty and Support section for assistance.

SPECIFICATIONS

Feature	Detail
Model Name	HP4076

Dimensions	16 inches x 30 inches
Material	Microplush
Wattage	75 watts
Heating Settings	6 levels (approx. 98°F to 140°F)
Timer Modes	4 levels (0.5, 1, 1.5, 2 hours)
Auto Shut-Off	Yes
Power Cord Length	7 feet
Special Features	Breathable, Lightweight, Skin Friendly

WARRANTY AND SUPPORT

Your toberto Electric Heating Pad comes with the following assurances:

- **60-Day Return:** We offer a free 60-day return policy from the date of purchase.
- **365-Day Replacement:** We provide a free 365-day replacement for any manufacturing defects or issues.

If you have any problems or concerns regarding your heating pad, please do not hesitate to contact us. Our customer support team is available to assist you.