

- Manuals.plus /
- › [PEPE Mobility](#) /
- › PEPE Mobility Adjustable Step Stool Model 2412 User Manual

PEPE Mobility 2412

PEPE Mobility Adjustable Step Stool Model 2412 User Manual



Brand: PEPE Mobility | Model: 2412

1. INTRODUCTION AND PRODUCT OVERVIEW

The PEPE Mobility Adjustable Step Stool is designed to provide a stable and secure elevated surface, assisting users with accessing bathtubs, high shelves, or vehicles like caravans. Its robust aluminum construction ensures durability and corrosion resistance, while the non-slip surface and rubber tips enhance user safety. This step stool features tool-free assembly and adjustable height settings to accommodate various needs.



Image 1.1: Overview of the PEPE Mobility Adjustable Step Stool, highlighting its 395 lbs (180 kg) weight capacity and adjustable height range of 8" to 10" (20.5 - 25.5 cm). The diagram also indicates tool-free assembly.



20,5 - 25,5 cm
(8" - 10")

Image 1.2: The PEPE Mobility Adjustable Step Stool fully assembled, showcasing its black non-slip platform and sturdy legs.

2. SAFETY INFORMATION

Please read and understand all safety instructions before using the PEPE Mobility Adjustable Step Stool. Failure to follow these instructions may result in injury.

- **Maximum Weight Capacity:** Do not exceed the maximum user weight of 395 lbs (180 kg).
- **Stable Surface:** Always place the step stool on a flat, level, and stable surface. Ensure all four legs are firmly on the ground before use.
- **Non-Slip Surface:** The step stool features a non-slip surface. However, exercise caution when stepping on or off, especially if the surface is wet.
- **Leg Adjustment:** Ensure all leg adjustments are securely locked in place before use. Uneven leg heights can cause instability.
- **Intended Use:** This product is designed as a step stool for assistance with elevation. Do not use it for purposes other than its intended use.
- **Children:** Supervise children when they are near or using the step stool.
- **Inspection:** Regularly inspect the step stool for any signs of damage, wear, or loose components. Do not use if damaged.

3. PACKAGE CONTENTS

Verify that all components are present and undamaged upon opening the package.

- 1x Step Stool Platform
- 4x Leg Assemblies (each with adjustable height extension and rubber tip)
- 4x Securing Knobs with Washers

4. SETUP AND ASSEMBLY

The PEPE Mobility Adjustable Step Stool is designed for tool-free assembly. Follow these steps to assemble your step stool:

1. **Attach Leg Assemblies:** Insert each leg assembly into the designated slots on the underside of the step stool platform. Ensure the metal pins on the leg assemblies align with the holes in the platform.
2. **Secure Legs:** From the underside of the platform, insert a securing knob with its washer through the hole and into the threaded opening of each leg assembly. Turn clockwise to tighten until snug. Do not overtighten.
3. **Adjust Height (Initial):** Before final use, adjust each leg to your desired height. Press the brass pin on each leg extension and slide the inner tube up or down until the pin locks into one of the three adjustment holes. Ensure all four legs are set to the same height for stability.
4. **Final Check:** Once assembled and adjusted, gently press down on the step stool to ensure it is stable and all components are securely fastened.



Video 4.1: This video demonstrates the tool-free assembly process of the PEPE Mobility Adjustable Step Stool, showing how to attach and secure the legs to the platform.



Image 4.2: An underside view of the step stool, illustrating the attachment points for the leg assemblies.



Image 4.3: A close-up view showing how to secure the leg assemblies to the platform using the provided knobs and washers.

5. OPERATING INSTRUCTIONS

To use the PEPE Mobility Adjustable Step Stool safely and effectively:

1. **Placement:** Position the step stool directly in front of the object or area you wish to access (e.g., bathtub, cabinet, vehicle entrance).
2. **Stability Check:** Before stepping on, ensure the step stool is stable and does not wobble. If it wobbles, re-check leg

adjustments and placement.

3. **Step On:** Carefully step onto the center of the step stool's non-slip surface. Use any available support (e.g., grab bar, wall) for balance if needed.
4. **Step Off:** When descending, step off the stool carefully, maintaining balance.



Video 5.1: This video demonstrates the practical use of the PEPE Mobility Adjustable Step Stool for safely entering and exiting a bathtub.



Image 5.2: A user demonstrating the use of the step stool to access a caravan, highlighting its versatility for various elevated access needs.

6. ADJUSTING HEIGHT

The step stool offers three height adjustment positions, ranging from 8 inches (20.5 cm) to 10 inches (25.5 cm).

1. **Unlock Leg:** Turn the step stool upside down or tilt it to access the legs. Locate the brass pin on each adjustable leg extension.
2. **Slide to Desired Height:** Press the brass pin inward and slide the inner leg tube up or down to the desired height. The pin will snap into place when it aligns with one of the adjustment holes.
3. **Ensure Evenness:** Repeat this process for all four legs, ensuring they are all set to the exact same height to maintain stability.
4. **Verify Lock:** Confirm that all four brass pins are fully engaged in their respective holes before using the step stool.



Image 6.1: A detailed view of the leg adjustment mechanism, showing how to press the pin to change the height setting.

7. MAINTENANCE AND CARE

- **Cleaning:** Clean the step stool with a damp cloth and mild soap. Avoid abrasive cleaners or solvents that could damage the material.
- **Drying:** Ensure the step stool is completely dry before storage or next use, especially the non-slip surface, to prevent slips.
- **Inspection:** Periodically check all screws, knobs, and leg pins for tightness and proper function. Inspect the rubber tips for wear and tear. Replace worn parts if necessary.
- **Storage:** Store the step stool in a dry place away from direct sunlight and extreme temperatures.

8. TROUBLESHOOTING

Problem	Possible Cause	Solution
Step stool wobbles or is unstable.	◦ Legs not set to the same height. ◦ Securing knobs are loose. ◦ Placed on an uneven surface.	◦ Ensure all four legs are adjusted to the same height and pins are locked. ◦ Tighten all securing knobs. ◦ Move the step stool to a flat, level surface.
Leg adjustment pins do not lock.	◦ Pin is obstructed or misaligned.	◦ Ensure the inner tube is correctly aligned with the adjustment holes. Wiggle slightly while pressing the pin.

9. SPECIFICATIONS

Brand	PEPE Mobility
Model Number	2412
Material	Aluminum
Color	Black
Product Dimensions (L x W x H)	18.7" x 12.4" x 8" - 10" (Adjustable)
Weight Capacity	395 lbs (180 kg)

Item Weight	1.89 kg
Assembly Required	Yes (Tool-Free)
Special Features	Adjustable Height, Non-Slip Surface, Rubber Tips



Image 9.1: Dimensional diagram of the step stool, illustrating its width, depth, and maximum height.

10. WARRANTY AND SUPPORT

For warranty information, technical support, or to inquire about replacement parts, please contact PEPE Mobility customer service. Refer to your purchase documentation or the retailer's website for specific contact details.

Manufacturer: IMD SL

Brand: PEPE Mobility

[Visit the PEPE Mobility Store on Amazon](#)

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B09DT9ZZW7>

More documents and user manuals: **<https://manuals.plus/>**