

[Manuals.plus](#) /

› [Starkey Products](#) /

› Starkey Hearing Aid Domes User Manual

Starkey Products B09H2NB1KH

Starkey Hearing Aid Domes User Manual

Comfort Ear Buds Open/Vented Variety Pack

1. INTRODUCTION

This manual provides essential information for the proper use and care of your Starkey Hearing Aid Domes Comfort Ear Buds Open/Vented Variety Pack. This pack contains 32 domes of various sizes (5mm, 6mm, 7mm, 9mm) to help you find the most comfortable and effective fit for your hearing aids. Proper selection and maintenance of domes are crucial for optimal sound quality and comfort.



Image 1.1: The Starkey Hearing Aid Domes Variety Pack (Kit 53264-000) in its closed packaging.

2. SETUP AND SELECTION

The variety pack is designed to assist you in determining the ideal dome size for your individual ear canal and hearing aid. Selecting the correct size is vital for comfort, sound delivery, and preventing feedback.

2.1 Identifying Dome Sizes

The variety pack includes open/vented domes in 5mm, 6mm, 7mm, and 9mm sizes. Each size is typically separated within the packaging for easy identification.



Image 2.1: Contents of the Starkey Hearing Aid Domes Variety Pack, showing different sizes.

2.2 Attaching the Dome

1. **Choose a Size:** Start with a size that feels comfortable and provides good sound. If you experience feedback (whistling) or the dome feels loose, try a larger size. If it feels too tight or causes discomfort, try a smaller size.
2. **Remove Old Dome (if applicable):** Gently pull the old dome off the receiver or thin tube of your hearing aid.
3. **Attach New Dome:** Align the opening of the new dome with the end of the hearing aid's receiver or thin tube. Push firmly until the dome is securely seated. Ensure it is fully attached and not loose.

3. OPERATING (WEARING)

Once the correct dome is attached, insert the hearing aid into your ear canal. The dome should sit comfortably and securely within the ear, creating a seal that allows sound to be delivered effectively without discomfort or excessive feedback.

- **Gentle Insertion:** Gently push the hearing aid into your ear canal until the dome is comfortably seated. Do not force it.
- **Check for Comfort:** The dome should feel snug but not painful. If there is discomfort, try adjusting the position or consider a different size.
- **Sound Quality:** Listen for clear sound and minimal feedback. An improperly seated or sized dome can lead to poor sound quality or whistling.

4. MAINTENANCE

Regular cleaning and timely replacement of your hearing aid domes are essential for hygiene, performance, and longevity of your hearing aids.

4.1 Cleaning

- **Daily Wipe:** After each use, gently wipe the dome with a clean, dry cloth to remove earwax and debris.
- **Avoid Liquids:** Do not use water, alcohol, or cleaning solutions directly on the domes while they are attached to the hearing aid, as this can damage the hearing aid components. If a more thorough cleaning is needed, remove the dome and clean it separately with a mild, damp cloth, then ensure it is completely dry before reattaching.

4.2 Replacement

Domes should be replaced regularly, typically every 2-4 weeks, or sooner if they become stiff, discolored, torn, or if sound quality degrades. Frequent replacement ensures optimal hygiene and acoustic performance.

5. TROUBLESHOOTING

Problem	Possible Cause	Solution
Hearing Aid Whistles (Feedback)	Dome is too small or not properly seated; earwax buildup.	Try a larger dome size. Ensure the dome is fully inserted into the ear canal and securely attached to the hearing aid. Clean the dome and ear canal.
Discomfort or Irritation	Dome is too large or improperly inserted.	Try a smaller dome size. Re-insert the hearing aid gently. If discomfort persists, consult an audiologist.
Poor Sound Quality / Muffled Sound	Dome is clogged with earwax or damaged.	Clean the dome thoroughly. If cleaning does not resolve the issue, replace the dome.
Dome Falls Off Easily	Dome is not securely attached or is worn out.	Ensure the dome is pushed firmly onto the receiver/tube. Replace the dome if it appears stretched or damaged.

6. SPECIFICATIONS

- **Product Name:** Starkey Hearing Aid Domes Comfort Ear Buds Open/Vented Variety Pack
- **Model (ASIN):** B09H2NB1KH
- **Kit Number:** 53264-000
- **Dome Types:** Open/Vented
- **Sizes Included:** 5mm, 6mm, 7mm, 9mm
- **Quantity:** 32 Domes (8 of each size)
- **Package Dimensions:** 2.99 x 2.2 x 0.67 inches
- **Package Weight:** 0.63 ounces
- **Manufacturer:** Starkey
- **First Available:** September 24, 2021

7. WARRANTY AND SUPPORT

As these are consumable accessory items, a specific manufacturer's warranty beyond initial defects may not apply. For any issues related to the fit, performance, or compatibility of these domes with your Starkey hearing aids, it is

recommended to consult with your audiologist or the hearing aid professional who fitted your device.

For general product inquiries or support, you may visit the[Starkey Products Store on Amazon](#) or refer to the documentation provided with your specific hearing aid device.

Legal Disclaimer: *Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. This product is an accessory for hearing aids and not a medical device itself.*

© 2024 Starkey Products. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question