

MEDISANA FS 300

MEDISANA FS 300 Foot Spa Instruction Manual

Model: FS 300

1. INTRODUCTION

Thank you for choosing the MEDISANA FS 300 Foot Spa. This device is designed to provide a relaxing foot massage experience, combining bubble massage, vibration massage, and a heat retention function. Please read this instruction manual carefully before initial use to ensure safe operation and optimal performance.

2. SAFETY INSTRUCTIONS

WARNING: Failure to follow these safety instructions may result in electric shock, fire, or serious injury.

- Always read the entire user manual before operating the device.
- Do not use the product if you suffer from a medical condition that may be aggravated by massage. Consult a physician if you have any concerns or pre-existing medical conditions, especially if you have diabetes, circulatory problems, open wounds, or skin infections.
- Keep the device away from water. Never immerse the unit in water or other liquids. Do not use near bathtubs, showers, or sinks.
- Ensure the device is plugged into a safe, grounded power outlet. Do not use extension cords unless absolutely necessary and ensure they are rated for the device's power requirements.
- Unplug the device from the power outlet immediately after use, before filling or emptying water, and before cleaning.
- Do not operate the device if the cord or plug is damaged, if it is not working properly, or if it has been dropped or damaged. Contact customer support for repair.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Do not stand in the foot spa. Use only while seated.
- Do not fill the foot spa above the maximum fill line.
- Adjust temperature and massage settings according to your personal preferences and needs. Discontinue use if you experience any discomfort.

3. PRODUCT OVERVIEW

The MEDISANA FS 300 Foot Spa is designed for ease of use and effective foot relaxation. It features a robust housing, integrated massage rollers, and a simple control panel.



Figure 1: MEDISANA FS 300 Foot Spa, top-down view with water and massage rollers.

Intensive foot reflexology massage




Figure 2: Close-up of the MEDISANA FS 300 Foot Spa highlighting key features.

Key Components:

- **Control Panel:** Located at the front, with a power button and function indicator.
- **Massage Rollers:** Two sets of removable rollers for intensive foot reflexology massage.
- **Bubble Jets:** Integrated jets for a soothing bubbling massage.
- **Vibration Function:** Provides a gentle vibrating massage.
- **Heat Retention Function:** Helps maintain the water temperature.
- **Water Fill Line:** Indicates the maximum safe water level.

4. SETUP

1. **Placement:** Place the foot spa on a stable, level, and non-slip surface. Ensure it is near a power outlet and away from other electrical appliances.
2. **Filling with Water:** Fill the foot spa with warm water up to the indicated maximum fill line. While the device has a heat retention function, it does not heat cold water. For best results, use pre-warmed water.
3. **Power Connection:** Plug the power cord into a suitable electrical outlet.

5. OPERATING INSTRUCTIONS

The MEDISANA FS 300 Foot Spa offers a 3-in-1 function for a comprehensive foot massage experience.



Figure 3: Control panel with function settings.

1. **Power On:** Press the power button on the control panel. The indicator light will illuminate.
2. **Function Selection:** The device typically operates in a combined mode: **1 = Bubble + Massage + Heat**. This activates the bubbling massage, vibration massage, and heat retention functions simultaneously.
3. **Using the Foot Spa:** Once the functions are active, carefully place your feet into the water. You can gently move your feet over the massage rollers for an enhanced reflexology massage.
4. **Using Without Water:** The vibration massage function can also be used without water for a gentle, dry foot massage.
5. **Duration:** Use the foot spa for a comfortable duration, typically 10-20 minutes.
6. **Power Off:** To turn off the device, press the power button again (**0 = OFF**). The indicator light will turn off.



Figure 4: User enjoying the foot spa with feet immersed.

6. MAINTENANCE AND CLEANING

Regular cleaning ensures hygiene and extends the lifespan of your foot spa.

1. **Unplug:** Always unplug the device from the power outlet before cleaning.
2. **Empty Water:** Carefully empty the water from the foot spa.
3. **Clean Interior:** Wipe the interior and exterior surfaces with a soft, damp cloth. A mild, non-abrasive cleaner can be used if necessary.
4. **Do Not Immerse:** Never immerse the entire unit in water or place it under running water.
5. **Storage:** Store the foot spa in a cool, dry place, away from direct sunlight and out of reach of children.

7. TROUBLESHOOTING

If you encounter any issues with your MEDISANA FS 300 Foot Spa, please refer to the following common solutions:

- **Device does not turn on:**

- Ensure the power cord is securely plugged into a working electrical outlet.
- Check if the power button is pressed to the 'ON' position.

- **No bubbles or vibration:**

- Verify that the device is powered on and the function is selected (Mode 1).
- Ensure the water level is not too high, which might impede bubble generation.

- **Water does not stay warm:**

- The device has a heat *retention* function, not a heating function. Always fill with pre-warmed water.
- Ensure the heat function is active (part of Mode 1).

If the problem persists after attempting these solutions, please contact MEDISANA customer support.

8. SPECIFICATIONS

Technical details for the MEDISANA FS 300 Foot Spa:




Figure 5: Dimensions of the MEDISANA FS 300 Foot Spa.

Specification	Detail
Brand	Medisana
Model Number	FS 300
Product Dimensions (L x W x H)	34 x 42 x 16 cm
Item Weight	1.24 Grams
Material	Plastic
Color	White
Power Source	Electricity
Special Feature	Heat (Retention)

Massager Form	Massage Roller
Suitable for Shoe Size	Up to 46 (EU)
Country of Origin	China

9. WARRANTY AND SUPPORT

For detailed warranty information, please refer to the warranty card included with your product packaging or visit the official MEDISANA website. In case of technical issues, questions, or to request service, please contact MEDISANA customer support directly. Keep your proof of purchase for warranty claims.