

[Manuals.plus](#) /

› [UPBRIGHT](#) /

› UpBright 5V AC/DC Adapter and Micro USB Charging Cable User Manual for Belifu TENS001, SM9126, AS8015, AS8012B Massagers

UPBRIGHT Belifu TENS001 SM9126 USB Cable Power

UpBright 5V AC/DC Adapter and Micro USB Charging Cable User Manual

Compatible with Belifu TENS001, SM9126, AS8015, AS8012B Electronic Pulse Massagers

1. INTRODUCTION

This manual provides essential information for the safe and effective use of your UpBright 5V AC/DC Adapter and Micro USB Charging Cable. This product is designed to provide power and charging capabilities for compatible devices, specifically Belifu TENS001, SM9126, AS8015, and AS8012B Dual Channel Electronic Pulse Massagers.

1.1 What's Included

- 1 x UpBright 5V AC/DC Adapter
- 1 x Micro USB Charging Cable



Figure 1.1: The UpBright 5V AC/DC Adapter and Micro USB Charging Cable. The image displays a white USB wall adapter and a black Micro USB cable, which are the primary components of this product.

2. SETUP INSTRUCTIONS

Follow these steps to properly connect and use your UpBright power supply.

1. **Connect the Micro USB Cable:** Insert the standard USB-A end of the Micro USB charging cable into the USB port of the UpBright 5V AC/DC Adapter. Ensure a secure connection.
2. **Connect to Device:** Insert the Micro USB end of the cable into the charging port of your compatible device (e.g., Belifu TENS001 massager).
3. **Plug into Power Outlet:** Plug the AC/DC Adapter into a standard wall power outlet (100-240VAC, 50/60Hz).
4. **Verify Connection:** Confirm that your device indicates it is receiving power or charging, typically shown by an indicator light or on-screen display.



Figure 2.1: The Micro USB Charging Cable. This cable connects the power adapter to your compatible electronic device for charging or power supply.

3. OPERATING INSTRUCTIONS

Once connected, the UpBright 5V AC/DC Adapter and Micro USB Charging Cable will provide continuous power or charge the internal battery of your compatible device.

- The adapter is designed for a wide input voltage range (100-240VAC), making it suitable for international travel with appropriate plug adapters (not included).
- The output is 5VDC, 500mA, which is standard for many small electronic devices, including the specified Belifu massagers.
- Always ensure the adapter and cable are fully inserted into their respective ports to prevent intermittent power supply.

4. MAINTENANCE

Proper care will extend the life of your power adapter and cable.

- **Cleaning:** Disconnect the adapter from the power outlet and device before cleaning. Use a dry, soft cloth to wipe the surfaces. Do not use liquid cleaners or solvents.
- **Storage:** Store the adapter and cable in a cool, dry place when not in use. Avoid extreme temperatures.
- **Handling:** Do not bend or crimp the cable excessively, especially near the connectors, as this can damage internal wires. Avoid dropping the adapter.

- **Inspection:** Periodically inspect the cable and adapter for any signs of damage, such as frayed wires, cracked casing, or bent pins. Discontinue use if damage is observed.

5. TROUBLESHOOTING

If you encounter issues with your UpBright power supply, refer to the following common solutions.

- **Device Not Charging/Powering On:**
 - Ensure the adapter is securely plugged into a working wall outlet. Test the outlet with another device.
 - Verify that both ends of the Micro USB cable are firmly connected to the adapter and the device.
 - Inspect the cable and adapter for visible damage. If damaged, discontinue use and replace.
 - If your device uses a rechargeable battery, note that batteries have a limited lifespan (typically 1 to 3 years). An old or faulty battery may not accept a charge, which is a device issue, not a charger issue.
- **Slow Charging:**
 - Ensure the device is not in heavy use while charging.
 - Confirm the cable is the original or a high-quality replacement.

6. SPECIFICATIONS

Input Voltage	100-240VAC, 50/60Hz
Output Voltage	5VDC
Output Current	500mA (0.5A)
Connector Type	Micro USB
Protection Features	OVP (Over Voltage Protection), OCP (Over Current Protection), SCP (Short Circuit Protection)
Product Dimensions	Approximately 4 x 2 x 2 inches (Adapter)
Item Weight	Approximately 3 ounces
Compatible Devices	Belifu TENS001, SM9126, AS8015, AS8012B Dual Channel Electronic Pulse Massagers

7. WARRANTY AND SUPPORT

UpBright products are designed for reliability and performance.

- **Refund Policy:** UpBright offers a 30-day refund policy from the date of purchase.
- **Exchange Policy:** The product is covered by a 24-month exchange policy for units found to be in great working condition upon testing.
- For support or warranty claims, please contact your retailer or the UpBright customer service channel provided at the point of purchase.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B09QKCCLSW>

More documents and user manuals: <https://manuals.plus/>