

[Manuals.plus](#) /

› [Veridian Healthcare](#) /

› FeverScan Forehead Strip Thermometer User Manual

Veridian Healthcare 09-116DS

FeverScan Forehead Strip Thermometer User Manual

Model: 09-116DS

INTRODUCTION

Monitor your health easily and accurately with the FeverScan Forehead Strip Thermometer. Designed for fast and convenient temperature checks, this innovative thermometer provides a 15-second readout of your body temperature without the need for batteries or complicated instructions. Featuring a liquid crystal display and light adhesive application, it ensures quick, non-invasive, and safe readings for both children and adults.

SAFETY INFORMATION

Please read all instructions carefully before use. This device is intended for personal use to indicate body temperature. It is not a substitute for professional medical advice, diagnosis, or treatment. If you have concerns about your health, consult a healthcare professional.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. While this specific disclaimer may not directly apply to this product, it is provided as general important information.

PACKAGE CONTENTS

- FeverScan Forehead Strip Thermometer (Model: 09-116DS)



Image: FeverScan Forehead Thermometer in its clear plastic packaging, showing the product name and basic information.

SETUP

The FeverScan Forehead Strip Thermometer requires no batteries or complex setup. It is ready for use directly out of the package.

1. Carefully remove the thermometer strip from its packaging.
2. Ensure the forehead area is clean and dry for optimal contact.

OPERATING INSTRUCTIONS

The FeverScan Forehead Strip Thermometer utilizes liquid crystal technology to provide a quick temperature reading.



Image: Close-up of the FeverScan strip thermometer displaying temperature ranges in both Fahrenheit and Celsius, with colored indicators for different temperature zones.

1. **Placement:** Gently press the thermometer strip flat against the center of the forehead, ensuring full contact with the skin. The strip uses a gentle adhesive for comfortable application.
2. **Reading Time:** Hold the strip in place for approximately 15 seconds.
3. **Interpreting Results:**
 - Observe the colored squares on the strip. The temperature is indicated by the square that changes color.
 - If two squares change color, the temperature is approximately halfway between the two indicated values.
 - The strip displays temperatures in both Fahrenheit (°F) and Celsius (°C).
4. **Removal:** After reading, gently peel the strip off the forehead.

Note: For best accuracy, ensure the thermometer is at room temperature before use. Avoid direct sunlight or extreme temperatures during measurement.

MAINTENANCE AND CARE

The FeverScan Forehead Strip Thermometer is designed to be reusable and durable.

- **Cleaning:** Wipe clean with a soft, damp cloth. Do not immerse in water or use harsh chemical cleaners.
- **Storage:** Store in a cool, dry place away from direct sunlight and extreme temperatures. Keep out of reach of children.
- **Durability:** The thermometer is unbreakable under normal use conditions.

TROUBLESHOOTING

Problem	Possible Cause	Solution
No color change on strip.	Insufficient contact with skin; thermometer too cold; temperature outside measurable range.	Ensure firm, flat contact with forehead. Allow thermometer to reach room temperature before use. Confirm temperature is within 95°F - 104°F (35°C - 40°C).
Inaccurate reading.	Movement during reading; improper placement; external factors (e.g., sweating, drafts).	Remain still during measurement. Reapply the strip ensuring full contact. Ensure forehead is dry. Avoid drafts.

SPECIFICATIONS

- **Model Number:** 09-116DS
- **Product Dimensions:** 3 x 0.39 x 1 inches
- **Weight:** 0.35 ounces
- **Display Type:** Analog (Liquid Crystal)
- **Readout Time:** Approximately 15 seconds
- **Material:** Rubber (Outer Material), Latex-Free
- **Age Range:** All ages
- **Manufacturer:** Veridian Healthcare
- **Country of Origin:** USA
- **Special Feature:** Fast Reading System

WARRANTY INFORMATION

Veridian Healthcare products are manufactured to exacting standards and inspected prior to shipment. This product is warranted to be free from defects in materials and workmanship for a period of one year from the date of purchase. This warranty does not cover damage caused by misuse, abuse, accident, or unauthorized alteration. For warranty service, please contact customer support.

CUSTOMER SUPPORT

For questions, assistance, or to report issues with your FeverScan Forehead Strip Thermometer, please contact Veridian Healthcare customer support:

- **Website:** [Veridian Healthcare Store on Amazon](#)
- Please refer to the product packaging or the manufacturer's official website for the most current contact information.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question