

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

› [Black Wolf](#) /

› [Wush Pro By Black Wolf - Deluxe Water Powered Ear Cleaner - User Manual](#)

## Black Wolf Wush Pro

# Wush Pro By Black Wolf - User Manual

DELUXE WATER POWERED EAR CLEANER

Model: Wush Pro

## INTRODUCTION

The Wush Pro By Black Wolf is an advanced water-powered ear cleaning device designed to safely and effectively remove ear wax buildup. Utilizing a triple jet stream and multiple pressure settings, it offers a gentle yet thorough cleaning experience. This manual provides detailed instructions for the setup, operation, maintenance, and troubleshooting of your Wush Pro device.

## WHAT'S INCLUDED

Your Wush Pro package includes the following components:

- WUSH Device (Ear Irrigator)
- Water Basin
- 6 Interchangeable Tips (color-coded for hygienic sharing)
- USB Charging Cord



*Image: The Wush Pro device, water basin, six ear tips (three black, three white), and a USB-C charging cable are displayed on a countertop.*

## SETUP

### Charging the Device

Before first use, fully charge your Wush Pro device using the provided USB charging cord. Connect the USB-C end to the device's charging port and the other end to a compatible USB power source. The battery indicator lights will show charging progress. A full charge provides approximately 3 hours of battery life.

### Filling the Water Basin

Fill the integrated water basin with warm water. It is crucial to use warm water only, as room temperature or cold water can cause dizziness or discomfort during use.



safe and effective ear  
cleaning solution



triple jet streams  
+ three pressure settings



loosen and prevent  
wax build-up

# Wush Ear Cleaner

**Note:** Use WUSH with warm water, as cold or room temperature water can cause dizziness.



*Image: The Wush Ear Cleaner device is shown with icons and text explaining its features: safe and effective cleaning, triple jet streams with three pressure settings, and its function to loosen and prevent wax build-up. A crucial note advises using warm water to avoid dizziness.*

## OPERATING INSTRUCTIONS

### Attaching the Nozzle

Select one of the interchangeable tips. The flexible, compact nozzle connects to a soft irrigation tip that fits comfortably in your ear. Ensure the tip is securely attached before use.

### Selecting Pressure Settings

The Wush Pro offers 3 pressure settings: Low, Medium, and High. Choose the setting that best suits your needs, from a gentle clean to more intensive wax removal. It is recommended to start with the lowest setting and gradually increase if needed.



### **INTENTIONAL COMFORT**

flexible and soft irrigation tip  
fits comfortably in your ear



### **COMPACT DESIGN**

foldable nozzle for sleek  
storage and device protection



### **USB RECHARGEABLE**

3 hours of battery life and  
water-resistant IPX-7 rating

*Image: A man is shown using the Wush Pro ear cleaner in the shower, illustrating its water-resistant design.*

*Accompanying icons emphasize its intentional comfort with a flexible tip, compact foldable nozzle for storage, and USB rechargeable battery with an IPX-7 water-resistant rating.*

## **Using the Device**

Once the basin is filled with warm water and a tip is attached, flip the nozzle into position. To begin cleaning, press the power button. When cleaning your right ear, tilt your head to the right to allow water and wax to flow out. Similarly, for your left ear, tilt your head to the left. The device can be used in the shower or with the provided basin placed below the ear to catch draining water.

# Triple Jet Stream Gentle Wax Remover

## SAFE AND EFFECTIVE

Gently clean your ears with a  
massaging water stream

---

## THREE PRESSURE SETTINGS

From gentle clean to wax  
removal, adjust to your  
needs and preferences

---

## UPGRADE YOUR ROUTINE

Water jets flush out wax  
without pushing it deeper

*Image: The Wush Pro ear cleaner is shown actively spraying water, demonstrating its triple jet stream. Text highlights its effectiveness for gentle wax removal and the availability of three pressure settings to suit individual needs.*

# Includes:

## WUSH Device

Safe ear wax remover

---

## 6 Interchangeable Tips

Color-coded for hygienic sharing

---

## Water Catch Cup

Over the sink use

---

## USB-C Charging Cable

Power up your routine



*Image: Two images side-by-side show a clear ear model. The first image, labeled 'BEFORE WUSH', shows visible earwax. The second image, labeled 'AFTER WUSH', shows the ear model completely clear of wax, illustrating the device's cleaning effectiveness.*

## MAINTENANCE

### Cleaning the Device

The Wush Pro is water resistant and can be safely rinsed under running water. After each use, detach the tip and rinse both the tip and the nozzle thoroughly. Empty and rinse the water basin. Ensure all components are dry before storing.

### Storing the Device

The compact nozzle folds down for easy storage, protecting the device when not in use. Store the Wush Pro in a clean, dry place away from direct sunlight and extreme temperatures.

## TROUBLESHOOTING

Issue	Possible Cause / Solution
Device not turning on or charging.	Ensure the device is fully charged. Check the USB charging cable and power source connection.
Water not spraying or weak spray.	Check if the water basin is sufficiently filled. Ensure the nozzle and tip are clear of any obstructions. Select a higher pressure setting if desired.
Experiencing dizziness during use.	This can occur if cold or room temperature water is used. Always use warm water for ear irrigation to prevent dizziness.

## SPECIFICATIONS

- **Product Dimensions:** 10.35 x 4.37 x 3.42 inches
- **Item Weight:** 1.06 Pounds
- **Battery Life:** Approximately 3 hours of continuous use
- **Water Resistance:** IPX-7 rating
- **Pressure Settings:** 3 (Low, Medium, High)
- **Charging:** USB Rechargeable (USB-C)

## SAFETY INFORMATION

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition. Consult a healthcare professional if you have any pre-existing ear conditions, ear infections, or concerns about ear wax removal.

## WARRANTY & SUPPORT

For information regarding product warranty, returns, or customer support, please refer to the official Black Wolf website or contact their customer service directly. Each item has a unique code that is verified before shipping.