

[Manuals.plus](#) /

› [QUINEAR](#) /

› QUINEAR Leg Massager QN-072A-BU Instruction Manual

QUINEAR QN-072A-BU

QUINEAR Leg Massager QN-072A-BU Instruction Manual

Model: QN-072A-BU

INTRODUCTION

The QUINEAR Leg Massager QN-072A-BU is designed to provide air compression massage and heat therapy for your feet, calves, and thighs. Utilizing airwave pressure therapy, it inflates and deflates 13 airbags to create a circulating pressure that promotes blood circulation, reduces restless leg syndrome (RLS), alleviates swelling, and relieves muscle soreness. This device is suitable for individuals seeking relaxation and relief, including sports enthusiasts, office workers, and the elderly.

SAFETY INFORMATION

- Use the massager responsibly, avoiding extended use.
- Adjust compression level to personal comfort.
- Be cautious when using the heat feature to prevent burns.
- Consult a doctor before use if you are receiving medical treatment or have pre-existing health conditions.

PACKAGE CONTENTS

- 2 Leg Massager Wraps (Left & Right)
- 1 Handheld Digital Controller
- 1 Power Adapter
- 2 Extender Wraps
- 1 Portable Storage Bag
- 1 User Manual

PRODUCT OVERVIEW



The QUINEAR Leg Massager system, including the wraps and handheld controller.

Better Compression System



3 Heating Mode



4 Heating Areas



10 mins Fast Heating



Insulated Cotton Material



20 mins Auto Off



Before



After



A user comfortably wearing the QUINEAR Leg Massager for a full leg compression experience.



360° Fully Wrap Compression



Flexible Velcro Adjustment



High Tenacity Silicone Air Hose



Finely Sewn Fabrics



Handheld Controller Connection

Close-up of the handheld digital controller, showing the display and control buttons for modes, intensity, and heat.

Air Compression Massager




13*2 Air Bags



4 Heating Area & 3 Heating Level




360° Fully Wrapped Compression



The leg massager wraps with the included extender pieces, allowing for adjustable sizing.

MORE SUITABLE SIZE DESIGN


2 extensions to expand the size for bigger legs



With 2 Extensions
Thighs up from 27" to 40"
Calves up from 19" to 28.5"




The QUINEAR Leg Massager neatly packed into its portable storage bag for convenience.



Visual representation of key features: 3 heating modes, 4 heating areas, 10 minutes fast heating, insulated cotton material, and 20 minutes auto-off timer.

Sequence Mode / Circulation Mode / Combination Mode



Detailed diagram illustrating the 13x2 airbags, 4 heating areas, 3 heating levels, and 360° fully wrapped compression design.



Illustrations explaining the three distinct massage modes: Sequence Mode, Circulation Mode, and Combination Mode.

Your browser does not support the video tag.

Official QUINEAR Leg Massager Overview: This video provides a general overview of the QUINEAR Leg Massager, highlighting its features and benefits for relaxation.

Your browser does not support the video tag.

QUINEAR Leg Recovery System Features: This video demonstrates the various features of the QUINEAR Leg Recovery System, showcasing its professional capabilities.

SETUP INSTRUCTIONS

1. Unpack the QUINEAR Leg Massager from its packaging.
2. Connect the massager to a power source.
3. Place the leg sleeves on your legs, ensuring a snug fit. Use the 2 extender wraps if needed for larger leg circumferences. The wraps are labeled 'L' for left and 'R' for right.
4. Connect the air hoses from the leg wraps to the handheld controller. Ensure the 'UP' indicators on the connectors align with the controller ports.

OPERATING INSTRUCTIONS

The handheld digital controller allows you to customize your massage experience. Follow these steps to operate the device:

1. **Power On:** Press the power button on the handheld controller to turn on the device. The display will illuminate.
2. **Select Mode:** Use the 'M' button to cycle through the three compression modes:
 - **Sequence Mode:** Airbags inflate and deflate sequentially from feet to thighs.
 - **Circulation Mode:** Airbags inflate and deflate in a circulating pattern.
 - **Combination Mode:** Alternates between Sequence and Circulation modes.
3. **Adjust Intensity:** Use the 'Strength' button (up arrow icon) to select one of three intensity levels: Minimum (Min), Medium (Mid), or Maximum (Max).
4. **Set Heat Level:** Use the 'Heat' button (wavy lines icon) to choose from three heating levels: Minimum (Min), Medium (Mid), or Maximum (Max). You can also select specific heating areas (feet, calves, thighs) if available on your model.
5. **Start Massage:** Once settings are selected, the massage session will begin automatically. The device has a built-in 20-minute auto-shutoff feature for safety.
6. **Power Off:** Press and hold the power button to turn off the device manually.

MAINTENANCE

To ensure the longevity and hygiene of your QUINEAR Leg Massager, follow these maintenance guidelines:

- After each use, clean the leg sleeves with a damp towel. The nylon waterproof fabric is designed for easy cleaning.
- Store the massager in a dry, cool place, away from direct sunlight.
- Do not immerse the controller or any electrical components in water.

TROUBLESHOOTING

If you encounter any issues with your QUINEAR Leg Massager, please refer to the following common troubleshooting tips:

- **Device not powering on:** Ensure the power adapter is securely connected to both the controller and a working power outlet.
- **No compression:** Check that all air hoses are firmly connected to both the leg wraps and the controller. Ensure the 'UP' indicators are correctly aligned.
- **No heat:** Verify that a heat level (Min, Mid, or Max) is selected on the controller. Allow a few minutes for the heating elements to warm up.
- **Uneven compression:** Ensure the leg wraps are applied snugly and evenly to both legs. Adjust the Velcro straps as needed.
- If problems persist, consult the full user manual or contact customer support.

SPECIFICATIONS

Model Number	QN-072A-BU
Package Dimensions	15.86 x 8.82 x 7.05 inches
Item Weight	4.85 Pounds
Power Source	Corded Electric
Material	Carbon Fiber
Manufacturer	QUINEAR

WARRANTY AND SUPPORT

For warranty information or customer support, please refer to the contact details provided in the product packaging or visit the official QUINEAR website. Keep your purchase receipt for warranty claims.

LEGAL DISCLAIMER

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question