

VAAGHANM NA-H1223C

VAAGHANM 33 x 17 inches Electric Heating Pad (Model NA-H1223C) Instruction Manual

Comprehensive guide for safe and effective use.

INTRODUCTION

Thank you for choosing the VAAGHANM 33 x 17 inches Electric Heating Pad. This manual provides essential information for the safe and proper operation, maintenance, and care of your new heating pad. Please read all instructions carefully before use and retain this manual for future reference.

IMPORTANT SAFETY INFORMATION

It is crucial to follow these safety guidelines to prevent injury or damage to the product:

- Avoid folding the heating pad while in use.
- Always detach and cool the controller before washing.
- Ensure the pad is completely dry before turning it on again.
- Keep out of reach of children and pets.
- Do not use if damaged or malfunctioning.
- Always follow the manufacturer's instructions for safe use and storage.

PRODUCT FEATURES

The VAAGHANM Electric Heating Pad is designed for comfort and ease of use, featuring:

- **Extra Large Size:** 33 x 17 inches (84cm x 43cm) for broad coverage.
- **Ultra-Soft Flannel Material:** Provides a comfortable and gentle touch.
- **6 Heat Settings:** Adjustable temperature from 105°F to 140°F (40°C to 60°C).
- **6 Timer Settings:** Automatic shut-off options (0.5H, 1H, 1.5H, 2H, 2.5H, 3H) for safety and convenience.
- **Automatic Shut-Off:** Enhances safety by turning off after the set time.
- **Overheat Protection:** Built-in safety mechanism to prevent overheating.

- **Long Detachable Cord:** A 9-foot (118-inch) cord for flexible use and easy storage.
- **Machine Washable:** Easy to clean after detaching the controller.



Image: The VAAGHANM 33x17 inch Electric Heating Pad in gray, showing the soft flannel material and the attached controller.

SETUP

1. **Unpack:** Carefully remove the heating pad and controller from its packaging.
2. **Inspect:** Check the heating pad and cord for any signs of damage. Do not use if damaged.
3. **Connect Controller:** Ensure the controller is securely connected to the heating pad. The connector should be fully inserted.
4. **Plug In:** Insert the power plug into a standard electrical outlet.

OPERATING INSTRUCTIONS

The heating pad features an intuitive LED controller for easy operation.

Multi-Functional LED Controller



Image: Close-up of the multi-functional LED controller, highlighting the display screen, timer button, temperature buttons, and power switch.

1. **Power On:** Press the **Power Switch** button on the controller to turn on the heating pad. The LED display will illuminate.
2. **Adjust Heat Setting:** Use the **Temperature buttons** ("**+**" or "**-**") to select your desired heat level from 1 to 6. Level 1 is the lowest (105°F/40°C) and Level 6 is the highest (140°F/60°C).
3. **Set Timer:** Press the **Timer Button** to cycle through the automatic shut-off options: 0.5H, 1H, 1.5H, 2H, 2.5H, or 3H. The pad will automatically turn off after the selected time. If no timer is set, the pad will default to a 2-hour auto-off for safety.
4. **Power Off:** Press the **Power Switch** button again to turn off the heating pad. Always unplug the pad from the electrical outlet when not in use.

Note: The heating pad heats up quickly. Adjust settings as needed for your comfort.

INDICATIONS FOR USE

This heating pad is intended for temporary relief of minor aches and pains associated with muscles and joints. It can be used for:

- Back pain

- Shoulder pain
- Cramp relief
- Menstrual pain
- Enhancing blood circulation
- General warmth and comfort



Image: A woman comfortably using the VAAGHANM heating pad on her lower back while seated on a couch, demonstrating its large size and flexible application.

CARE AND MAINTENANCE

Proper care ensures the longevity and safe operation of your heating pad:

1. **Disconnect:** Always detach the controller from the heating pad before cleaning.
2. **Machine Washable:** The heating pad is machine washable. Use cold water on a gentle cycle.
3. **Drying:** Tumble dry on low heat or air dry. Ensure the heating pad is **completely dry** before reconnecting the controller and using it again.
4. **Storage:** Roll up the heating pad with the controller inside and secure it with the elastic straps for easy and compact storage.



Image: A woman demonstrating how to place the VAAGHANM heating pad into a washing machine after detaching the controller, highlighting its machine-washable feature.

PRODUCT SPECIFICATIONS

| Feature | Detail |
|-------------------|--|
| Model Number | NA-H1223C |
| Dimensions | 33 x 17 inches (84cm x 43cm) |
| Material | Flannel |
| Heat Settings | 6 levels (105°F - 140°F / 40°C - 60°C) |
| Timer Settings | 6 options (0.5H, 1H, 1.5H, 2H, 2.5H, 3H) |
| Cord Length | 9 feet (118 inches) |
| Special Features | Automatic Shut Off, Overheat Protection, Portable, Timer |
| Power Consumption | 120 Watts (for Large and Double XL sizes, current model is XL) |

Weight

1.8 Pounds

Extra large: **84cm*43cm** in




Image: The VAAGHANM heating pad laid flat, illustrating its 33x17 inch dimensions and the length of its detachable power cord.

OFFICIAL PRODUCT VIDEOS

Watch these official videos from VAAGHANM for a visual guide on using your heating pad.

VAAGHANM Large Heating Pad Overview



Video: A demonstration of the VAAGHANM large heating pad, showcasing its features and size. Duration: 0:38.

VAAGHANM Heating Pad for Back Use



Video: This video illustrates the VAAGHANM heating pad being used for back relief, demonstrating its flexibility and coverage. Duration: 0:37.

LEGAL DISCLAIMER

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.