

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

> [veicomtech](#) /

> [veicomtech VEICOMTECH-# Ionic Foot Bath Detox Machine User Manual](#)

veicomtech VEICOMTECH-#

veicomtech VEICOMTECH-# Ionic Foot Bath Detox Machine User Manual

Model: VEICOMTECH-#

[Contents](#) [Setup](#) [Introduction](#) [Safety Information](#) [Package](#)
[Operation](#) [Maintenance](#) [Troubleshooting](#) [Specifications](#) [Warranty & Support](#)

1. INTRODUCTION

This manual provides comprehensive instructions for the safe and effective use of your veicomtech VEICOMTECH-# Ionic Foot Bath Detox Machine. This 3-in-1 system includes a foot spa, a far infrared waist belt, and a wrist strap, designed for personal relaxation and wellness. Please read this manual thoroughly before initial use and retain it for future reference.

2. SAFETY INFORMATION

Important Safety Precautions:

- Use the massager strictly according to the instructions provided in this manual.
- Do not use any substances other than water and a small amount of salt in the foot bath.
- Handle the machine and its components with care to prevent damage.
- Discontinue use immediately if you experience any discomfort, pain, or adverse reactions during operation.
- Keep the device out of reach of children.
- Do not immerse the main control unit in water or any other liquid.
- Ensure your hands are dry when plugging in or unplugging the power cord.
- Consult a healthcare professional before use if you have any medical conditions, are pregnant, or have electronic implants.

3. PACKAGE CONTENTS

Please verify that all items listed below are included in your package:

- 1 x Ionic Detox Foot Bath Machine (Main Control Unit)

- 1 x Far Infrared Waist Belt
- 2 x Detox Arrays
- 1 x Wrist Strap
- 1 x Power Cord
- 2 x TENS Pads
- 1 x Foot Tub
- 1 x Blue Carry Bag
- 10 x Tub Liners
- 1 x Instruction Manual (this document)



Image 3.1: Complete package contents of the veicomtech Ionic Foot Bath Detox Machine.

4. SETUP INSTRUCTIONS

Follow these steps to set up your Ionic Foot Bath Detox Machine:

1. **Placement:** Place the foot tub on a flat, stable surface. Ensure the main control unit is placed nearby on a dry, stable surface, away from water splashes.
2. **Fill Tub:** Fill the foot tub with warm water, ensuring it is deep enough to cover your feet comfortably.
3. **Add Salt:** Add a small amount of non-iodized salt (approximately 1/4 to 1/2 teaspoon) to the water. Stir until dissolved. The salt helps to conduct the electrical current.
4. **Connect Array:** Connect one or two detox arrays to the designated ARRAY port on the main control unit. Place the array(s) into the water in the foot tub.
5. **Connect Wrist Strap/Waist Belt:** If using, connect the wrist strap to the 'Wrist Strap' port or the far infrared waist belt to the 'Belt' port on the main control unit. Secure the wrist strap around your wrist or the waist belt around your waist.
6. **Connect Power:** Plug the power cord into the main control unit's power input, then plug the other end into a standard electrical outlet.



SHARE FOOT SPA LEISURE WITH FAMILY

Image 4.1: Proper setup with feet in the tub, array, and wrist strap connected.



Image 4.2: Rear connections of the main control unit.

5. OPERATION

Once the machine is set up, follow these steps to begin your session:

1. **Power On:** Press the 'ON/OFF' button on the control panel to turn on the machine. The display will illuminate.
2. **Select Mode:** Press the 'Detox Mode' button to cycle through the available modes. The VEICOMTECH-# typically offers five modes:
 - **Ion Detox Mode:** Focuses on ionizing the water.
 - **Hydrogen Spa Mode:** Combines ionization with hydrogen therapy.
 - **Relaxation Mode:** Provides a gentle experience without specific detoxification.
 - **Infrared Therapy Mode:** Activates the far infrared waist belt for warmth and circulation.
 - **Custom Mode:** Allows for personalized settings.
3. **Set Time:** Press the 'Detox Time' button to adjust the session duration. A typical session lasts 30 minutes.
4. **Start Session:** Once the mode and time are set, place your feet into the water. The machine will begin

operation.

5. **Monitor Amperage:** Observe the amperage (A) reading on the display screen. It should ideally stay within the recommended range (e.g., 1.7 to 2.2 amps). If the amperage is too low, you may need to add a small pinch more salt. If it's too high, dilute the water slightly.
6. **Session Completion:** The machine will automatically turn off when the set time expires. Remove your feet from the tub and dry them thoroughly.




Image 5.1: Control panel for selecting mode and time.



Image 5.2: Steps for setting time, mode, and starting the machine.

Understanding Water Color Change

During a session, the water in the foot tub may change color. This phenomenon is primarily due to the interaction between the array's metal, the electrical current, and the minerals present in the water itself, as well as any impurities from your feet. The color change is a natural part of the electrochemical process and does not solely indicate the release of toxins from the body. The primary benefit of the ionic foot bath is believed to be through the ionization process in the water.

PROCESS OF IONIC DETOXIFICATION



Image 5.3: Visual representation of water color change during a session.

6. MAINTENANCE

Proper maintenance ensures the longevity and optimal performance of your device.

Cleaning the Array

- After each use, disconnect the array from the main unit.
- Rinse the array thoroughly with clean water to remove any residue.
- For deeper cleaning, soak the array in a solution of 50% white vinegar and 50% water for 10-15 minutes. Use a brush to gently scrub away any buildup.
- Rinse again with clean water and allow it to air dry completely before storage.

Cleaning the Foot Tub

- Empty the water from the tub after each session.
- Wash the tub with mild soap and water.

- Rinse thoroughly and dry completely before storing. Using tub liners can simplify cleaning.

Cleaning the Main Control Unit

- Ensure the unit is unplugged before cleaning.
- Wipe the exterior of the unit with a soft, damp cloth.
- Do not use abrasive cleaners or immerse the unit in water.

Storage

Store all components in a cool, dry place, away from direct sunlight and extreme temperatures. Use the provided blue carry bag for organized storage and transport.

7. TROUBLESHOOTING

If you encounter issues with your veicomtech Ionic Foot Bath Detox Machine, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Machine does not power on.	Power cord not properly connected; power outlet issue.	Ensure power cord is securely plugged into both the machine and a working electrical outlet. Test the outlet with another device.
Low or no amperage reading.	Insufficient salt in water; faulty array connection.	Add a small pinch more non-iodized salt to the water and stir. Ensure the array is securely connected to the ARRAY port.
High amperage reading.	Too much salt in water.	Remove a small amount of water and add fresh water to dilute the salt concentration.
Water does not change color.	Insufficient salt; array not functioning; water mineral content.	Ensure adequate salt is added. Verify the array is connected and functioning. Note that water color change is primarily due to the interaction of the array with water minerals, not solely body detox.
Screen is cracked or unreadable.	Physical damage.	Contact customer support for assistance.
Discomfort during use.	Individual sensitivity; incorrect settings.	Discontinue use immediately. Consult a healthcare professional if discomfort persists. Review settings and ensure proper usage.

8. SPECIFICATIONS

Feature	Detail
Brand	veicomtech
Model Number	VEICOMTECH-#
Product Dimensions	14 x 8 x 4 inches
Item Weight	10 Pounds

Feature	Detail
Material	Aluminum (main unit casing)
Power Source	Corded Electric
UPC	751271408155
ASIN	B0B5KGDS7W

9. WARRANTY & SUPPORT

For warranty information, technical support, or any questions regarding your veicomtech Ionic Foot Bath Detox Machine, please refer to the contact information provided on the product packaging or visit the official veicomtech website. Please retain your proof of purchase for warranty claims.