

## FlyingJoy 68-black

# FlyingJoy Heavy-Duty Rollator Walker with Seat

MODEL: 68-BLACK

## Introduction

This instruction manual provides essential information for the safe and effective use, assembly, and maintenance of your FlyingJoy Heavy-Duty Rollator Walker with Seat, Model 68-black. Please read this manual thoroughly before initial use and retain it for future reference. This rollator is designed to provide mobility assistance and support for individuals requiring walking aid, offering a stable frame, comfortable seating, and convenient storage.

## Safety Information

Your safety is paramount. Adhere to the following guidelines to prevent injury and ensure proper operation:

- Always ensure the rollator is fully unfolded and locked before use.
- Before sitting, engage the parking brakes to prevent the rollator from moving.
- Do not push the rollator while sitting on it. It is intended for walking assistance and temporary seating.
- Regularly inspect all components, including brakes, wheels, and frame, for wear or damage. Do not use if any part is compromised.
- After adjusting the handle or seat height, ensure that the adjustment pins are securely engaged in the corresponding holes and the knobs are tightened.
- Be aware of your surroundings, especially in low-light conditions. The rollator features high-visibility reflective strips for added safety.
- The maximum weight capacity for this rollator is 450 lbs (204 kg). Do not exceed this limit.
- Use the handgrips to assist walking; do not lean excessively on the frame.

## Components Overview

The FlyingJoy Heavy-Duty Rollator Walker consists of the following main components:

- **Frame:** Reinforced aluminum frame with a triangular design for stability.
- **Handles:** Ergonomically designed handgrips with integrated brake levers.

- **Braking System:** Dual-brake system for controlled stopping and parking.
- **Wheels:** Four 8-inch airless tires, with 360-degree rotatable front wheels for maneuverability.
- **Seat:** Extra-wide (22 inches / 55 cm) padded seat with an air-permeable cushion.
- **Backrest:** Padded backrest for comfort when seated.
- **Storage Bag:** Large capacity under-seat storage bag (up to 11 lbs / 5 kg).



Image: The FlyingJoy Heavy-Duty Rollator Walker, Model 68-black, showcasing its main components including the frame, wheels, seat, backrest, and handles.

## Setup and Assembly

---

The rollator is designed for easy assembly. Follow these steps for initial setup:

1. **Unpacking:** Carefully remove all components from the packaging. Inspect for any visible damage or missing parts.
2. **Unfolding the Frame:** Place the rollator on a flat surface. Lift the seat upwards, then pull the folding lever to fully extend the frame until it locks into place. Ensure all locking mechanisms are engaged.
3. **Handlebar Installation/Adjustment:** Insert the handlebars into the frame. Adjust the height by pressing the push-button (marble) and sliding the handlebar to your desired height. There are 6 levels for handlebar adjustment (33.27"-41.49"). Ensure the adjustment pins snap securely into place and the knobs are tightened.
4. **Seat Height Adjustment (if applicable):** The seat height can be adjusted to 5 levels (19.48"-23.03"). Press the push-button and adjust the seat to a comfortable position, ensuring it locks securely.
5. **Brake Check:** Test the brakes by pushing down on the brake levers to engage the parking brake. The wheels should lock. Release the levers to disengage.

SPACIOUS **26"** HEAVY DUTY WALKERS

## EASY TO FOLD AND STORAGE

STORAGE DOES NOT TAKE UP SPACE



Image: The rollator in its folded state, illustrating its compact dimensions for storage and transport. It also shows the under-seat storage bag.

## FOLDABLE AND EASY TO ASSEMBLE

Convenient to carry and store,  
does not take up space



Image: A close-up view of the one-touch height adjustment mechanism, showing how to press the button to adjust the handgrip and seat height.

## Operating Instructions

Proper operation ensures safety and comfort.

### Walking with the Rollator

- Stand upright between the handles, maintaining a comfortable posture.
- Push the rollator forward a short distance, then step forward with your weaker leg, followed by your stronger leg.
- Keep your weight centered over the rollator for stability.
- The 8-inch airless tires and 360-degree rotatable front wheels provide smooth movement across various terrains, including hardwood, carpet, gravel, grass, dirt, and sidewalks.

THE STRONGEST CONSTRUCTION

## ALL-ROUND PROTECTION FOR SENIORS SAFETY

FLYINGJOY 20 YEARS OF EXPERIENCE

\*SUPPORTS UP TO  
450 LBS/177 KG

\*DOUBLE SAFETY,  
DOUBLE COMFORT

\*RIGOROUS US CERTIFICATIONS  
& QUALITY ASSURANCE



Image: A senior man confidently walking outdoors with the FlyingJoy Rollator Walker, demonstrating its use in a natural environment.

# ALL TERRAIN STABILITY

Enjoy smooth, stable, and quiet rides everywhere

Safe to Go Anywhere!

8-inch airless tire  
shock absorption function




Gravel



Uneven



Grass



Sand



Image: The rollator being used on different surfaces like gravel, uneven ground, grass, and sand, highlighting its all-terrain stability with 8-inch airless tires.

## Using the Brakes

- **To Slow Down or Stop:** Squeeze the brake levers gently towards the handles.
- **To Engage Parking Brakes:** Push the brake levers down until they lock into place. This is crucial before sitting on the seat.
- **To Release Parking Brakes:** Pull the brake levers upwards until they disengage.



Image: A visual demonstration of the dual-brake system, showing how to engage the hand brakes for slowing down and how to push them down to lock for parking.

### Using the Seat

- Ensure the rollator is on a level surface and the parking brakes are fully engaged.
- Turn around and slowly lower yourself onto the padded seat.
- The extra-wide 22-inch (55 cm) seat and padded backrest provide comfort for resting.

# DOUBLE COMFORT

NO MORE  
CRAMPED SEATS!

EXTRA-WIDE  
SEAT & BACKREST

*Reducing Pressure  
on Hips and Back*



Image: A detailed view of the rollator's extra-wide seat and backrest, highlighting the 22-inch backrest width and 23.62-inch seat depth for enhanced comfort.



## DUAL-HEIGHT ADJUSTMENT MODE

Different heights, same comfort, creating a personalized travel experience for you.



Adjustable handle



Adjustable seat

Image: A close-up of the air-permeable cushion design on the seat, indicating its breathable and quick-drying properties for extended comfort.

### Using the Storage Bag

- Lift the seat to access the spacious under-seat storage bag.
- The bag can hold up to 11 lbs (5 kg) of personal items.



Image: The super large capacity storage bag located under the seat, shown filled with various items, demonstrating its practical use for carrying personal belongings.

### Folding for Transport and Storage

- Ensure the storage bag is empty.
- Lift the seat upwards.
- Pull the folding lever (typically located under the seat or on the frame) to collapse the rollator.
- The rollator folds compactly for easy transport and storage.

### Maintenance

Regular maintenance helps prolong the life of your rollator:

- **Cleaning:** Wipe the frame and seat with a damp cloth and mild detergent. Avoid abrasive cleaners.
- **Inspection:** Periodically check all screws, bolts, and locking mechanisms to ensure they are tight and secure.
- **Brakes:** Test the brakes regularly to ensure they are functioning correctly. If they feel loose or ineffective, consult a qualified technician.
- **Wheels:** Keep the wheels free of debris. The airless tires require no inflation.

- **Storage:** Store the rollator in a dry place away from extreme temperatures.

## Troubleshooting

If you encounter issues with your rollator, consider the following:

- **Rollator feels unstable:** Ensure all frame components are fully extended and locked. Check that height adjustment pins are securely engaged.
- **Brakes are not engaging/disengaging properly:** Inspect the brake cables for damage or looseness. Ensure the brake levers are fully pushed down to lock or pulled up to release. If issues persist, contact customer support.
- **Difficulty folding/unfolding:** Ensure the seat is lifted and the folding lever is correctly actuated. Check for any obstructions.

## Specifications

| Feature                       | Specification                                              |
|-------------------------------|------------------------------------------------------------|
| Product Dimensions (Unfolded) | 15.94"D x 25.98"W x 46.6"H (40.5cm D x 66cm W x 118.4cm H) |
| Item Weight                   | 23.5 Pounds (10.66 kg)                                     |
| Material                      | Aluminum                                                   |
| Color                         | 68-black                                                   |
| Maximum Weight Capacity       | 450 lbs (204 kg)                                           |
| Adjustable Handle Height      | 33.27"-41.49" (84.5cm-105.4cm)                             |
| Adjustable Seat Height        | 19.48"-23.03" (49.5cm-58.5cm)                              |
| Seat Width                    | 22" (55cm)                                                 |
| Wheel Size                    | 8-inch airless tires                                       |

SPACIOUS **26"** HEAVY DUTY WALKERS

## 26" EXTRA-WIDE SIZE

PLEASE ALLOW 1-2CM ERROR AS FOR MANUAL MEASUREMENT

### Notice

**NOT FOR  
STANDARD  
DOOR**

Measure  
Before Purchase



Image: A diagram illustrating the key dimensions of the FlyingJoy Rollator Walker, including height, width, and depth, along with its weight.

### Warranty and Support

The FlyingJoy Heavy-Duty Rollator Walker is US certified and comes with a lifetime warranty, reflecting our commitment to quality and safety. Each unit undergoes over 100 quality checks during manufacturing.

For any questions, concerns, or support needs, please contact our customer service team. We are committed to resolving any issues within 24 hours.

**Contact Information:** Please refer to the product packaging or the official FlyingJoy website for the most current customer support contact details.

### Official Product Videos

No official seller-created videos were found for this product in the provided data.

© 2025 FlyingJoy. All rights reserved.

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

**Ask a question about this manual**

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

**Name**

Anonymous

**Email**

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

**Question**

Example: How do I reset this device or fix this error?

**Details**

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question