

[Manuals.plus](#) /

› [ALLJOY](#) /

› ALLJOY Leg Air Compression Massager User Manual

ALLJOY AJFLM22

ALLJOY Leg Air Compression Massager User Manual

Model: AJFLM22 | Brand: ALLJOY

INTRODUCTION

The ALLJOY Leg Air Compression Massager is designed to provide comprehensive air pressure massage, heating, and vibration for your feet, calves, and thighs. It aims to relieve leg fatigue, improve blood circulation, and offer a relaxing massage experience. This manual provides detailed instructions for safe and effective use of your new massager.



Figure 1: ALLJOY Leg Air Compression Massager and included accessories.

Who Can Benefit?

For Each One of You



Workers

The Elders



Athletes



Leg diseases



Figure 2: The massager is suitable for various users, including those who stand for long periods, the elderly, and athletes.

IMPORTANT SAFETY INFORMATION

- Consult a doctor before use if you have a medical condition, are pregnant, or have sensitive skin.
- Do not use if you have a pacemaker, diabetes, or an artificial joint.
- Do not use on open wounds or swollen areas.
- Keep out of reach of children.
- Do not immerse the device in water.

- The device is equipped with a 30-minute auto-shut-off timer to prevent overheating and overuse.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

PRODUCT COMPONENTS

Your ALLJOY Leg Massager comes with the following components:

- Full Leg Massager Sleeves (Left and Right)
- Handheld LCD Controller
- Power Adapter
- 2 Extension Wraps
- Storage Bag

Product Details



Figure 3: Various parts of the ALLJOY Leg Massager, including sleeves, controller, and connections.

Handheld LCD Controller

The intuitive LCD controller allows you to easily select massage modes, adjust intensity, and control heat and vibration functions.

Upgraded Remote Controller

High quality & Easy to use



3 Vibration



2 Heat Levels



6 Air compression Mode



Long press:
Turn on/off
Short press:
Adjust intensity

Figure 4: The upgraded remote controller features touch buttons for various settings.

Adjustable Sleeves and Extensions

The massager sleeves feature adjustable velcro wraps to fit various leg sizes. Two extension wraps are included to accommodate larger calves and thighs, ensuring a comfortable fit for the whole family.



How to connect the controller with air hose?

Insert the words "R" and "L" on the hose into the corresponding "R" and "L" holes in the controller respectively.

Figure 5: Max calves size can be 25 inches, and thighs size can be up to 35 inches with extensions.



Figure 6: Strong velcro ensures a secure and adjustable fit.

Power Adapter

The massager is powered by a safety low voltage AC adapter.




Figure 7: The power adapter ensures safe operation.

SETUP GUIDE

Follow these steps to set up your ALLJOY Leg Massager:

- 1. **Wrap your foot:** Place your foot into the foot section of the massager sleeve.
- 2. **Wrap your calf:** Secure the calf section around your calf using the velcro straps.
- 3. **Wrap your thigh:** Secure the thigh section around your thigh using the velcro straps.
- 4. **Wrap your knee:** Ensure the knee section is positioned correctly over your knee.
- 5. **Plug in controller:** Connect the air hoses from the massager sleeves to the handheld controller. Ensure the 'R' and 'L' indicators on the hose match the corresponding holes on the controller.
- 6. **Connect the adapter:** Plug the power adapter into the controller and then into a power outlet.
- 7. **Press the power button:** Long press the power button on the controller to turn on the device.



Figure 8: Step-by-step visual guide for installation.



Figure 9: Proper connection of air hoses to the controller.

OPERATING INSTRUCTIONS

Once set up, use the handheld controller to customize your massage experience.

Massage Modes and Intensities

The massager offers 6 different massage modes and 3 intensities (Min, Mid, Max) to suit your preferences. Modes include Auto, Leg, Foot, Pinch, Press, and Relax.

6 Massage Modes & 2 Heat Levels

Relieving leg fatigue and improving blood circulation



Figure 10: Choose from various modes and intensity levels for personalized relief.

Heat and Vibration Functions

Enjoy soothing warmth with 2 heating levels (42°C, 52°C) for your soles and knees. The massager also provides a vibration function for your calves.

Video 1: Demonstrates the ALLJOY Leg Massager with heat and compression features. (Source: Alljoy Direct Store)

Auto-Shut-Off Function

For your safety and to prevent overheating, the massager is equipped with a 30-minute auto-shut-off timer. If you wish to continue massaging, simply restart the device.

MAINTENANCE

- **Cleaning:** Wipe the massager sleeves with a soft, damp cloth. Do not use abrasive cleaners or immerse in water.
- **Storage:** After use, fold the massager neatly and store it in the provided storage bag. Keep it in a cool, dry place away from direct sunlight.




Figure 11: Follow the steps for proper storage of the massager.

TROUBLESHOOTING

Problem	Possible Cause	Solution
Device does not turn on	Not properly plugged in; Power adapter issue	Ensure all connections are secure. Check power outlet.
Weak compression	Sleeves too loose; Low intensity setting	Tighten velcro straps. Increase intensity level on controller.
No heat	Heat function off; Device not fully warmed up	Activate heat function on controller. Allow a few minutes for heating elements to warm up.
Device shuts off unexpectedly	Auto-shut-off timer activated	This is a safety feature. Restart the device if you wish to continue.

SPECIFICATIONS

Feature	Detail
---------	--------

Feature	Detail
Product Dimensions	16.14 x 9.84 x 8.27 inches
Item Weight	4.41 Pounds (2 Kilograms)
Model Number	AJFLM22
Manufacturer	ALLJOY
Material	Nylon
Power Source	Corded Electric
Auto-Shut-Off	30 minutes
Heat Levels	2 (42°C, 52°C)
Massage Modes	6
Intensity Levels	3

WARRANTY AND SUPPORT

ALLJOY provides a 365-day worry-free warranty for this product. If you encounter any problems or have questions regarding further use, please contact our customer support team.

Customer Support: For assistance, please refer to the contact information provided on the product packaging or the official ALLJOY website. Our team aims to provide professional answers within 24 hours.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question