

helavo H1050-Serie

Helavo Replacement Bag Instruction Manual

For Helavo H1050 and H1150 Rolling Walker Series

1. INTRODUCTION

This manual provides essential information for the proper installation, use, and maintenance of your Helavo Replacement Bag. This bag is designed as a direct replacement part for specific Helavo rolling walker models, ensuring a perfect fit and continued functionality.



Image 1: The Helavo Replacement Bag. This image shows the black fabric bag, designed to attach to a rolling walker frame. It features multiple fabric loops along the top edge for secure attachment.

2. COMPATIBILITY

The Helavo Replacement Bag is specifically designed to fit the following Helavo rolling walker models:

- Helavo H1050 Rolling Walker Series
- Helavo H1150 Rolling Walker Series

To confirm your walker's model number, please check the sticker located on the frame of your walker. Using this replacement bag with walkers from other manufacturers is not recommended as it may not fit correctly or securely.

3. SETUP AND INSTALLATION

Follow these steps to properly install the replacement bag onto your Helavo rolling walker:

1. **Prepare the Walker:** Ensure your rolling walker is in a stable, upright position. If replacing an old bag, remove it completely.
2. **Orient the Bag:** Position the new replacement bag so that the opening faces upwards and the attachment loops align with the designated attachment points on your walker's frame.
3. **Attach the Loops:** Carefully thread each fabric loop on the bag through the corresponding attachment points (e.g., hooks, bars, or clips) on the walker frame. Ensure all loops are securely fastened.
4. **Verify Security:** Gently tug on the bag to ensure it is firmly attached and does not sag or come loose. The bag should sit snugly within the walker's frame.

If you encounter any difficulties during installation, please refer to your original walker's instruction manual or contact Helavo customer service.

4. OPERATING AND USAGE

The Helavo Replacement Bag is designed to provide convenient storage for personal items while using your rolling walker. It can hold items such as groceries, books, or other small belongings.

- **Weight Capacity:** Avoid overloading the bag. While durable, excessive weight can affect the walker's stability and the bag's longevity. Refer to your walker's manual for overall weight limits.
- **Even Distribution:** Distribute items evenly within the bag to maintain balance and prevent tipping.
- **Secure Contents:** Ensure items placed in the bag are secure and will not fall out during movement.

5. MAINTENANCE

Proper care will extend the life of your Helavo Replacement Bag:

- **Cleaning:** Wipe the bag clean with a damp cloth and mild soap. Do not use harsh chemicals or abrasive cleaners. Allow the bag to air dry completely before reattaching or storing.
- **Inspection:** Periodically inspect the bag for any signs of wear, tear, or damage to the fabric or attachment loops.
- **Storage:** When not in use, store the bag in a clean, dry place away from direct sunlight and extreme temperatures.

6. TROUBLESHOOTING

If you encounter issues with your replacement bag, consider the following:

- **Bag Does Not Fit:** Double-check your walker's model number against the compatible models listed in Section 2. This bag is specifically designed for Helavo H1050 and H1150 series walkers.
- **Bag Sags or Feels Loose:** Ensure all attachment loops are correctly and securely fastened to the walker's frame. Re-thread them if necessary. Avoid overloading the bag.
- **Damage/Wear:** If the bag is torn or severely worn, it may need to be replaced to ensure safe and effective use.

7. SPECIFICATIONS

Product Name	Helavo Replacement Bag
Model Number	H1050-Serie (compatible with H1050 & H1150 walker series)
Dimensions	17.7 x 13.8 x 5.1 inches
Weight	7.41 ounces
Material	Alloy Steel
Color	Black
Manufacturer	Helavo

Note: The "Alloy Steel" material specification from the product data refers to the walker frame, not the bag itself, which is made of fabric.

8. WARRANTY AND SUPPORT

Helavo products are manufactured with German engineering design and service, ensuring quality and reliability. For specific warranty information regarding your replacement bag or your Helavo rolling walker, please refer to the documentation provided with your original walker purchase or visit the official Helavo website. If you require further assistance, have questions about compatibility, or need to contact customer service, please reach out to Helavo support. Contact details can typically be found on the Helavo website or through your purchase platform.