

Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

[manuals.plus](#) /

> [veicomtech](#) /

> [Veicomtech Ionic Foot Detox Machine User Manual](#)

veicomtech Veicomtech Ionic Foot Bath

Ionic Foot Detox Machine User Manual

Model: Veicomtech Ionic Foot Bath

PRODUCT OVERVIEW

The veicomtech Ionic Foot Detox Machine is designed to provide a relaxing and beneficial foot spa experience in the comfort of your home or salon. This system utilizes ionic technology to facilitate foot detoxification, aiming to enhance overall well-being. It includes a foot basin, controller unit, arrays, wrist band, power supply, FIR heating belt, and massage pads for a comprehensive experience.



Figure 1: Complete Veicomtech Ionic Foot Detox Machine with all included accessories.

SAFETY INFORMATION

For safe and effective operation, please adhere to the following safety guidelines:

- Ensure the metal plate of the wrist strap maintains constant contact with your skin during use.
- Avoid direct contact with the array or power supply during operation.
- Maintain an appropriate water level in the basin, ensuring it covers the array adequately but does not overflow.
- Do not use this product if you have certain health conditions, including but not limited to: pregnancy, epilepsy, pacemakers, organ transplants, or open wounds on the feet. Consult a healthcare professional if you have any concerns.
- Always follow the instructions provided in this user manual for safe and effective use.

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

PRODUCT COMPONENTS

The Veicomtech Ionic Foot Detox Machine package includes the following items:



Figure 2: All components included in the Veicomtech Ionic Foot Detox Machine package.

- 1 x Controller Unit
- 1 x Foot Basin
- 2 x Arrays
- 1 x Wrist Band
- 1 x Power Supply
- 1 x FIR Belt
- 2 x Massage Pads
- 10 x Ginger Tablets
- 25 x Basin Liners

- 1 x Instruction Manual
- 1 x Connect Line
- 1 x Gift Box

SETUP INSTRUCTIONS

Before operating your ionic foot detox machine, ensure all components are present and undamaged.

Controller Unit Overview




Figure 3: Detailed view of the controller unit's ports and buttons.

Familiarize yourself with the controller unit's interface:

- **Power ON/OFF:** Main power switch.
- **Foot Spa Timer:** Adjusts the duration of the foot spa session.
- **FIR Belt ON/OFF & Level Choose:** Activates and sets intensity for the FIR belt.
- **Message Intensity Choose:** Adjusts intensity for the massage pads.

- **Wrist Band Port:** Connects the wrist band.
- **Array Port:** Connects the ionic array.
- **LF Port:** Connects the massage pads.
- **DC Power Supply Feed In:** Connects the power adapter.
- **FIR Belt Port:** Connects the FIR heating belt.

OPERATING INSTRUCTIONS

Ionic Foot Spa Operation

Follow these steps to begin your ionic foot spa session:

1. **Prepare the Basin:** Add warm water into the foot basin. Ensure the water level is sufficient to immerse the array and your feet comfortably.
2. **Wear Wrist Strap:** Place the wrist strap on your wrist, ensuring the metal plate makes good contact with your skin.
3. **Connect Components:** Connect the wrist strap, array, and power supply to their respective ports on the controller panel.
4. **Position Feet and Array:** Place your feet and the array into the basin.
5. **Add Salt:** Add approximately 5g of sea salt or bath salt to the water.
6. **Power On:** Turn on the main power button. The display will show the default session time (e.g., 30 minutes).
7. **Set Time:** Press the "TIME" button to adjust the desired foot spa duration.
8. **Start Session:** Press the "DETOX" button to start the foot spa session. Press again to pause or stop.
9. **Add Ginger Tablet (Optional):** For enhanced benefits, you may add one ginger tablet to the basin after a few minutes. Observe the water changing color as the session progresses.



Figure 4: Step-by-step visual guide for operating the ionic foot spa.



Figure 5: Example of water changing color during an ionic foot spa session.

FIR Belt Operation

The FIR heating belt can be used concurrently with the foot spa for targeted relief:

1. **Wear Belt:** Securely wear the FIR belt on the desired body area (e.g., lower back).
2. **Connect Belt:** Connect the FIR belt to its designated port on the controller panel.
3. **Power On:** Turn on the main power button.
4. **Activate Belt:** Press the "BELT" button to activate the FIR heating function.
5. **Adjust Level:** Press the "BELT" button again to cycle through Low, Medium, and High heating levels.



Figure 6: Proper placement of the FIR heating belt for use.

Massage Pads Operation

The massage pads provide targeted electrical stimulation:

1. **Prepare Pads:** Detach the protective film from the white massage pads.
2. **Attach Pads:** Apply the pads to the target body position, ensuring good skin contact.
3. **Connect Pads:** Connect the pads to the "TENS" port on the controller panel.
4. **Power On:** Turn on the main power button.
5. **Activate Pads:** Press the "TENS" button to start or stop the massage function.
6. **Adjust Intensity:** Use the "+" and "-" buttons to adjust the intensity level of the massage.




Figure 7: User demonstrating the application of massage pads during a foot spa session.

MAINTENANCE

Proper maintenance ensures the longevity and optimal performance of your device:

- **Cleaning the Basin:** After each use, empty the basin and rinse it thoroughly with clean water. Wipe it dry with a soft cloth.
- **Cleaning the Array:** After each session, rinse the array under running water to remove any residue. For deeper cleaning, soak the array in a mixture of 50% white vinegar and 50% water for 10-15 minutes, then scrub gently with a brush and rinse. Ensure it is completely dry before storage.
- **Cleaning the Controller Unit:** Wipe the controller unit with a dry or slightly damp cloth. Do not immerse it in water or use abrasive cleaners.
- **Storage:** Store all components in a cool, dry place away from direct sunlight and extreme temperatures.
- **Liner Replacement:** Use a new basin liner for each session to maintain hygiene and simplify cleanup.

TROUBLESHOOTING

If you encounter issues with your Veicomtech Ionic Foot Detox Machine, refer to the following common problems and solutions:

Problem	Possible Cause	Solution
Machine does not power on.	Power cable not connected properly; Power outlet issue; Unit malfunction.	Ensure power cable is securely plugged into the unit and a working outlet. Try a different outlet. If problem persists, contact customer support.
No bubbles or activity from the array.	Insufficient salt in water; Array not properly connected; Array worn out; Wrist strap not making good contact.	Add more salt to the water (ensure it's dissolved). Check array connection. Ensure wrist strap metal plate is firmly against skin. Replace array if it's old or damaged.
Water does not change color.	Insufficient salt; Array not functioning; Short session time.	Ensure adequate salt is added. Verify array is bubbling. Extend session time if needed. Note that water color change is also influenced by water impurities and array oxidation.
FIR Belt or Massage Pads not working.	Not properly connected; Not activated on controller; Low intensity setting.	Check connections to the controller unit. Press the respective "BELT" or "TENS" button to activate. Increase intensity using the level/intensity buttons.

SPECIFICATIONS

Key technical specifications for the Veicomtech Ionic Foot Detox Machine:

Feature	Detail
Model Number	Veicomtech Ionic Foot Bath
Manufacturer	Veicomtech
Power Source	Corded Electric
Material	Metal (for certain components)
Color	Red 801m (Controller Unit)
Package Dimensions	17.48 x 13.7 x 6.81 inches
Item Weight	8.04 Pounds

WHAT IS IN THE BOX

10 x Ginger Tablets

1 x Gift Box

1 x Controller Unit

1 x 25 Tub Liners

1 x Wrist Strap

2 x White Pads

1 x Instruction

1 x Basin

1 x Power Supply

2 x Arrays

1 x FIR Waist Belt

1 x Connect Line



Figure 8: Approximate dimensions of the foot spa basin and array.

WARRANTY AND SUPPORT

Veicomtech stands by the quality of its products. The controller unit of this Ionic Foot Detox Machine comes with a **12-month warranty** from the date of purchase.

For any questions, technical assistance, or warranty claims, please contact Veicomtech customer support. Refer to the contact information provided on the product packaging or the official Veicomtech website for the most up-to-date support details.