

JTF P200BH

JTF P200 Sonic Toothbrush Replacement Heads User Manual

Model: P200BH

INTRODUCTION

This manual provides essential information for the proper use and maintenance of your JTF P200 Sonic Toothbrush Replacement Heads. These standard brush heads are designed for compatibility with JTF P200 series rechargeable electric toothbrush handles, ensuring effective oral hygiene.

Regular replacement of toothbrush heads is crucial for maintaining optimal cleaning performance and hygiene. Dentists recommend replacing brush heads every 3 months.

COMPATIBILITY

These replacement heads are specifically designed for use with the JTF P200 Sonic Toothbrush handle. Ensure your toothbrush handle is a P200 series model for proper fit and function.

Compatible with JTF Brand P200 Sonic Handle Only

White / Black / Pink



Image: JTF P200 Sonic Toothbrush handles in various colors, illustrating compatibility with the replacement heads.

SETUP AND INSTALLATION

Follow these steps to replace your JTF P200 Sonic Toothbrush head:

1. **Detach the Old Brush Head:** Gently pull the old brush head straight off the toothbrush handle. Do not twist.
2. **Install the New Brush Head:** Align the new replacement brush head with the metal shaft on the toothbrush handle.
3. **Secure the Head:** Press the new brush head down firmly onto the shaft until it clicks into place. Ensure it is securely attached and there is no gap between the head and the handle.



Image: Step-by-step guide for detaching and installing a JTF P200 replacement brush head.

OPERATING RECOMMENDATIONS

For optimal cleaning performance and oral health, follow these recommendations:

- **Brush Head Replacement:** Dentists recommend replacing your toothbrush head every 3 months, or sooner if the bristles become frayed or worn.
- **Brushing Technique:** Use the JTF P200 sonic toothbrush as directed in its main instruction manual. The replacement heads feature W-shaped bristles designed to adapt to the shape of your teeth for effective cleaning.
- **Bristle Hardness:** These are medium hardness brush heads, providing effective plaque control and removal.

W-shaped Bristles better adapt to the shape of your teeth



Image: Illustration of the W-shaped bristles on the JTF P200 replacement head, designed for thorough tooth coverage.

40000 VPM

Sonic Tech remove up to 100 % more plaque



Image: Visual representation of the sonic technology, demonstrating 40,000 vibrations per minute for plaque removal.

MAINTENANCE AND CARE

- **Rinsing:** After each use, rinse the brush head thoroughly under running water to remove toothpaste residue and debris.
- **Drying:** Allow the brush head to air dry.
- **Hygienic Cap:** Each replacement head comes with a hygienic cap. Use this cap for storage and travel to protect the bristles from dust and bacteria.
- **Cleaning the Handle:** Refer to your JTF P200 sonic toothbrush manual for instructions on cleaning the handle.



Image: Detail of the medium soft bristles, emphasizing their design for effective cleaning.

TROUBLESHOOTING

- **Brush Head Does Not Fit:** Ensure you have a genuine JTF P200 series toothbrush handle. These replacement heads are designed specifically for the P200 model. Do not force the head onto the handle.
- **Reduced Cleaning Performance:** If you notice reduced cleaning effectiveness, it may be time to replace your brush head. Bristles can wear down over time, typically after 3 months of use.
- **Loose Brush Head:** Ensure the brush head is pressed down firmly onto the metal shaft until it clicks into place. If it remains loose, inspect both the brush head and the handle for any damage.

SPECIFICATIONS

Brand	JTF
Model Number	P200BH
Compatibility	JTF P200 Sonic Toothbrush Handles
Bristle Firmness	Medium
Color	White
Number of Pieces	3 per pack
Recommended Replacement	Every 3 months
Package Dimensions	3.74 x 2.32 x 1.14 inches
Item Weight	2.4 ounces
UPC	761736230443



Image: Detailed dimensions of the JTF P200 replacement brush head.

OFFICIAL PRODUCT VIDEO

Your browser does not support the video tag.

Video: An overview of the JTF P200 Sonic Electric Toothbrush, demonstrating its features and various color options. This video provides context for the compatible replacement heads.

WARRANTY AND SUPPORT

For information regarding warranty coverage or customer support, please refer to the documentation included with your JTF P200 Sonic Toothbrush handle or contact JTF customer service directly. The standard return policy for this product is 30 days for refund or replacement.

Manufacturer: JTF

Seller: Greek-US

