

[Manuals.plus](#) /

› [babytone](#) /

› Babytone Baby Sleep Monitor S1 & S2 Spare Fabric Sock Instruction Manual

babytone S1, S2

Babytone Baby Sleep Monitor S1 & S2 Spare Fabric Sock Instruction Manual

Model: S1, S2 | Brand: babytone

1. PRODUCT OVERVIEW

The Babytone Baby Sleep Monitor Spare Fabric Socks are designed as a comfortable and skin-friendly replacement for the original babyO2 straps. These socks ensure stable and accurate data measurement for your baby's sleep monitor. They are made from soft nylon fabric and come in two sizes to accommodate different foot circumferences.




Image 1.1: Overview of the Babytone S1 and S2 spare fabric socks, illustrating the two sizes (S and L) and their design.

2. KEY FEATURES

- **Comfortable and Skin-Friendly:** Made from soft nylon fabric, ensuring comfort and safety for your baby's delicate skin.
- **Enhanced Stability:** Designed to securely hold the baby sleep monitor's sensors in place, providing more stable data measurement compared to traditional straps.
- **Two Sizes Included:** Each package contains two large (L) and two small (S) sizes to fit a range of baby foot sizes.
- **Size Specifications:**
 - Size S: Suitable for foot width of 80 to 120mm (3-4.7 inches), typically for babies weighing less than 12 lbs.
 - Size L: Suitable for foot width of 100 to 153mm (3.9-6 inches), typically for babies weighing 12-30 lbs.
- **Recommended Foot Circumference:** Designed for feet with a circumference of 3-6 inches (80-150mm).

Material Comfortable and Size Fit

The sock is with soft nylon material and two kinds of sizes, perfect designed for babies from little to big



Image 2.1: Illustration of the comfortable material and size fit, showing a baby wearing the sock and the size guide.

3. SETUP AND INSTALLATION

Follow these steps to properly set up the Babytone spare fabric sock with your baby sleep monitor:

1. **Choose the Correct Size:** Select either the Small (S) or Large (L) sock based on your baby's foot width and weight. Refer to the size specifications in Section 2.

2. **Measure Foot Circumference:** For optimal fit, measure your baby's foot circumference to ensure it falls within the recommended 3-6 inch range.
3. **Insert the Monitor Sensor:** Carefully place the Babytone sleep monitor sensor into the designated pocket within the fabric sock. Ensure it is snug and secure.
4. **Position on Baby's Foot:** Gently slide the sock onto your baby's foot, ensuring the sensor is positioned correctly on the top or side of the foot as per your monitor's instructions.
5. **Secure the Sock:** Use the integrated hook-and-loop fasteners to secure the sock comfortably around your baby's foot. Do not overtighten; it should be snug but not restrictive.

Fit 3-6 inch foot circumference





Image 3.1: Guide for measuring foot circumference and ensuring proper fit of the sock.



Image 3.2: A baby's foot with the Babytone sleep monitor sock and sensor properly installed.

4. OPERATING INSTRUCTIONS

The fabric socks are a component of the Babytone Baby Sleep Monitor system. Once the sock is properly installed with the monitor sensor, operate your Babytone Sleep Monitor according to its specific user manual. The sock's primary function is to provide a secure and comfortable fit for the sensor, ensuring continuous and accurate data collection.

- Ensure the monitor is powered on and connected to its respective application or display unit.
- Verify that the sensor is making proper contact with the baby's skin through the sock.
- Monitor data as usual through the Babytone app or display.

5. MAINTENANCE AND CARE

Proper care of your Babytone fabric socks will extend their lifespan and maintain hygiene:

- **Cleaning:** Hand wash the fabric socks gently with mild soap and cold water. Do not use harsh detergents or bleach.
- **Drying:** Air dry completely. Do not tumble dry or iron, as high heat can damage the fabric and hook-and-loop fasteners.
- **Storage:** Store in a clean, dry place away from direct sunlight when not in use.
- **Inspection:** Regularly inspect the socks for any signs of wear, tear, or damage. Replace if necessary to ensure optimal performance and safety.

6. TROUBLESHOOTING

If you encounter issues with the fabric socks, consider the following:

Problem	Possible Cause	Solution
Sensor data is inconsistent or inaccurate.	Sock is too loose or too tight; sensor not properly seated; incorrect sock size.	Adjust the sock for a snug but comfortable fit. Ensure the sensor is fully inserted. Try a different size sock if available.
Baby shows discomfort or irritation.	Sock is too tight; material irritation (rare); prolonged wear without repositioning.	Loosen the sock slightly. Ensure the sock is clean. Reposition the sock periodically. If irritation persists, discontinue use and consult a pediatrician.
Sock appears worn or damaged.	Normal wear and tear; improper cleaning/drying.	Replace the sock with a new one. Follow maintenance instructions carefully.

7. SPECIFICATIONS

Brand: babytone

Model Compatibility: Babytone Baby Sleep Monitor S1 & S2

Material: Nylon fabric

Color: Blue

Sizes Included: 2 x Small (S), 2 x Large (L)

Size S Foot Width: 80-120mm (3-4.7 inches)

Size L Foot Width: 100-153mm (3.9-6 inches)

Recommended Foot Circumference: 80-150mm (3-6 inches)

Item Package Dimensions: 4.09 x 2.8 x 1.26 inches

Package Weight: 0.04 Kilograms

UPC: 850040046724

ASIN: B0BJDRHDRB

8. WARRANTY AND SUPPORT

For warranty information and customer support regarding your Babytone products, please refer to the original Babytone Baby Sleep Monitor manual or visit the official Babytone website. You can also contact Babytone customer service directly for assistance with spare parts or product inquiries.

Babytone Official Store: [Visit the babytone Store on Amazon](#)



© 2024 babytone. All rights reserved.

This manual is for informational purposes only. Specifications are subject to change without notice.