

## Manuals+

[Q & A](#) | [Deep Search](#) | [Upload](#)

manuals.plus /

> [Pride Mobility](#) /

> [Pride Mobility Go-Go Elite Traveller 3-Wheel SC40E 18AH Mobility Scooter User Manual](#)

## Pride Mobility SC40E

# Pride Mobility Go-Go Elite Traveller 3-Wheel SC40E 18AH Mobility Scooter User Manual

Model: SC40E | Brand: Pride Mobility

## 1. INTRODUCTION

This manual provides essential information for the safe and efficient operation, maintenance, and care of your Pride Mobility Go-Go Elite Traveller 3-Wheel SC40E 18AH mobility scooter. The compact design allows for easy maneuverability in tight spaces while providing stable outdoor performance. Key features include feather-touch disassembly for transport and a delta tiller with a convenient charger port.

The Go-Go Elite Traveller is a FDA Class II Medical Device, designed to assist individuals with mobility impairments. Please read this manual thoroughly before operating the scooter to ensure proper usage and safety.

## 2. SAFETY INFORMATION

Always operate the scooter in a safe manner. Ensure you are familiar with all controls and functions before use. Avoid operating the scooter in conditions that may compromise stability or control. Always use caution when navigating inclines or uneven terrain.

## 3. PRODUCT OVERVIEW

The Go-Go Elite Traveller 3-Wheel SC40E 18AH mobility scooter is designed for both indoor and outdoor use, offering a blend of portability and performance. It features a 300 lb. weight capacity and a maximum speed of up to 4 mph.

- **Feather-Touch Disassembly:** Allows for simple frame separation into five lightweight pieces for easy transport and storage.
- **Delta Tiller with Charger Port:** Provides convenient charging options, including on-board or off-board charging with the dual voltage charger.
- **Microprocessor-Based Controller:** Optimizes power usage and enhances safety during operation.
- **Non-Scuffing Pride Tires:** Designed for durability and a long-lasting ride.
- **Front Frame-Mounted Seat Post:** Enhances stability during rides.



Image: The Pride Mobility Go-Go Elite Traveller 3-Wheel SC40E 18AH mobility scooter in blue, showcasing its compact design and front basket.

# THE QUICKEST WAY TO TRAVEL



Image: A detailed diagram highlighting key features such as the compact, foldable seat with armrests, feather-touch disassembly, easy-to-grip tiller adjustment knob, standard front basket, and interchangeable shroud panels.

## 4. SETUP

The Go-Go Elite Traveller is designed for easy assembly and disassembly. The feather-touch disassembly feature allows the scooter to break down into five super lightweight pieces for convenient transport and storage. To assemble, simply connect the main body sections, attach the seat, and secure the tiller. Ensure all connections are firm before use.

Video: This video demonstrates the assembly and key features of a Go-Go Ultra X Scooter, which shares similar design principles and functionality with the Elite Traveller model, including feather-touch disassembly and compact storage.

## 5. OPERATING INSTRUCTIONS

---

To operate your scooter, ensure it is fully charged and the key is in the 'ON' position. Adjust the seat and tiller for optimal fit and comfort. The scooter features intuitive controls for forward and reverse movement, typically located on the tiller. The speed can be adjusted up to 4 mph. The compact turning radius allows for easy navigation in various environments.

### Charging

The delta tiller includes a charger port for convenient charging. The dual voltage charger allows for both on-board and off-board charging. Ensure the scooter is turned off during charging. Refer to the user manual for detailed instructions on battery charging procedures.

Video: This video provides an overview of mobility scooter features, including controls, charging, and general operation, which are applicable to the Go-Go Elite Traveller model.

## 6. MAINTENANCE

---

Regular maintenance is crucial for the longevity and performance of your scooter. Clean the scooter regularly to prevent dirt and debris buildup. Conduct routine checks on the battery and tires to ensure they are in good condition. For detailed maintenance schedules and procedures, please refer to the comprehensive user manual provided with your product.

## 7. TROUBLESHOOTING

---

If you encounter any issues with your scooter, please refer to the troubleshooting section in the full user manual. Common issues may include battery performance or control responsiveness. If problems persist, contact customer support for assistance.

## 8. SPECIFICATIONS

---

Feature	Specification
Model Name	Go-Go Elite Traveller - 3 Wheel Travel Scooter
Material	Aluminum
Style	Compact, Electric
Color	Red, Blue (Shroud Inserts)
Brand	Pride Mobility
ASIN	B0BNSMSZMK
Manufacturer	Pride Mobility Products Corporation
Date First Available	December 1, 2022
Weight Capacity	300 lbs.
Maximum Speed	Up to 4 MPH
Range Per Charge	Up to 10.5 miles (18AH)

## 9. WARRANTY AND SUPPORT

---

The Pride Mobility Go-Go Elite Traveller 3-Wheel SC40E 18AH mobility scooter comes with a 2-Year Warranty on the Frame, 2-Year Warranty on Electronics & Drivetrain, and a 6-Month Warranty on the Battery. For specific warranty details, questions, or exclusions, please contact the seller or Pride Mobility Products Corporation directly.

For further assistance or support, please contact Pride Mobility Products Corporation.