

[Manuals.plus](#) /

› [FACEMOON](#) /

› FACEMOON KN95 Face Mask Instruction Manual

FACEMOON KN95 Face Mask

FACEMOON KN95 Face Mask Instruction Manual

1. INTRODUCTION TO YOUR FACEMOON KN95 FACE MASK

Thank you for choosing the FACEMOON KN95 Face Mask. This manual provides essential information for the proper use, care, and understanding of your disposable protective mask. Designed for adults, these masks offer 5-layer protection and are engineered for breathability and comfort.



Each FACEMOON KN95 mask features a robust 5-layer design, providing enhanced filtration. The construction includes two outer layers of non-woven fabric, two inner layers of melt-blown filters, and a middle layer of non-woven cotton. This design aims to maximize protection while maintaining optimal breathability.

30% INCREASE IN THICKNESS



Image 1.2: A visual comparison highlighting the increased thickness of FACEMOON KN95 masks for improved protection.

2. SETUP: HOW TO PROPERLY WEAR YOUR KN95 MASK

For effective protection and comfort, it is crucial to wear your KN95 mask correctly. Follow these steps to ensure a secure fit:

1. **Unfold the Mask:** Gently unfold the mask to its full shape.
2. **Position the Mask:** Place the mask over your nose and mouth, ensuring the nose clip is at the top.
3. **Secure Ear Loops:** Place the ear loops around each ear.
4. **Adjust Nose Clip:** Press the flexible nose clip firmly against the bridge of your nose to form a tight seal.
5. **Check for Gaps:** Ensure there are no gaps around the edges of the mask. Adjust the mask and ear loops as needed.
6. **Utilize Adjustable Clips (Optional):** If the ear loops cause discomfort or the fit is not snug enough, use the provided adjustable clips. Attach the ear loops to the clips behind your head and select one of the four settings on each side to achieve a customized, comfortable, and secure fit. This helps reduce ear pressure.

FACEMOON KN95 MASK



Image 2.1: A woman demonstrating the proper fit of a FACEMOON KN95 mask, with key features labeled.

GB2626-2019 TEST RESULTS

KN95 MEETS INTERNATIONAL STANDARDS

Tests Performed	Test Stanaaras	Test Results	Comment
Particle Filtration Efficiency	GB2626-2019	KN95 $\geq$ 95%	Pass

Image 2.2: Diagram illustrating the use of adjustable ear clips to achieve a personalized and comfortable mask fit.

3. OPERATING: INTENDED USE AND PROTECTION

The FACEMOON KN95 mask is designed to provide respiratory protection in various environments. Its 5-layer filtration system is GB2626-certified, ensuring a particle filtration efficiency of 95% or greater.

5-LAYERS FILTER PROTECTION

≥ 95%PFE
Filtration




Image 3.1: Detailed diagram of the 5-layer filtration system, showcasing each protective layer.

This mask is suitable for use in conditions involving:

- Wildfire smoke
- Smog
- Dust (e.g., during grinding, sanding, sawing, sweeping)
- Urban pollution

PERSONAL PROTECTIVE EQUIPMENT



SMOG WEATHER



DUSTWORKING



URBAN POLLUTION



WILDFIRE SMOKE

Image 3.2: Examples of environments where the KN95 mask offers personal protective equipment.

The mask's design prioritizes breathability, allowing for comfortable wear over extended periods.

MORE BREATHABLE

FACEMOON KN95 MASK



OTHERS KN95 MASK



Image 3.3: A visual demonstration of the mask's breathability.

4. MAINTENANCE: CARE AND STORAGE

FACEMOON KN95 masks are designed for single use and are individually packaged to ensure hygiene and maintain cleanliness until use. They are odorless and compact, making them convenient for storage and transport.



INDIVIDUALLY PACKAGED CLEAN AND ODORLESS

Image 4.1: An individually packaged FACEMOON KN95 mask, highlighting its hygienic storage.

- **Storage:** Store unused masks in their individual packaging in a cool, dry place away from direct sunlight and contaminants.
- **Disposal:** Dispose of used masks responsibly according to local regulations. Do not reuse disposable masks.

5. TROUBLESHOOTING

While FACEMOON KN95 masks are designed for comfort and a secure fit, some common concerns may arise:

- **Discomfort Around Ears:** If the ear loops cause discomfort, utilize the adjustable clips provided with your masks. These clips allow you to adjust the tension and shift pressure from behind the ears to the back of the head.
- **Mask Feels Loose/Tight:** Ensure the nose clip is properly molded to your nose. If the mask still feels too loose, adjust the ear loops using the adjustable clips. If it feels too tight, try a looser setting on the clips.
- **Glasses Fogging:** A common issue, often indicating a poor seal around the nose. Re-adjust the nose clip to ensure it is pressed firmly against your nose, preventing air from escaping upwards.

6. SPECIFICATIONS

Feature	Detail
Brand	FACEMOON
Filter Class	KN95 (GB2626-2019 certified)
Layers	5-Ply (Non-woven fabric, Melt-blown filters, Non-woven cotton)
Quantity	100 Count (Pack of 100)
Color	Colored assortment
Department	Unisex Adult
Manufacturer	Changzhoushi Kangmin Medical Equipment Co.,Ltd.
ASIN	B0BW41M4CJ

Image 6.1: Official GB2626-2019 test results confirming the mask's filtration performance.

7. IMPORTANT SAFETY INFORMATION

Please note the following critical information regarding your FACEMOON KN95 Face Mask:

- This KN95 mask is **NOT an N95 mask**. While both offer filtration, N95 masks are certified by NIOSH in the United States, and KN95 masks are certified under Chinese standards (GB2626-2019).
- For more information about KN95 masks and their differences from N95 masks, please refer to [amazon.com/AboutKN95s](https://www.amazon.com/AboutKN95s).
- Do not use this mask in environments with insufficient oxygen, unknown contaminants, or immediately dangerous to life or health (IDLH) conditions.
- Do not use this mask for gas or vapor protection.
- If the mask becomes damaged, soiled, or breathing becomes difficult, discard it and replace it with a new one.

8. WARRANTY AND SUPPORT

Information regarding product warranty and customer support is not available in the provided product data. For any inquiries, please refer to the seller's contact information on the purchase platform.

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question