

AIFINI X2

AIFINI Ear Wax Removal Tool Camera X2 User Manual

Model: X2 | Brand: AIFINI

INTRODUCTION

This manual provides essential information for the safe and effective use of your AIFINI Ear Wax Removal Tool Camera X2. Please read these instructions carefully before operation and retain for future reference.

SAFETY INFORMATION

- Do not use the product while charging.
- Ensure the silicone sleeve on the ear spoon is tightened before each use to prevent scratching the ear canal.
- Move the otoscope gently and slowly within the ear canal.
- Consult a medical professional if you experience pain, discomfort, or have pre-existing ear conditions.
- Keep out of reach of children.

PACKAGE CONTENTS

The AIFINI Ear Wax Removal Tool Camera X2 package includes:

- Otoscope with light (1)
- Type-C charging cable (1)
- Silicone ear spoons (6)
- User Manual (1)
- 7-piece ear wax removal tool kit (1)



Figure 1: Complete Package Contents

PRODUCT OVERVIEW

The AIFINI Ear Wax Removal Tool Camera X2 features a 1080P HD camera and bright LED lights for clear visualization of the ear canal. Its all-in-one design securely integrates the camera and ear spoon, enhancing safety during use. The device is equipped with a Type-C charging port and offers stable Wi-Fi connectivity for real-time viewing on your smartphone or tablet.

1080P HD Ear Wax Removal Tool

Opening up new horizons in ear
wax removal



Easy to fall off and
fall in the ear



Ear canal fluid
accumulation, prone
to infection



Camera and ear
spoon in one

Figure 2: 1080P HD Camera and LED Lights

Earwax Removal Tool Features and Details



Figure 3: Key Features of the Earwax Removal Tool

SETUP INSTRUCTIONS

1. **Download the App:** Scan the QR code provided in the manual or visit <http://qr19.cn/CRUIhc> to download the "Suear" app. The app is compatible with both Android and iOS devices.
2. **Power On:** Press and hold the ear wax cleaner on/off button for 1 to 2 seconds to turn on the device.
3. **Connect to Wi-Fi:** On your smartphone or tablet, go to your Wi-Fi settings and connect to the Wi-Fi network named "Suear -xxxx" (where xxxx is a unique identifier). It is normal for a "Network unavailable" prompt to appear, as the device does not support external internet networking.
4. **Start Using:** Return to the "Suear" app. The real-time camera feed from the otoscope should now be displayed.



Figure 4: Product Connection Guide

OPERATING INSTRUCTIONS

Once connected, you can view the inside of your ear canal in real-time on your device's screen. The otoscope features a six-axis gravity sensor and screen reverse design, allowing for flexible viewing from various positions (left ear, right ear, sitting, reclining).

- Gently insert the ear spoon into the ear canal.
- Use the real-time view on your phone to guide the ear spoon and carefully remove earwax.
- The kit includes five highly flexible interchangeable silicone ear spoons with three kinds of curved designs (low arc, middle arc, high arc) to comfortably fit and clean the ear.
- Additional 7-piece ear wax removal tools are provided for various cleaning needs, including spring-loaded ear picks for cleaning and massaging, and tools for rotating to clean.

Six-Axis Gravity Sensor & Screen Reverse Design

Left/right ear view switch at will



Figure 5: Flexible Viewing Angles

Ear Cleaning Kit with Six Silicone Ear Spoons

(3 kinds of curves)



Figure 6: Silicone Ear Spoons and Design

MAINTENANCE

- After each use, clean the ear spoons and the camera tip thoroughly.
- Use the provided cleaning brush to remove any debris from the tools.
- The camera lens is waterproof (IP67), allowing for easy cleaning.
- Store all components in a clean, dry place.

TROUBLESHOOTING

- **Device not turning on:** Long press the ON/OFF button for 1-2 seconds. Ensure the device is charged.
- **No image on app:** Confirm the device is powered on and connected to the "Suear -xxxx" Wi-Fi network. Restart the app and device if necessary.
- **"Network unavailable" prompt:** This is normal. The device creates its own Wi-Fi network for direct connection to your phone and does not provide internet access.

SPECIFICATIONS

Feature	Detail
Item Model Number	X2
Product Dimensions	6.69 x 2.55 x 0.7 inches
Item Weight	2.4 ounces
Batteries	1 Lithium Polymer battery (included)
Manufacturer	AIFINI
Camera Resolution	1080P HD
Waterproof Rating	IP67

OFFICIAL PRODUCT VIDEOS

AIFINI Ear Wax Removal Tool Review

Your browser does not support the video tag.

This video provides a detailed review of the AIFINI Ear Wax Removal Tool, demonstrating its features and functionality for effective ear cleaning.

AIFINI Ear Wax Removal Tool Comparison

Your browser does not support the video tag.

This video offers a comparison between the AIFINI Ear Wax Removal Tool and a similar product, highlighting differences in features and user experience.

WARRANTY AND SUPPORT

For warranty information and customer support, please refer to the contact details provided in the original product packaging or visit the official AIFINI website.

© 2023 AIFINI. All rights reserved.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page’s manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for

conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question